RIDA[®] C18 column

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RIDA[®] C18 column

1. General

RIDA[®] C18 column (R2002) are C18 columns for solid phase extraction. The solid phase extraction (SPE) is an efficient technique for a fast and selective sample preparation, especially for the RIDASCREEN[®] anabolic and hormone enzyme immunoassays.

The solid phase extraction offers a lot of advantages compared with other sample preparation techniques:

- high recovery rates
- sample concentration
- high clean up
- easy automation
- compatibility with instrumental analysis
- reduction of solvent volumes

2. Reagents provided

100 x C18 columns (100 mg/ml column reservoir)

The columns are packed with 100 mg end-capped sorbent (octadecyl-sorbent).

This sorbent is manufactured using trifunctional silane under carefully controlled conditions.

The columns have been subjected to a number of quality control tests (see quality assurance certificate).

3. Principle of solid phase extraction

Sample analytes are retained, when the fluid sample passes the column.

In the interference elution step an appropriate wash solution is indicated.

For the following selective elution of the analyte an appropriate elution solvent is needed.

Concentration of the analytes is often possible by the extraction procedure.

A solid phase extraction procedure normally consists of 5 steps (see Fig. 1):

- sample pre-treatment (extraction and conditioning of the extract)
- column solvation and pre-equilibration (1)
- sample application (2)
- elution of interfering matrix components / rinsing of the column (3)
- elution of the analyte (4)

The columns are packed with sorbent of an average particle size of 140 μm .

Most solvents pass the column without technical help.

Viscous samples and solvents can be passed through the column using either a vacuum, centrifuge or positive pressure, shown in Fig. 2.


1. Column solvation / pre-equilibration

2. Apply sample

3. Elution of interfering matrix components (interferences) rinsing step


4. Analyte elution


● = Analyte
▲ }
□ } Interferences
△ }


Fig. 1: Typical steps for solid phase extraction

Vacuum Manifold

Centrifugation

Positive pressure
(manual)

Fig. 2: Possible techniques for solid phase extraction

4. Handling of RIDA® C18 column

The flow rate control during sample loading, interference and analyte elution steps are important.

Too fast flow rates can result in low recoveries due to breakthrough during the retention step or inadequate elution during the elution step.

Too high flow rates during the interference elution step can result in dirty extracts.

A flow rate of 1 - 2 ml per min is recommended. (Usually working steps are possible under hydrostatic conditions.)

Flow rates should be specific in all protocols to ensure trouble free method transfers from one analyst or laboratory to another and reproducibility.

Please note, that the presence of significant concentrations of particulate matter can result in clogging of the inlet frit. Filtration or centrifugation of a sample to remove particulate matter can often enhance the reliability of an SPE procedure. However, care must be taken to ensure the analyte is not removed with the particulate matter. Pre-treatment of the sample (addition of methanol, pH adjustment etc.) may be required to minimize analyte-matrix interaction.

Caution:

Drying the column prior to application of the sample can result in analyte breakthrough and variable recoveries.

If the column dries out before loading the sample, simply repeat the solvation and pre-equilibration procedure.

5. Sample preparation

See product information of the respective RIDASCREEN® enzyme immunoassay.

Several other SPE procedures are possible with the RIDA® C18 columns. Please contact your local distributor for more information.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG