


UV method for approx. 16 assays each

 For laboratory use only  
 Store between +2 and +8°C

The method is contained in the Austrian, Dutch, German, Swiss food laws. Recommended e. g. by IFU, AIJN, MEBAK, OICCC, VDLUFA. Standardized by DIN, EN, GOST, NEN, NF.

## Principle



Ref.: Bergmeyer, H.U. & Bernt, E. (1974) in Methods of Enzymatic Analysis (Bergmeyer, H.U., ed.) 2nd ed., vol. 3, pp. 1176-1179; Bergmeyer, H.U., Bernt, E., Schmidt, F. & Stork, H. (1974) in Methods of Enzymatic Analysis (Bergmeyer, H.U., ed.) 2nd ed., vol. 3, pp. 1196-1201; Verlag Chemie Weinheim / Academic Press Inc., New York and London.

## Assay performance

Wavelength: 340 nm (NADPH),  $\epsilon = 6.3 \text{ l} \times \text{mmol}^{-1} \times \text{cm}^{-1}$   
 Light path: 1.00 cm (glass or plastic cuvettes)  
 Temperature: +20 to +25°C  
 Assay volume: 3.020 ml  
 Measurement: against air or against water  
 Sample solution: 4 to 150 µg sucrose + D-glucose in 0.100 to 1.800 ml (sucrose), resp. 2.000 ml (D-glucose) sample solution.

## Reagents

- # S: Lyophilizate with citrate buffer, pH approx. 4.6, approx. 510 U of  $\beta$ -fructosidase (for stability see pack label). Dissolve contents of bottle #S with 7 ml redist. water. The solution is stable for 1 month at +2 to +8°C, resp. for 2 months at -15 to -25°C.
- # 1: Powder mixture with triethanolamine buffer, pH approx. 7.6, approx. 80 mg NADP, approx. 190 mg ATP, magnesium sulfate (for stability see pack label). Dissolve contents of bottle # 1 with approx. 31 ml redist. water. The solution is stable for 1 month at +2 to +8°C, resp. for 2 months at -15 to -25°C.
- # 2: Approx. 0.7 ml hexokinase (HK) / glucose-6-phosphat dehydrogenase (G6P-DH) suspension (approx. 200 U / 100 U) in ammonium sulfate (for stability see pack label). The suspension is ready for use. Swirl bottle carefully before the suspension is pipetted.

### In addition (not contained in the kit):

Standard solution sucrose, ultrapure, 0.8 g/l for test control only; standard solution D-glucose, anhydrous, ultrapure, 0.5 g/l for test control only.

The reagents for the determination of sucrose and D-glucose are not hazardous. The general safety rules for the work in chemical laboratories should be applied. After use the reagents can be disposed of with the laboratory waste. Packaging materials may be recycled.

## Procedure

Pipette into cuvettes:	Blank sucrose	Sucrose standard assay <sup>1</sup>	Sucrose sample assay <sup>2</sup>	Blank D-glucose	D-glucose sample assay	Assay with internal standard <sup>3</sup>
Citrate buffer, $\beta$ -fructosidase, solution # S	0.200 ml	0.200 ml	0.200 ml	-	-	0.200 ml
<b>Sample solution<sup>4</sup></b> (e.g. 0.08 to 0.8 g sucrose/l)	-	-	<b>0.100 ml</b>	-	<b>0.100 ml</b>	<b>0.100 ml</b>
Standard solution <sup>4</sup> (e.g. 0.8 g sucrose/l)	-	0.100 ml	-	-	-	0.100 ml
<b>Mix, e.g. by gentle swirling of the cuvette. Incubate at +20°C to +25°C (+37°C) for 15 min (5 min). Add :</b>						
Tea buffer, NADP and ATP solution # 1	1.000 ml	1.000 ml	1.000 ml	1.000 ml	1.000 ml	1.000 ml
Redist. Water	1.800 ml	1.700 ml	1.700 ml	2.000 ml	1.900 ml	1.600 ml
<b>Mix<sup>5</sup>, after approx. 3 min read the absorbances (A<sub>1</sub>). Add:</b>						
HK/G6P-DH suspension # 2	0.020 ml	0.020 ml	0.020 ml	0.020 ml	0.020 ml	0.020 ml
<b>Mix<sup>5</sup>, after approx. 10 to 15 min read the absorbances (A<sub>2</sub>). Repeat absorbance reading after another 2 min<sup>6</sup>.</b>						

## Notes

- Run a "standard" to see "accidents" in analysis. The measurement of the standard is not necessary for calculating results.
- This assay together with the blank is a single determination. In the case of a double determination, run two assays with different sample volumes. The absorbance differences measured have to be proportional to the sample volumes (calculate with resp. volumes).
- Recovery =  $[(\Delta A_{\text{sample+standard}} - \Delta A_{\text{sample}}) / \Delta A_{\text{standard}}] \times 100 [\%]$ .
- Before dispensing, rinse the enzyme pipette, resp. the tip of the piston pipette with sample resp. with standard solution.
- e.g. with a plastic spatula, or after closing the cuvette with Parafilm (trademark of American Can Co., Greenwich Ct., USA)
- The reaction has stopped when the absorbance is constant. If the reaction has not stopped, continue to read absorbances until they increase constantly over e. g. 2 min. Extrapolate absorbances to the time of the addition of HK/G6P-DH (suspension # 2).

**Calculation**

$$\Delta A_{D\text{-glucose}} = (A_2 - A_1)_{\text{sample, resp. standard}} - (A_2 - A_1)_{\text{blank}}$$

$$\Delta A_{\text{sucrose}} = \Delta A_{\text{test sucrose}} - \Delta A_{\text{test D-glucose}}$$

$$= [(A_2 - A_1)_{\text{assay sucrose}} - (A_2 - A_1)_{\text{blank sucrose}}] - [(A_2 - A_1)_{\text{assay D-glucose}} - (A_2 - A_1)_{\text{blank D-glucose}}]$$

$$c = (V \times MW \times \Delta A) / (\epsilon \times d \times v \times 1000) \text{ [g D-glucose, resp. sucrose /l sample solution]}$$

$$c = (3.020 \times 180.16 \times \Delta A) / (6.3 \times 1.00 \times 0.100 \times 1000) = \mathbf{0.8636 \times \Delta A \text{ [g D-glucose/l sample solution]}}$$

$$c = (3.020 \times 342.3 \times \Delta A) / (6.3 \times 1.00 \times 0.100 \times 1000) = \mathbf{1.641 \times \Delta A \text{ [g sucrose/l sample solution]}}$$

If the sample has been diluted during preparation, multiply the result with the dilution factor F.

When analyzing samples which are weighed out for sample preparation, calculate the content from the amount weighed:

$$\text{Content}_{\text{sucrose/D-glucose}} = \frac{C_{\text{sucrose/D-glucose}} \text{ [g/l sample solution]}}{\text{weight}_{\text{sample}} \text{ [in g/l sample solution]}} \times 100 \text{ [g/100 g]}$$

**Sample preparation**

If the sample has one of the characteristics below, which hamper the test, please follow the corresponding sample preparation procedure.

1. Dilute *clear, colourless and almost neutral liquid samples* to get a sample solution with 0.08 to 0.8 g sucrose + D-glucose/l.
2. Filter or centrifuge *turbid solutions*, dilute (see pt. 1).
3. Degas *samples containing carbon dioxide*, e.g. by filtration, or add NaHCO<sub>3</sub> until the solution is slightly alkaline, dilute (see pt. 1).
4. Crush (corn size < 0.3 mm) or homogenize *solid or semi-solid (pasty) samples*, extract with water or dissolve in water, filter and dilute (see pt. 1) if necessary.
5. Extract *fat containing samples* with hot water at a temperature above the melting point of fat, e.g. in a 100 ml volumetric flask. Adjust to +20°C, fill volumetric flask to the mark. Store in ice or in refrigerator for approx. 15 resp. 30 min, filter. Alternatively, clarify with Carrez reagents (which can be recommended).
6. Clarify *samples containing protein* with Carrez reagents  
Weigh sufficient quantity of solid or pasty sample into 100 ml volumetric flask, add approx. 60 ml water. Or pipette liquid sample into 100 ml volumetric flask containing approx. 60 ml water. Add, and mix after each addition, 5 ml Carrez-I-solution (3.60 g K<sub>4</sub>[Fe(CN)<sub>6</sub>] x 3H<sub>2</sub>O = potassium hexacyanoferrate(II)/100 ml), 5 ml Carrez-II-solution (7.20 g ZnSO<sub>4</sub> x 7 H<sub>2</sub>O = zinc sulfate hepta-hydrate/100 ml). Adjust to pH 7.5 to 8.5 by the addition of e.g. 10 ml NaOH (0.1 M). Fill the flask to the mark, mix and filter.
7. Do not deproteinize samples with perchloric acid because sucrose is hydrolyzed.

**Assay characteristics**

1. **Specificity:** Specific for D-glucose. Relatively specific for sucrose in the absence of 2-β-fructosans. (2-β-fructosans, if present in the sample, react slower than sucrose.) In the analysis of commercial sucrose results of 100 % have to be expected, in the analysis of D-glucose and D-glucose monohydrate results of < 100 % because the materials absorb moisture.
2. **Sensitivity:** 0.2 mg D-glucose/l (Δ A = 0.005; v = 2.000 ml; V = 3.020 ml)  
1 mg sucrose/l (Δ A = 0.010; v = 1.800 ml; V = 3.020 ml)
3. **Detection limit:** 0.4 mg D-glucose/l (Δ A = 0.010; v = 2.000 ml; V = 3.020 ml)  
2 mg sucrose/l (Δ A = 0.020; v = 1.800 ml; V = 3.020 ml)
4. **Linearity:** from 4 µg sucrose + D-glucose /assay (v = 1.800 ml; V = 3.020 ml)  
to 150 µg sucrose + D-glucose /assay (v = 0.100 ml; V = 3.020 ml)
5. **Precision:** Δ A = +/- 0.005 absorbance units (D-glucose)  
Δ A = +/- 0.010 absorbance units (sucrose)  
CV = approx. 1 to 2 % (D-glucose)  
CV = approx. 1 to 3% (sucrose)  
Fruit juice:  
r = 1.9 + 0.033 x C<sub>sucrose</sub> [g/l]  
R = 3.3 + 0.061 x C<sub>sucrose</sub> [g/l]  
r = 0.42 + 0.027 x C<sub>D-glucose</sub> [g/l]  
R = 1.0 + 0.042 x C<sub>D-glucose</sub> [g/l]
6. **Interferences:** none known.
7. **Technical Information:** the reagents can also be used for the determination of D-fructose (with additional PGI).