**CONGEN**

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## 4 General information

### 4.1 Description

The multiplex test detects DNA of peanut, hazelnut and walnut according to directive (EC) 1169/2011. Each reaction contains an internal amplification control. If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). Also spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be used with established real-time PCR instruments, equipped for the detection of four fluorescence emissions at 510 nm, 580 nm, 610 nm and 660 nm (FAM, VIC/HEX, ROX and Cy5) at the same time. Internal validation was performed on Roche LightCycler® 480 II<sup>2</sup>, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96 and Agilent Mx3005P.

#### Note:

The SureFood®ALLERGEN Peanut/Hazelnut/Walnut+IAC PCR using SureFood® PREP Advanced, protocol 1 has a limit of detection of:

- ≤ 1 mg / kg Peanut
- ≤ 0.4 mg / kg Hazelnut
- ≤ 0.4 mg / kg Walnut

The limit of detection of the complete quantitative determination (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

### 4.2 DNA-preparation

For DNA-preparation the use of SureFood® PREP Advanced, protocol 1 is recommended.

### 4.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple use on the same day.**

### 4.4 Additionally required equipment and materials

- real-time PCR instrument, equipped with four detection channels (510 nm, 580 nm, 610 nm, 660 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer

<sup>2</sup> Note: For users of the Roche LightCycler® 480 II instrument a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

# SureFood® ALLERGEN 4plex Peanut/Hazelnut/Walnut+IAC (100 React.)

Art. No. S3402

May 2018

- micro centrifuge with a rotor for the reaction tubes

## 4.5 Setup

	Blockcycler / Bio Molecular Systems MIC	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	35	35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

## 4.6 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Agilent Mx3005P</b>	Peanut	FAM	+	
	IAC	HEX	+	
	Hazelnut	ROX	+	
	Walnut	Cy5	+	
<b>Agilent AriaDx</b>	Peanut	FAM	+	
	IAC	HEX	+	
	Hazelnut	ROX	+	
	Walnut	Cy5	+	
<b>Applied Biosystems 7500</b>	Peanut	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
	Hazelnut	ROX	None	
	Walnut	Cy5	None	
<b>Bio-Rad CFX96</b>	Peanut	FAM	+	
	IAC	VIC	+	
	Hazelnut	ROX	+	
	Walnut	Cy5	+	
<b>BioMolecular Systems MIC</b>	Peanut	green	+	
	IAC	yellow	+	
	Hazelnut	ROX	+	

**SureFood® ALLERGEN 4plex  
Peanut/Hazelnut/Walnut+IAC (100 React.)**

**Art. No. S3402**

May 2018

	Walnut	Cy5	+	
--	--------	-----	---	--

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Qiagen Rotor-Gene Q</b>	Peanut	green	+	
	IAC	yellow	+	
	Hazelnut	orange	+	
	Walnut	red	+	
<b>Roche LightCycler® 480 II</b>	Peanut	465-510	+	
	IAC	533-580	+	
	Hazelnut	533-610	+	
	Walnut	618-660	+	
<b>Roche cobas® z 480 Analyzer</b>	Peanut	465-510	+	
	IAC	540-580	+	
	Hazelnut	540-610	+	
	Walnut	610-670	+	

## 5 Qualitative Analysis

### 5.1 Protocol

#### 5.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative peanut/hazelnut/walnut detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1 positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19,3 µl	212,3 µl
Taq Polymerase	0,7 µl	7,7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

**Note: The Taq Polymerase can be stored at +2 to +8 °C for multiple use on the same day.**

### 5.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

### 5.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

A sample is stated **positive** for the respective parameter (peanut/hazelnut/walnut), if the sample DNA shows amplification in the respective channel (FAM/ROX/Cy5).

A sample is stated **negative** for the respective parameter (peanut/hazelnut/walnut), if the sample DNA shows no amplification in the respective channel (FAM/ROX/Cy5). The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value  $\leq 2$  compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific soya PCR assay.

#### In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.

result in the respective channel				
FAM-channel peanut	ROX-channel walnut	Cy5-channel hazelnut	VIC/HEX-channel amplification control	result
<b>positive</b>	negative	negative	<b>positive/ negative</b>	peanut DNA detected
negative	<b>positive</b>	negative	<b>positive/ negative</b>	walnut DNA detected
negative	negative	<b>positive</b>	<b>positive/ negative</b>	hazelnut DNA detected
negative	negative	negative	negative	invalid

## 6 Further information

### 6.1 Product Information

- **Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices**  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report
- Product Information

### 6.2 Technical Support

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

### 6.3 Distribution and ordering

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)