RIDASCREEN® T-2 Toxin

Art. No. R3801

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RIDASCREEN® T-2 Toxin

Brief information

RIDASCREEN® T-2 Toxin (Art. No. R3801) is a competitive enzyme immunoassay for the quantitative determination of T-2 toxin in cereals and feed.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction, centrifugation and dilution

Time requirement: sample preparation
(for 10 samples) approx. 30 min
sample preparation with column
(for 10 samples) approx. 1 h
test implementation
(incubation time) approx. 1 h 30 min

Detection limit: Measuring range: 3.5 – 56 µg/kg (ppb)
(corresponding to standard substance) Barley, rye, corn (maize), wheat approx. 7 µg/kg (ppb)
oats approx. 11 µg/kg (ppb)

Measuring range: 35 – 560 µg/kg (ppb)
Corn (maize), wheat, oats approx. 30 µg/kg (ppb)

Recovery rate: in spiked cereal samples approx. 90 % ± 10 %

Specificity: T-2 toxin 100 %
Acetyl T-2 toxin approx. 114 %
HT-2 toxin approx. 7 %
Iso T-2 toxin approx. 2 %

The specificity of the RIDASCREEN® T-2 Toxin test was established by analyzing the cross-reactivity to corresponding mycotoxins in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the

respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related T-2 and T-2/HT-2 Toxin test systems

RIDA[®] QUICK T-2/HT-2 RQS (R5304)

RIDASCREEN[®] FAST T-2 Toxin (R3801)

RIDASCREEN[®] T-2/HT-2 Toxin (R3805)

TRILOGY[®] Liquid Standard T-2 Toxin (TSL-314)

TRILOGY[®] Dried Standard T-2 Toxin (TS-314)

TRILOGY[®] Liquid Standard HT-2 Toxin (TSL-333)

TRILOGY[®] Dried Standard HT-2 Toxin (TS-333)

1. Intended use

RIDASCREEN[®] T-2 Toxin is a competitive enzyme immunoassay for the quantitative determination of T-2 toxin in cereals and feed.

2. General

T-2 toxin belongs to the trichothecene group of mycotoxins and is formed by fungi of the genus *Fusarium*. T-2 toxin is often found in agricultural commodities, although the incidence and the concentrations found show a broad regional variation. Due to its cytotoxic and immunosuppressive mode of action T-2 toxin is a threat for human and animal health.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-T-2 toxin antibodies. T-2 toxin standards or sample solutions, T-2 toxin enzyme conjugate and anti-T-2 toxin antibodies are added. Free T-2 toxin and T-2 toxin enzyme conjugate compete for the T-2 toxin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-T-2 toxin antibodies are also bound by the immobilized

capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Enzyme substrate and chromogen are added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the T-2 toxin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
Microtiter plate M		Ready to use	96 wells
Sample buffer	White	Ready to use	60 ml
Standard 1	White	Ready to use	0 µg/L 1.3 ml
Standard 2	White	Ready to use	0.1 µg/L 1.3 ml
Standard 3	White	Ready to use	0.2 µg/L 1.3 ml
Standard 4	White	Ready to use	0.4 µg/L 1.3 ml
Standard 5	White	Ready to use	0.8 µg/L 1.3 ml
Standard 6	White	Ready to use	1.6 µg/L 1.3 ml
Conjugate	Red	Ready to use	6 ml
Antibody	Black	Ready to use	6 ml
Substrate	Green	Ready to use	7 ml
Chromogen	Blue	Ready to use	7 ml
Stop solution	Yellow	Ready to use	14 ml

5. Materials required but not provided

5.1. Equipment

- grinder (mill)
- magnetic stirrer
- blender
- centrifuge or paper filter
- evaporator
- glas tubes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- microtiter plate spectrophotometer (450 nm)

5.2. Reagents

- Methanol
- Acetonitrile
- PuriTox Trichothecene; Product code: TC-T220
- for dilution of samples (for dilutions > dilution factor of 35):
PBS buffer, pH 7.2: (0.55 g $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$ + 2.85 g $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$ + 9 g NaCl; fill up to 1000 ml with methanol/distilled water 10/90 (v/v)

This buffer has to contain 10 % methanol in order to keep a 10 % methanol concentration of the samples used in the assay (also see section 9.1./9.2.).

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

The standards contain T-2 toxin, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and T-2 toxin solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

T-2 toxin is light sensitive, therefore, avoid exposure to direct light.

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the chromogen prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

9.1. Cereals: Measuring range 3.5 – 56 ppb (sample preparation with solid phase column: PuriTox Trichothecene; TC-T220)

- weigh 5 g of ground sample and add 20 ml acetonitrile/dist. water 84/16 (v/v)
- blend at high speed for 3 min
- centrifuge: 10 min / ≥ 3000 rpm / room temperature (20 - 25 °C / 68 – 77 °F)
- pass 6 ml of the supernatant through the PuriTox Trichothecene (TC-T220) column by applying pressure with a plunger and collect the cleaned-up filtrate in a glass tube; pass air through the column to remove residual liquid
- evaporate 4 ml of the cleaned-up filtrate to dryness under air at 60 - 70 °C
- dissolve the dried residue in 5 ml methanol/dist. water 70/30 (v/v)
- dilute the sample 1:7 (1+6) with sample buffer (see 4.), (e.g. 100 μ l sample + 600 μ l sample buffer)
- use 50 μ l per well in the assay

Note: The clean-up method was originally developed for GC analysis by R-Biopharm Rhone (see instruction for use TC-T220). The method described above for the ELISA analysis was modified by R-Biopharm.

Remark:

At high T-2 toxin concentrations (> 56 ppb), the 1:7 diluted extract solution must be further diluted e. g. 1:10 (1+9) with PBS buffer containing 10 % methanol (see 5.2.), e. g. 50 μ l of the diluted extract solution + 450 μ l PBS buffer with 10 % methanol.

9.2. Cereals and feed: Measuring range 35 – 560 ppb

- weigh 5 g of sample and dissolve with 25 ml of methanol/distilled water 70/30 (v/v) for extraction*)
- blend at high speed for 3 min
- filter the extract over paper filter or centrifuge the extract (10 min / ≥ 3000 rpm / room temperature (20 - 25 °C / 68 – 77 °F))
- dilute the filtrate or the supernatant 1:7 (1+6) with sample buffer (see 4.), e. g. 100 μ l of filtrate or supernatant + 600 μ l of sample buffer
- the diluted extract solution must be further diluted; e.g. 50 μ l diluted extract + 450 μ l PBS buffer containing 10 % methanol (see 5.2.)
- use 50 μ l per well in the assay

*) sample size may be increased if required, but the volume of methanol/distilled water 70/30 (v/v) must be adapted accordingly, e.g.: 25 g in 125 ml methanol/distilled water 70/30 (v/v) or 50 g in 250 ml methanol/distilled water 70/30 (v/v)

Remark:

At high T-2 toxin concentrations (> 560 ppb), the 1:7 diluted extract solution must be further diluted e. g. 1:10 (1+9) with PBS buffer containing 10 % methanol (see 5.2.), e. g. 50 μ l of the diluted extract solution + 450 μ l PBS buffer with 10 % methanol.

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.

2. Add 50 µl of standard or prepared sample to separate duplicate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of conjugate to each well.
4. Add 50 µl of the antibody to each well. Mix gently by shaking the plate manually and incubate for 1 h at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl distilled water and pour out the liquid again. Repeat the washing procedure two times.
6. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

11. Results

A special software, the RIDASOFT® Win.NET (Art. No. Z9996), is available for evaluation of the RIDASCREEN® enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the T-2 toxin concentration [µg/kg].

In order to obtain the T-2 toxin concentration in µg/kg actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

Cereals (9.1):	35 (or 350)
Ceraels, feed (9.2):	350 (or 3500)

Therefore, the standard curve is in the range of 3.5 to 56 µg/kg (ppb) T-2 toxin when using sample preparation 9.1. or in the range of 35 to 560 µg/kg (ppb) T-2 toxin when using sample preparation 9.2..

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-mail: info@r-biopharm.de

www.r-biopharm.com

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are registered trademarks of R-Biopharm AG