r-biopharm®

RIDASCREEN® FAST Peanut

Art. No. R6202

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RIDASCREEN®FAST Peanut

Brief information

The RIDASCREEN®FAST Peanut (Art. No.: R6202) is a sandwich enzyme immunoassay developed for the quantitative analysis of peanut and parts of peanut in food. Exemplarily for the food groups pastries and sweets the following commodities have been tested during the assay's development: breakfast cereals, cookies, ice cream, and milk chocolate. The ELISA was approved by the AOAC *Performance Tested Methods*SM Program and was assigned PTM Certification No. 030404 by the AOAC Research Institute. The assay can be likely used for analysis of other food samples too; this has to be checked by the user. Limit of detection and limit of quantification depend on the tested food matrix, the level of food processing, and the extraction method. Detailed results and further validation information about other sample matrices can be found in the brochure *Product Information*. Additionally, results from proficiency tests and inter laboratory trials may exist for different food matrices. Further applications are frequently validated in our laboratory and provided in the *Application Notes*.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction, and centrifugation

Time requirement: sample preparation (for 10 samples) approx. 20 min
test implementation (incubation time) 30 min

Limit of detection: 0.13 mg/kg (ppm) peanut
0.03 mg/kg (ppm) peanut protein

Limit of quantification: 2.5 mg/kg (ppm) peanut
0.6 mg/kg (ppm) peanut protein

Standard material: reference material NIST SRM 2387 peanut butter

Specificity: The antibodies specifically detect peanut proteins, including the peanut allergen Ara h 1 and Ara h 2.

There is a cross reactivity to green pea, lentils, wheat semolina and fenugreek.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments. During the development of the test kit more than 70 potential cross reactants are tested. Further information is described in the updated validation report.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under from the website <http://www.r-biopharm.com/products/food-feed-analysis>.

Related products

- RIDASOFT® Win.NET (Z9996)
- Bioavid Lateral Flow Erdnuss/Peanut (BL606-10 / BL606-25)
- SureFood® ALLERGEN Peanut (S3603)
- SureFood® ALLERGEN 4 Plex Peanut/Hazelnut/Walnut + IAC (S3402)

1. Intended use

The RIDASCREEN®FAST Peanut (Art. No.: R6202) is a sandwich enzyme immunoassay developed for the quantitative analysis of peanut and parts of peanut in food. Exemplarily for different food groups, breakfast cereals, cookies, ice cream, and milk chocolate have been approved. The assay can be likely used for analysis of other food samples too; this has to be checked by the user. Limit of detection and limit of quantification depend on the tested food matrix, the level of food processing, and the extraction method. Detailed results and further validation information about other sample matrices can be found in the brochure *Product Information*. Additionally, results from proficiency tests and inter laboratory trials may exist for different food matrices. Further applications are frequently validated in our laboratory and provided in the *Application Notes*.

2. General

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, peanut and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

RIDASCREEN®FAST Peanut was validated with breakfast cereals, cookies, ice cream, and milk chocolate in an independent study by three independent US laboratories under the supervision of the AOAC Research Institute. Each laboratory analyzed 20 blank samples and 20 samples spiked at 5 mg/kg peanut of each matrix. In total, two of 240 blank samples were measured positive (> 2.5 mg/kg peanut) corresponding to a specificity of 99.2 % (238/240). All of the spiked samples (240/240) were found positive (sensitivity = 100 %).

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to peanut. By adding standards and samples to the wells, peanut present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of peanut takes place by adding substrate/chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color

change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the peanut concentration of the sample. The result is expressed in mg/kg peanut.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0 mg/kg	2.6 ml
Standard 2*	Transparent	Ready to use	2.5 mg/kg	2.6 ml
Standard 3*	Transparent	Ready to use	5.0 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	10.0 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	20.0 mg/kg	1.3 ml
Wash buffer**	Brown	Concentrate	10x	100 ml
Conjugate	Red	Ready to use		6 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

- *) The dilution factor 20 for the sample has already been considered when labeling. Therefore, the peanut concentration of a sample can directly be read from the standard curve.
- ***) The peanut washing buffer has to be used in the RIDASCREEN®FAST Peanut only.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- shaker
- water bath
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax, or homogenisator
- graduated pipettes
- variable 20 - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents:

- distilled or deionized water
- skim milk powder

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any of the test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained substrate/chromogen prior to test implementation
- a value of less than 1.2 absorbance units ($A_{450\text{ nm}}$) for standard 5

9. Preparation of Samples

Working devices such as a mill, glass vials, or spatulas must be cleaned before and after each sample preparation to remove any remains of peanut and to avoid contamination.

The **Allergen Extraction buffer** is provided as a **10fold concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate **1:10 (1+9) with distilled water** before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is either stable at 20 - 25 °C (68 - 77 °F) for four weeks or at 2 - 8 °C (35 - 46 °F) for 12 weeks.

9.1. Sample preparation for all samples not mentioned in 9.2.

- grind 5 g of the sample carefully and mix thoroughly
- weigh 1 g of sample and add 20 ml diluted extraction buffer (the extraction buffer should have already been heated to approx. 60 °C (140 °F))
- mix intensively and extract for 10 min at 60 °C (140 °F) by shaking casually
- centrifuge: 10 min, 2500 x *g*, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged for 10 minutes at high speed in reaction caps by using a micro centrifuge)
- use 100 µl of the supernatant or filtrate per well in the assay

9.2. Sample preparation **with** skim milk powder for tannin or polyphenol containing samples (e.g. spices, pepper, paprika, chocolate, caramel, and ice cream)

- 5 g of the sample should be ground well and thoroughly mixed
- Take hereof 1 g sample and add 1 g skim milk powder; mix properly
- add 20 ml diluted extraction buffer (the extraction buffer should have already been heated to approx. 60 °C (140 °F))
- mix intensively and extract for 10 min at 60 °C (140 °F) by shaking casually
- centrifuge: 10 min, 2500 x *g*, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged for 10 minutes at high speed in reaction caps by using a micro centrifuge)

- use 100 µl of the supernatant or filtrate per well in the assay

Remark:

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 5 days. The extracts can be stored at -20 °C (-4 °F) for several months.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at room temperature for approx. four weeks. The peanut washing buffer has to be used in the RIDASCREEN®FAST Peanut only.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.

4. Add 100 µl of the enzyme conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of the reddish substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win.NET (Art. No. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. The document *Compliance Criteria* provides criteria for evaluating standard curves.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the standard curve, especially for the zero standard, may be a result of insufficient washing or peanut contamination.

Please note:

When working according to the test kit instructions, the dilution factor is 20. The allergen concentration can be directly read from the standard curve (see 4. *) - the sample dilution factor of 20 is already taken into account).

For sample dilutions of more than 1:20, the further dilution factor must be considered for the calculation of the peanut concentration.

The test is calibrated against NIST Reference Material 2387 (peanut butter). The result is expressed as mg/kg peanut. The reference material 2387 contains 22.2 % +/- 1% total protein. If the result is multiplied by 0.22, then the result can be expressed as mg/kg (ppm) total protein.

In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

Allergen containing samples that have been heat treated show a reduced recovery because the proteins denature and are no longer recognized by the antibody. The reduction in recovery depends strongly on the temperature and the duration of the heat treatment. If samples are heat treated at high temperature the recovery can be significantly reduced.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report.

The protein content and the protein composition may vary considerably between different peanut species. Therefore, different varieties may produce different results, since exemplary varieties were used for calibration.

Recommendation:

In order to ensure a high analytical performance we recommend:

- Adjust the pH to a neutral value for extremely acidic or alkaline samples
- Use allergen-free and allergen containing (spiked) samples as test controls
- Carry out spiking experiments for an accurate and correct procedure
- Analyze each sample material in duplicates
- Perform SureFood[®] PCR to confirm results
- Contact sales@r-biopharm.de if automates (e.g. ThunderBolt[®] / Bolt[®]) are used

Further application notes:

- More sensitive analysis with RIDASCREEN®FAST Peanut (dilution of standards down to 0.625 mg/kg peanut)

For further product information and application notes, please contact sales@r-biopharm.de.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG