CONGEN

**SureFood® ANIMAL ID 4plex
Beef/Sheep/Goat+IAAC**

Art. No. S6121

100 rxn

User Manual

November 2018

 Inhalt /  Content

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1 General Information

1.1 Description

The test detects beef (*Bos taurus*), sheep (*Ovis aries*) and goat (*Capra hircus*) DNA. Each reaction contains an internal amplification control and an internal detection assay for vertebrates DNA (IAAC). The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of four fluorescence emissions at 522 nm, 553 nm, 610 nm and 670 nm (FAM, VIC, ROX and Cy5) at the same time. The technical validation of instruments was performed on Roche LightCycler® 480 II², Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96 and Agilent Mx3005P.

1.2 Limit of Detection

The SureFood® ANIMAL ID 4plex Beef/Sheep/Goat+IAAC real-time PCR is developed for the detection of beef, sheep and goat DNA in muscle meat mixture at a relative amount of 0.1 %. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

Note: In mixed samples, inconsistent mixing ratios (e.g. 99.9 % beef and 0.1 % sheep) may cause a loss of sensitivity in the low concentration channel.

1.3 DNA-preparation

For DNA-preparation the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -20°C and protected from light.

1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic or SureFood® PREP Advanced)
- real-time PCR instrument with four detection channels (522 nm, 553 nm, 610 nm and 670 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps) pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

² note: For the use of the Roche LightCycler® 480 II a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

1.6 Setup

	Blockcyler/MIC/LightCycler® 480 II	Rotorgene Q
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	35	35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of Extension Phase Detection System beef: Various devices FAM -channel, Quencher: BHQ LightCycler 480 II 465 nm - 510 nm Detection System sheep: Various devices VIC/HEX -channel, Quencher: BHQ LightCycler 480 II 533 nm - 580 nm Detection System goat: Various devices ROX -channel, Quencher: BHQ LightCycler 480 II 533 nm - 610 nm Detection System vertebrates and internal amplification control (IAAC): Various devices Cy5 -channel, Quencher: BHQ LightCycler 480 II 618 nm - 660 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The results of the control reactions have to be correct.

Beef DNA is detected in the FAM-channel, sheep DNA is detected in the VIC-channel and goat DNA is detected in the ROX-channel. In the Cy5-channel it is possible to detect animal DNA in the sample as well as the amplification control (IAC) in a sample with no animal DNA inside.

A sample is stated **positive** for the respective parameter (beef/sheep/goat), if the sample DNA shows amplification in the respective channel. A sample is stated **negative** for the respective parameter (beef/sheep/goat), if the sample DNA shows no amplification in the respective channel and the internal signal (Cy5-channel) of the sample is **positive** (see also table below).

result in the respective channel				result
FAM-channel beef	VIC/HEX-channel sheep	ROX-channel goat	Cy5-channel ¹ animal + IAC	
positive	negative	negative	positive	beef DNA detected
negative	positive	negative	positive	sheep DNA detected
negative	negative	positive	positive	goat DNA detected
negative	negative	negative	positive	vertebrates DNA² detected
negative	negative	negative	negative	Invalid ³

¹ If the signal of the internal control (Cy5-channel) of the sample DNA is detected significantly before the signal of the negative control (master-mix without DNA) the sample contains animal DNA.

² Is the Cp-value of the internal control (Cy5-channel) in the range of the negative control the sample contains no PCR-inhibiting substances but only a low amount or no animal DNA.

³ In the case of a negative result in all channels including the Cy5-channel the sample contains PCR inhibiting substances. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

Note: If the sample contains more than one animal species the DNA mixture can have a competitive influence on the absolute fluorescence. The lower the concentration of the determinant DNA is in a mixture of animal DNAs the lower is the fluorescence level of the amplification curve.

November 2018

3 Further Information

3.1 Product Information

- Validation Report

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com