Premi[®]Test Urin

Art. No. R3921

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Premi[®]Test Urine

Brief information

Premi[®]Test Urine (Art. No. R3921) is used for sample preparation of bovine and porcine urine for the analysis with Premi[®]Test. Premi[®]Test is a microbial inhibition test for screening of antibiotic residues in food and products of animal origin. Premi[®]Test (Art. No. R3925 / R3900) is required for analysis. The number of samples which can be prepared with Premi[®]Test urine is dependent on the volume of the samples.

Sample preparation	Dilution with Premi [®] Test Urine, mixing
Time requirement:	Urine (10 samples) approx. 15 min test implementation (incubation time) 3 – 4.5 h
Limit of Detection:	
β -Lactams	Amoxicillin 5 μ g/kg (ppb) Ampicillin 5 μ g/kg (ppb) Penicillin G 5 μ g/kg (ppb) Cloxacillin 50 μ g/kg (ppb)
Cephalosporins	Ceftiofur 100 μ g/kg (ppb) Cephapirin 100 μ g/kg (ppb)
Macrolide	Tylosin 50 μ g/kg (ppb) Erythromycin 100 μ g/kg (ppb) Lincomycin 100 μ g/kg (ppb) Tilmicosin 50 μ g/kg (ppb) Tulathromycin 18,000 μ g/kg (ppb)
Tetracycline	Chlortetracycline 50 μ g/kg (ppb) Oxytetracycline 50 μ g/kg (ppb)
Sulfonamide	Sulfamethazine 100 μ g/kg (ppb) Sulfadimethoxin 100 μ g/kg (ppb) Sulfachloropyridazine 100 μ g/kg (ppb)
Aminoglykoside	Gentamicin 100 μ g/kg (ppb) Neomycin 100 μ g/kg (ppb) Dihydrostreptomycin 3,000 μ g/kg (ppb)
Quinolones	Enrofloxacin 600 μ g/kg (ppb) Danofloxacin 600 μ g/kg (ppb)
Others	Florfenicol 100 μ g/kg (ppb)

Related products

Premi[®]Test (R3925/R3900)

Premi[®]Test Starter Kit (ZPT-2000)

Premi[®]Test Multipress (ZPT-2012)

1. Intended use

Premi[®]Test and Premi[®]Test urine are used in combination for the screening of antibiotic residues in bovine and porcine urine.

2. General

Antibiotics are frequently used for the treatment of bacterial infections and diseases in animal breeding. Improper use or failure to comply with the statutory waiting period may result in antibiotic residues in food of animal origin. In addition, antibiotics are illegally abused as antimicrobial growth promoters. Because of the potential carcinogenic, toxic and allergic properties of antibiotic residues, contaminated food may pose a direct health risk to consumers. In addition, the improper or excessive use of antibiotics in food production is suspected to promote the formation of multi-resistant pathogens. For reasons of consumer protection, the legislative authorities have established Maximum Residue Limits (MRLs) for antibiotic residues in food, e.g. by EU Regulation 37/2010.

3. Test principle

Premi[®]Test is a microbial inhibition test in ampoule form which uses the natural growth inhibitory function of antibiotic residues on the thermophilic test organism *Bacillus stearothermophilus*. Each Premi[®]Test ampoule contains a standardized number of spores which are embedded together with nutrients and the redox indicator bromocresol in an agar matrix. In the absence of antibiotics, the growth of the microorganism leads to an acid production and to a color change of the pH indicator from its initial purple to yellow. When samples which are contaminated with antibiotic residues are tested, the color of the ampoule remains purple. The comparison of the samples with a negative control which verifiable contains no antibiotic residues, allows a qualitative evaluation of the sample (negative/positive).

4. Reagents provided

Depending on the sample volume, a different number of samples may be prepared with the reagents, e.g. 30 samples, at a volume of 30 ml.

Component	Cap color	Format	Volume
Sample buffer	Blue	Ready to use	15 ml

5. Materials required but not provided

5.1. Reagents:

- Premi[®]Test (R3925/R3900)
- Negative control (urine of the same animal species), which verifiable contains no antibiotic residues

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain further hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 3 - 10 °C. Do not freeze any test kit components.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Turbidity of the reagent may be an indication for a microbial contamination.

9. Preparation of Samples

During sampling of urine, a contamination of the urine samples with faeces has to be avoided because urease and *enterobacteria* influence the test system.

Urine has to be tested within the day of sampling. Alternatively, samples can be frozen at -20 °C and stored for 3 months.

- Prewarm Premi[®]Test Urine to room temperature before use (20 – 25 °C)
- Transfer 10, 20, 30, 40 or 50 ml of bovine or porcine urine into a sample vial
- In accordance to the table, add the appropriate number of drops of Premi[®]Test urine to the urine sample:

urine volume (ml)	Number of drops of Premi [®] Test Urine
10	5
20	10
30	15
40	20
50	25

- Mix urine sample and Premi[®]Test Urine well
- Allow the samples to stand for 10 min until eventual sediment has settled and a clear supernatant can be seen. Optional samples can be centrifuged (5 min; 3000g; room temperature (20 – 25 °C)).

10. Test implementation

10.1. Test preparation

Wash hands and clean the working surface before starting the procedure.

Prewarm Premi[®]Test to room temperature (20 – 25 °C) before use.

10.2. Test procedure

1. Cut off the required number of Premi[®] Test ampoules and remove the aluminum foil from the ampoules.
2. Transfer 100 µl of the clear supernatant of the prepared urine sample onto the agar in the ampoule. Use a new clean pipette tip on the syringe for every sample and avoid to touch or damage the agar during pipetting.
3. Transfer 100 µl of the clear supernatant of the prepared negative control, which verifiable contains no antibiotic residues, into one ampoule.
4. Allow to stand at room temperature for 20 minutes.
5. Flush the urine away by gently filling and emptying the ampoules twice with distilled water.
6. Remove residual liquid by placing the ampoule upside-down on an absorbent paper towel. Do not remove residual water by tapping the ampoules.
7. Seal the ampoules with the perforated plastic foil supplied.
8. Incubate the ampoules at $64 \pm 0.5^{\circ}\text{C}$ in the Premi[®]Test incubator (supplied with the Starter Kit (Art. No. ZPT-2000)) or in a tempered water bath.

11. Results

- When the negative control changes its color from purple to a distinctive yellow, the results of the samples can be interpreted. The incubation time for bovine urine is approx. 3 h 40 min to 4 h 30 min and approx. 2 h 40 min to 3 h 30 min for porcine urine, but may vary slightly between runs. Therefore, it is recommended to check the color at least 10 min earlier and then intervals of 5 min until the color of the negative control changes. Read the results from the bottom 2/3 part of the agar only.
- If the color of the sample remains purple or if the color of the sample is clearly different to the color of the negative control, the sample contains antibiotic residues at a concentration above the Premi[®]Test detection limits. A positive result can also be caused by residues of 2 or more antibiotics where the respective individual concentrations are below the limit of detection or the MRL, but where the summarized effect inhibits the growth of the testing microorganism.
- A clear color change from purple to yellow indicates that the sample contains no antibiotic residues or that the concentration of antibiotic residues is below the Premi[®]Test detection limits.

More applications for Premi[®]Test for the screening of antibiotic residues are:

- Meat
- Liver
- Kidney
- Fish
- Shrimp

For the product information folder with further information and applications please contact: info@r-biopharm.de.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG
R-Biopharm AG