# **RIDASCREEN<sup>®</sup> Testosteron**

**Art. No. R2401**

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# RIDASCREEN<sup>®</sup> Testosteron

## Brief informatin

RIDASCREEN<sup>®</sup> Testosteron (Art. No.: R2401) is a competitive enzyme immunoassay for the quantitative analysis of testosterone in bovine plasma.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction, freezing and evaporation

Time requirement: sample preparation (for 10 samples) ..... approx. 2 h  
test implementation (incubation time) .....2.5 h

Detection limit: 20 ng/l (ppt)

(corresponding to the standard substance)

Recovery rate: in spiked plasma samples .....85 % ± 12 %

(corresponding to the standard substance)

The specificity of the RIDASCREEN<sup>®</sup> Testosteron-test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:	Testosterone .....	100 %
	17 $\beta$ -Nortestosterone .....	10 %
	17 $\beta$ -Trenbolone.....	1 %
	17 $\beta$ -Estradiol.....	< 0.1 %
	Cortisol .....	< 0.1 %
	Corticosterone.....	< 0.1 %
	Progesterone.....	< 0.1 %

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## Related products

RIDA <sup>®</sup> Testosteron Spiking Soltuion.....	(R2499)
RIDASCREEN <sup>®</sup> 17 $\beta$ -Östradiol .....	(R2301)
RIDA <sup>®</sup> 17 $\beta$ -Östradiol Spiking Solution .....	(R2399)

### 1. Intended use

RIDASCREEN<sup>®</sup> Testosteron is a competitive enzyme immunoassay for the quantitative analysis of testosterone in bovine plasma.

### 2. General

Testosterone is a natural male steroid hormone. Because of the illegal use as resource in cattle fattening, in the German residue control plan the legally demanded detection limit in bovine plasma for HPLC-SP and GC-MS is < 0.1  $\mu\text{g}/\text{kg}$  (ppb).

### 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-testosterone antibodies. Standards or sample, testosterone conjugate and anti-testosterone antibodies are added. Free and enzyme conjugated testosterone compete for the antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-testosterone antibodies are also bound by the immobilized capture antibodies. Any unbound conjugate is then removed in a washing step. Substrate and chromogen are added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the testosterone concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses).

Componet	Cap Colour	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Standard 1	White	Ready to use	0 ng/l	1.3 ml
Standard 2	White	Ready to use	50 ng/l	1.3 ml
Standard 3	White	Ready to use	200 ng/l	1.3 ml
Standard 4	White	Ready to use	800 ng/l	1.3 ml
Standard 5	White	Ready to use	3200 ng/l	1.3 ml
Standard 6	White	Ready to use	12800 ng/l	1.3 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Antibody	Black	Concentrate	11x	0.7 ml
Conjugate/antibody/sample buffer	White	Ready to use		50 ml
Substrate	Green	Ready to use		7 ml
Chromogen	Blue	Ready to use		7 ml
Stop solution	Yellow	Ready to use		14 ml

## 5. Materials required but not provided

### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- freezer (-25 °C or -60 °C / -13 °F or -76 °F)
- water bath (60 °C / 140 °F)
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

### 5.2. Reagents:

- tert. butylmethylether p. a.
- petroleum ether

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- any coloration of the colorless chromogen solution prior to test implementation
- a value of less than 0.6 absorbance units ( $A_{450\text{ nm}} < 0.6$ ) for standard 1

## 9. Preparation of Samples

### 9.1. Bovine plasma

- extract 1 ml bovine plasma with 5 ml ether mixture (tert. butylmethylether / petroleum ether (30:70, v/v) in a glass tube by shaking thoroughly at room temperature (20 - 25 °C / 68 - 77 °F) for 20 min
- freeze the extracted solution at -25 °C (-13 °F) for 60 min (or at -60 °C (-76 °F) for 30 min) and decant the supernatant into a glass tube
- evaporate the ether supernatant at 60 °C (e. g. in a water bath)
- redissolve the dried residue in 400 µl sample dilution buffer and mix thoroughly
- employ 20 µl per well in the assay

### Remark:

**Plasma samples, as well as sample extracts obtained, can be stored at 2 - 8 °C (36 - 46 °F) for up to one week.**

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **testosterone conjugate** is provided as a concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in conjugate buffer (e. g. 200 µl concentrate + 2 ml conjugate buffer, ready to use sufficient for 4 microtiter strips).

The **anti-testosterone antibody** is provided as a concentrate. Since the diluted antibody has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the antibody concentrate should be shaken thoroughly. For reconstitution, the antibody concentrate is diluted 1:11 (1+10) in antibody buffer (e. g. 200 µl concentrate + 2 ml antibody buffer, ready to use sufficient for 4 microtiter strips).

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 20 µl of each standard or prepared sample to separate duplicate wells.
3. Add 50 µl of diluted conjugate to each well.
4. Add 50 µl of the diluted antibody to each well. Mix gently by shaking the plate manually and incubate 2 h at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all wells with 250 µl of distilled water and pour out the liquid again. Repeat two more times.
6. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

## 11. Results

A special software, the RIDA®SOFT Win / RIDA®SOFT Win.net (Art. No. Z9996), is available for evaluation of the RIDASCREEN® enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate, enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the testosterone concentration [ng/kg].

In order to obtain the testosterone concentration in ng/kg actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factor is as follows:

bovine plasma .....0.4

**For further information or applications please contact your local distributor or [info@r-biopharm.de](mailto:info@r-biopharm.de)**

### **Further application notes:**

–Sample preparation for meat

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG