RIDASCREEN[®] Gliadin competitive

Art. No. R7021

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Symbol	Beer			Starch syrup		Sourdough	
	gluten-free	15 mg/kg prolamin	50 mg/kg prolamin	gluten-free	naturally contaminated	35 mg/kg prolamin	75 mg/kg prolamin
Total number of laboratories	13	12	11	13	13	13	13
Total number of replicates	26	24	22	26	26	26	26
Overall mean	1.2	13.1	59.7	0.7	5.3	24.2	72.8
Recovery (%)		87	119			69	97
Repeatability SD* s(r) (mg/kg)	1.2	4.0	18.6	1.0	0.9	5.6	14.2
Reproducibility SD* s(R) (mg/kg)	1.5	4.8	18.6	1.5	1.8	6.3	20.0
Repeatability rel. SD* (%) RSD(r)	97.9	30.2	31.2	157.2	16.3	23.1	19.5
Reproducibility rel. SD* (%) RSD(R)	126.1	36.9	31.2	236.1	34.4	25.9	27.5
Horrat value	8.1	3.4	3.6	13.8	2.8	2.6	3.3

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RIDASCREEN® Gliadin competitive

Brief information

RIDASCREEN® Gliadin competitive (Art. No. R7021) is used for the analysis of fermented and hydrolyzed food (e.g. beer, starch syrup, starch, malt extract, sourdough, soy sauce) which are declared as “gluten-free”.

The competitive enzyme immunoassay quantitates peptide fragments of prolamins from wheat (gliadins), rye (secalin) and barley (hordein). The used R5 monoclonal antibody recognizes among others the potentially toxic sequence QQPFP, which occurs repeatedly in the prolamins molecules.

Der R5 ELISA RIDASCREEN® Gliadin competitive is

- is accepted as AOAC-OMA (2015-05), First Action Status
- validated by the AACCI (AACCI 38-55.01).

All reagents required for the enzyme immunoassay are contained in the test kit. The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization, extraction
Time requirement:	sample preparation (for 10 samples).....approx. 30 min test implementation (incubation time).....40 min
Limit of detection:	2.3 mg Gliadin / kg food (approx. 4.6 mg gluten / kg food) (depending on the matrix)
Limit of quantification:	5 mg gliadin / kg food (approx. 10 mg gluten / kg food)
Standard material:	The RIDASCREEN® standard material is a hydrolysate (mixture of wheat, rye and barley).
Specificity:	The R5 monoclonal antibody recognizes potentially toxic peptide sequences of gliadins from wheat and related prolamins from rye and barley.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN® Gliadin (Art. No. R7001)

RIDASCREEN® FAST Gliadin (Art. No. R7002)

RIDA® QUICK Gliadin (Art. No. R7003, R7004, R7005)

RIDA® Extraction Solution (Art. No. R7098)

Cocktail (patented) (Art. No. R7006, R7016)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

SureFood® ALLERGEN (real time PCR) Gluten (Art. No. S3106)

SureFood® ALLERGEN QUANT (real time PCR) Gluten (Art. No. S3206)

1. Intended use

RIDASCREEN® Gliadin competitive (Art. No. R7021) is used for the analysis of fermented and hydrolyzed food (e.g. beer, starch syrup, starch, malt extract, sourdough, soy sauce) which are declared as “gluten-free”.

2. General

The use of wheat flour and gluten in foodstuffs is extremely common to improve e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley.

Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (Alinorm 08/31/26) two categories for labeling of food according to the gluten content now exist:

1. Food products which contain less than 20 mg/kg (ppm) can be labeled as **"gluten-free"**.

2. Food products labeled as "**very low gluten**" can have a gluten content above 20 and up to 100 mg/kg (ppm).

During food processing like e.g. fermentation or hydrolysis intact prolamin molecules are partly or completely degraded to small peptide fragments. For Coeliac patients these fragments still remain dangerous even after digestion in the stomach.

Single, small peptide sequences (motifs) cannot be detected by a sandwich ELISA format, because at least two epitopes are necessary for a sandwich ELISA. However, single peptide fragments can be detected with a competitive format. RIDASCREEN® Gliadin competitive (Art. No. R7021) contains a new standard material. For this, wheat, rye and barley were digested by pepsin and trypsin, the peptide fragments were mixed and lyophilized after protein determination. The protein content was determined according to Dumas.

The assay can be related to the prolamin concentration and therefore to the limit values stipulated by the Codex Alimentarius. The result is expressed in mg/kg gliadin. The hydrolyzed standard was produced by the working group of Prof. Dr. Köhler (German Research Centre for Food Chemistry).

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with gliadin. Standards (hydrolysate mixture of wheat, rye and barley), respectively sample solutions and peroxidase labeled anti-gliadin antibodies (conjugate with monoclonal R5-antibodies) are added. Free and immobilized gliadin compete for the gliadin antibody binding sites (competitive enzyme immunoassay). Any unbound enzyme conjugate is then removed by a washing step. Substrate/Chromogen is added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product. Addition of the stop solution causes a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the gliadin respectively prolamin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Buffer	White	Concentrate	5x	60 ml
Standard 1*	Transparent	Ready to use	0 ng/ ml	1.3 ml
Standard 2*	Transparent	Ready to use	10 ng/ ml	1.3 ml
Standard 3*	Transparent	Ready to use	30 ng/ ml	1.3 ml
Standard 4*	Transparent	Ready to use	90 ng/ ml	1.3 ml
Standard 5*	Transparent	Ready to use	270 ng/ ml	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

*The standards contain a hydrolysate of wheat, rye and barley prolamins. Contrary to the AOAC OMA Approval 2015.05 the concentrations of the standards are given in gliadin. To convert gliadin into gluten the factor of 2 used (definition Codex Alimentarius).

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials with a screw top
- shaker
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenizer
- graduated pipettes
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents:

- distilled or deionized water
- for the extraction of samples like starch and syrup an **ethanol solution (60 %)** is required: e.g. add 150 ml ethanol p.a. to 100 ml distilled water and shake well
- for the extraction of polyphenol containing samples like **beer, malt and hops** an **ethanol solution (60 %) containing 10 % liquid fish gelatin** is required: e.g. add 30 ml of distilled water into a 100 ml graduated cylinder; add 10 g of liquid fish gelatin (Serva, Art. No. 22156 or Sigma Art. No. G-7765, both 45% solid content) and mix well; add 60 ml ethanol p.a., mix and adjust pH with a pH electrode to 8.5 and fill up to 100 ml with distilled water. The fish gelatin cannot be dissolved completely in 60% ethanol, therefore take the required volume for

the sample preparation while stirring. The fish gelatin/ethanol solution can be stored for 2 weeks at room temperature.

–1 M sodium hydroxide

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag and reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The red stained substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

–any bluish coloration of the reddish substrate/chromogen prior to test implementation

–a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the zero standard

9. Preparation of Samples

9.1. Preliminary comments

Airborne cereal dust and used laboratory equipment may lead to gliadin contamination of the assay. Therefore, wear gloves during the assay and before starting with the assay:

–clean surfaces, glass vials, mincers and other equipment with 60 % ethanol (see 5.2.)

- carry out the sample preparation in a room isolated from the ELISA procedure
- check for gliadin contamination of reagents and equipment with the test strips RIDA[®]QUICK Gliadin (Art. No. R7003, R7004, R7005)

9.2. Extraction of samples

For the RIDASCREEN[®] Gliadin competitive (R7021) only samples extracted with ethanol should be tested (the Cocktail (patented) should not be used).

The yellow stained **buffer** is provided as a concentrate (5fold). Only the amount which is actually needed should be diluted with distilled water (e. g. 3 ml concentrate + 12 ml distilled water, sufficient for the dilution of 10 samples). Make sure that the buffer is not contaminated with gliadin from the sample.

Homogenize a representative amount of food material (5 - 50 g).

- Solid food (e.g. starch and starch syrup):** weigh 1 g of a representative, homogeneous sample and add 10 ml 60 % ethanol solution (see 5.2.)
- Polyphenol-containing solid food (e.g. malt and hops):** weigh 1 g of a representative, homogeneous sample and add 10 ml 60 % ethanol solution, containing 10 % liquid fish gelatin (see 5.2.)
- Liquid food (e.g. soy sauce):** mix 1 ml of the sample with 9 ml 60 % ethanol solution (see 5.2.)
- Beer:** mix 1 ml of the sample with 9 ml 60 % ethanol solution, containing 10 % fish gelatine (see 5.2.)

Further procedure for all samples

- mix thoroughly for at least 30 sec (vortex)
- shake well up side down or rotate on a rotator for further 10 min
- centrifuge: 10 min / at least 2500 g / room temperature (20 - 25 °C / 68 - 77 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the supernatant (or filtrate) 1:50 (1+49) with diluted sample diluent (i. e. 20 µl supernatant + 980 µl diluted buffer)
- use 50 µl per well in the assay

Remark:

All supernatants obtained after centrifugation can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °F) up to four weeks.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **antibody enzyme conjugate** (bottle with red cap) is provided as a concentrate (11fold). Since the diluted enzyme conjugate solution has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the conjugate concentrate should be shaken. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with distilled water (e. g. 100 µl conjugate concentrate + 1 ml water, sufficient for 2 microtiter strips) mix well. Take care that the water is not contaminated with gliadin.

The **wash buffer** is provided as a 10fold concentrate. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard or prepared sample to separate duplicate wells.
3. Add 50 µl of diluted conjugate (see 10.1.) to each well, mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.

6. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate, enclosed in the test kit.

The gliadin concentration in ng/ml is read from the standard curve and must be further multiplied by the dilution factor of 500. This factor of 500 is set by default in the RIDA[®]SOFT Win.

The RIDA[®]SOFT Win (Version 1.93 or newer) indicates the results in gliadin and gluten. The R5 antibody detects gliadin, secalin and hordein.

To convert gliadin into gluten the factor of 2 used (definition Codex Alimentarius). This factor depends strongly on the analyzed sample, especially washed starches may contain a higher proportion of glutelins, in this case the factor is bigger than 2 and the gluten concentration accordingly higher.

Example:

From the standard curve a concentration of 50 ng/ml (ppb) gluten is obtained. Multiplying by the recommended dilution factor 500 leads to 25000 ng/kg gliadin, corresponding to 25 mg/kg (ppm) gliadin. To calculate the gluten content, it is necessary to multiply by factor 2 which results in 50 mg/kg (ppm) gluten).

Background information:

The table below shows the range of variation for the analysis of hydrolyzed samples (further method performance results are published in J. AOAC. Int. 98:1346-1354, 2015).

Symbol	Beer			Starch syrup		Sourdough	
	gluten-free	15 mg/kg prolamin	50 mg/kg prolamin	gluten-free	naturally contaminated	35 mg/kg prolamin	75 mg/kg prolamin
Total number of laboratories	13	12	11	13	13	13	13
Total number of replicates	26	24	22	26	26	26	26
Overall mean	1.2	13.1	59.7	0.7	5.3	24.2	72.8
Recovery (%)		87	119			69	97
Repeatability SD* (mg/kg) s(r)	1.2	4.0	18.6	1.0	0.9	5.6	14.2
Reproducibility SD* (mg/kg) s(R)	1.5	4.8	18.6	1.5	1.8	6.3	20.0
Repeatability rel. SD* (%) RSD(r)	97.9	30.2	31.2	157.2	16.3	23.1	19.5
Reproducibility rel. SD* (%) RSD(R)	126.1	36.9	31.2	236.1	34.4	25.9	27.5
Horrat value	8.1	3.4	3.6	13.8	2.8	2.6	3.3

- The pepsin/trypsin digest used as a calibrator in this assay does not represent all hydrolysis processes. It is recommended that users confirm method performance for their specific processes.
- Samples tested negative still could contain a gluten contamination below the limit of detection of the assay.
- Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.
- Heat treated samples that are extracted with ethanol show a reduced recovery. For such sample the extraction should be carried out with Cocktail (patented) followed by analysis in the sandwich ELISA RIDASCREEN® Gliadin.

Recommendation:

In order to ensure a high analytical performance we recommend:

- to analyze each sample material in duplicates
- to use gluten-free and gluten-containing (spiked) samples as test controls
- to carry out spiking experiments for an accurate and correct procedure
- to perform SureFood® PCR for confirmation of the result

Further application notes:

- Analysis of Gluten in supplemental enzymes

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG