Inositol

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VitaFast® Inositol

Brief information

Easy to use microbiological microtiter plate test for the quantification of total inositol (added and native inositol) in food, animal feed and pharmaceutical products. All required reagents and standard are contained in the test kit. The kit is sufficient for 96 determinations incl. standard. Evaluation requires a microtiter plate reader (610 - 630 nm alternatively at 540 - 550 nm).

Sample preparation

Liquid and solid samples (native and added inositol):

Sample homogenization, enzymatic extraction,
centrifugation and dilution

Time requirement: Test conductionapprox. 60 min
Evaluation.....2 min

Incubation: 44 - 48 h in the dark at 30 °C (86 °F)

Standard range: 0.5 - 5.0 mg / 100 g

Recovery: 80 - 120 %

Intra-assay CV for standards: < 10 %

Inter-assay CV for standards: < 10 %

1. Principle of the test

The VitaFast[®] Inositol microtiter plate test is a microbiological method for the quantitative determination of total inositol (added and native inositol) in food, animal feed and in pharmaceutical products. The microbiological test system is in accordance with the international norms.

Inositol is extracted from the sample and the extract is diluted. The diluted extract and the inositol assay medium are pipetted into the wells of a microtiter plate which is coated with *Saccharomyces cerevisiae*. The growth of *Saccharomyces cerevisiae* is dependent on the supply of inositol. Following the addition of inositol as a standard or as a compound of the sample, the bacteria grow until the inositol is consumed. The incubation is done in the dark at 30 °C (86 °F) for 44 - 48 h. The intensity of metabolism or growth of *Saccharomyces cerevisiae* in relation to the extracted inositol is measured as turbidity and compared to a standard curve. The measurement is done using a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm).

2. Reagents provided

Each test kit contains sufficient materials for 96 determinations incl. standards. Each test kit contains:

- 1 x microtiter plate with 96 wells, coated with *Saccharomyces cerevisiae*
- 3 x redist., sterile water (30 ml) for the preparation of assay medium, standards and dilution of extracted samples
- 3 x inositol assay medium (solid)
- 3 x inositol standard (solid)
- 3 x adhesive foil (1 complete and 2 halves of adhesive foils, sufficient for 3 test runs)
- 1 x additional holder for microtiter strips

Remark:

No quality guarantee can be assumed after expiration of the shelf life.

3. Required reagents and instruments, not provided

- sterile bench (recommended for sterile working)
- microtiter plate reader 610 - 630 nm (540 - 550 nm)
- incubator with dark chamber, 30 °C (86 °F)
- incubator with dark chamber, 37 °C (98 °F)
- water bath heatable 95 °C (203 °F)
- pH meter
- centrifuge > 8000 x g (if the sample cannot be filtrated)
- sterile tips for graduated micropipette 20 - 200 µl and 100 - 1000 µl
- sterile graduated centrifuge vials with screw cap (15 and 50 ml) and sterile reaction vials 1.5 or 2.0 ml
- sterile filters polyethersulfon 0.2 µm with syringe
- redist. or deionized water for sample extraction
- HCl 1.0 mol/l and 0.1 mol/ l
- NaOH 1.0 mol/ l and 0.1 mol/ l
- NaOH concentrated (40 g NaOH ad 100 ml redist. or deionized water)
- NaOH 2 mol/ l (8 g NaOH ad 100 ml redist. or deionized water)

Reagents for the sample extraction

- redist. or deionized water for sample extraction
- taka diastase, aspergillus oryzae; (e.g. Fluka 86250)
- Pepsin from pork (Merck 1.07185.0100)
- citrate buffer pH 4.5 (weigh 1.5 g citric acid monohydrate (e.g. Roth 5110.3) in a 100 ml beaker with magnetic stirrer; solve the citric acid with about 50 ml redist. or. deionised water under stirring; thereafter add 12 ml NaOH 1 mol / l (or 0.48 g NaOH); the pH should be 4.5 (correct with HCl 0.1 mol / l); transfer the solution quantitatively with redist. or deionised water in a 100 ml volumetric flask and fill up with redist. or deionised water to the mark; the buffer can be stored 3 days at 2 - 8 °C (35.6 - 46.4 °F)

4. Warning and precautions for the user

- the assay medium could evoke irritations of mucosa, eyes and skin
- after running the test the strips used must be disposed of according to regulations (e.g. autoclaved)

5. Storage instructions

Store the kit / reagents at 2 - 8 °C (35.6 - 46.4 °F).

Use the prepared reagents (standard, medium) directly and reject them after the assay.

6. Sample preparation

For the determination of the **total inositol** content, the sample has to be treated by enzymatic extraction.

Samples should be stored protected from light at 4 °C (39.2 °F). Standard and samples should be run in triplicate (note: official norms recommend triplicate determination). **Unknown sample matrices** should be analyzed **with two dilutions** of the sample extract. Sample extracts have to be used within one day and should be stored in the dark until analyzing.

Sample extraction is carried out with 1 g homogenized sample in 40 ml extraction solution. This equals a **sample extraction dilution factor of 40**. This factor is already included in the standard curve (see Quality Assurance Certificate). **For low inositol concentrations a sample weight of up to 5 g can be used** (this has to be considered in the evaluation).

Only sterile sample extracts or sterile dilutions thereof should be pipetted onto the microtiter plate. Dilutions have to be prepared with sterile water from the test kit. **Therefore after the sample extraction sterile working conditions and sterile consumables are necessary.** A sterile filtration of the sample or the sample extract is **always** necessary for:

- samples containing herbs and spices
- samples with low concentrations of inositol that are highly coloured; the filtration step eliminates the colouring
- if filtration is not possible due to solid particles or due to cloudiness, centrifugation should be carried out before the sterile filtration step (greater than 8000 x g for 5 min)

The sterile filtration of the sample is not necessary if the sample extraction is carried out at 95 °C (203 °F) for 30 min. Nevertheless the extracted samples have to be diluted with sterile water from the test kit (the assay medium has to be always filtered).

Example for the dilution factors of the sample extract

Solid sample with a labeled concentration of 50 mg / 100 g

The dilution of the extract should be in the middle of the standard curve. Therefore, the concentration is divided by standard 2.

Calculation:

50 mg / 1.0 mg = 50
→ dilution factor 50 (1 : 50)

Dilution steps:

- a) 1 : 10 → 0.1 ml sample extraction solution + 0.9 ml sterile water from the test kit
- b) 1 : 5 → 0.2 ml from a) + 0.8 ml sterile water from the test kit

6.1. Total inositol content (native and added inositol)

Weigh exactly 1 g (ml) homogenized sample into a 50 ml centrifuge vial, add 20 ml redist. or deionized water. Shake and adjust pH to 4.5 with HCl.

Alternatively instead of water a citrate buffer can be used for the extraction (no pH adjustment is necessary); to 1 g sample add 20 ml citrate buffer pH 4.5 and shake.

Add 100 mg taka diastase and 100 mg pepsin, shake well and incubate for 1 h at 37 °C (98.6 °F) in the dark (shake at times). Fill up exactly to 40 ml with redist. or deionized water and extract for 30 minutes in a water bath at 95 °C (203 °F). During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vial is tightly closed. Chill down quickly to below 30 °C (86 °F). Thereafter centrifuge and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

7. Test implementation

7.1. Test preparation

The bottle with sterile water: push the coloured lid up, pull off right up to the glass rim and turn entire lid to remove it.

Inositol standards should be dissolved and diluted freshly. Each dilution is sufficient for three wells.

- open the inositol standard bottle, place the lid down with the opening facing upwards
- add **x ml (x = see quality assurance certificate and label standard bottle)** sterile water (from the test kit) to the standard bottle. Close the standard bottle with the lid and dissolve the standard by shaking = **standard concentrate**
- take 6 sterile vials (1.5 or 2.0 ml) and prepare from the dissolved standard concentrate a **standard curve** according to the following scheme:

standard curve* in mg / 100 g (ml)	sterile water in µl		standard concentrate in µl		total volume in µl
blank: 0	950	+	0	=	950
standard 1: 0.5	950	+	50	=	1000
standard 2: 1.0	900	+	100	=	1000
standard 3: 2.0	400	+	100	=	500
standard 4: 3.0	350	+	150	=	500
standard 5: 5.0	250	+	250	=	500

* the sample extraction dilution of 1 : 40 is already included in the standard curve

The **inositol assay medium** is sufficient for 6 microtiter strips. Open the assay medium bottle and remove the desiccant using tweezers (discard the desiccant).

- add 10 ml sterile water from the test kit to the Inositol assay medium bottle
- close the assay medium bottle carefully and shake well
- heat the bottle in a water bath to 95 °C (203 °F) for 5 min while shaking at least twice; always make sure that the bottle is tightly closed
- quickly chill down to room temperature below 30 °C (86 °F)
- filter the medium through a 0.2 µm filter into a sterile 15 ml centrifuge vial

7.2. Test procedure

Only sterile samples which are diluted with sterile water from the test kit should be pipetted onto the microtiter plate.

- remove the **required strips** of the microtiter plate and place them into the additional holder; return the unused strips together with the desiccant to the foil bag and seal it well, store at 2 - 8 °C (35.6 - 46.4 °F)

Pipette first the assay medium and then standard or diluted samples, as followed:

- pipette 150 µl inositol assay medium into the wells
- pipette 150 µl standard or diluted sample into the assigned wells (flush the pipette tip with standard or sample solution)
- cover the strips / cavities with adhesive foil: pull off the protective layer of the foil, place the foil flat onto the strips, smoothing it down by hand, press the foil firmly onto the strips
important: make sure the cavities are sealed airtight by smoothing down the foil over the cavities, take special care with the wells around the edges
- incubate at **30 °C (86 °F)** in the dark for **44 - 48 h** in an incubator

7.3. Measurement

- press down the adhesive foil once more, place the microtiter plate upside down on a table and dissolve the microorganisms thoroughly by shaking the plate on the surface of the desk
- invert the plate to the regular position and remove the adhesive foil diagonally, 180 degrees backwards, starting from the upper right; **the foil is strongly adhesive, so pulling it off the microtiter plate must be done with great care: hold the strips firmly in the frame with one hand while you pull off the foil diagonally from the top right to the back**
- destroy any bubbles on the surface of liquid in the wells (by means of a pipette tip or a needle)
- measure the turbidity with a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm)

Note:

- after 44 - 48 h of incubation, the microtiter plate can be stored for max. 48 h in the refrigerator, thereafter the turbidity should be measured
- to avoid any time losses due to weekends or bank holidays, the microtiter plate can be evaluated after 60 h. It is recommended to use a timer to turn off the incubator after 44 - 48 h

8. Evaluation

A **4-parameter** evaluation is recommended, e. g. the RIDA® SOFT Win from R-Biopharm.

The test evaluation is correct on condition that

- OD blank < OD standard 1
- OD standard 5 > 0.6 OD

The sample dilution factor of 40 is already included in the standard curve (see Quality Assurance Certificate). In the below formula merely the further dilution factor of the extract and a differing sample weight need to be taken into consideration.

$$\text{Inositol (in mg / 100 g)} = \frac{\text{conc. standard curve} \times \text{dilution factor}}{\text{amount of sample in g}}$$

Example:

Sample weight:	1 g
Sample extraction dilution:	1 : 40 (must not be considered)
Dilution of the sample extraction:	1 : 50 (has to be considered)
Measured concentration from the standard curve:	1.0 mg inositol / 100 g

$$1.0 * 50 / 1 = 50 \text{ mg inositol / 100 g}$$

Unknown sample matrices should be analyzed with two dilutions of the extract (the deviation should be lower than 10 %). If the Inositol content of the higher dilution is elevated, inhibitors such as heavy metals and antibiotics could be present.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

9. Literature

- Official Methods of Analysis (AOAC) 960.46
- Schweizer Lebensmittelbuch SLMB Kapitel 62 / 16 - 16.1
- Atkin, Schultz, Williams and Frey. 1943. End. & Eng. Chem., Ann. Ed. 15 : 141.
- Vitamin-Bestimmungen, Verlag Chemie GmbH Weinheim/Bergstr.; Strohecker und Henning
- The vitamins, 2. Auflage, Vol. VII. Academic Press, New York/London (1967): V. Herbert und J.R.Bertino

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is a registered trademark of the ifp Institut für Produktqualität GmbH.