RIDASCREEN[®] FAST Casein

Art. Nr. R4612

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RIDASCREEN[®]FAST Casein

Brief information

RIDASCREEN[®]FAST Casein (Art. No.: R4612) is a sandwich enzyme immunoassay for the quantitative analysis of casein in food, like bakery goods, sausages and baby food.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation (for 10 samples) approx. 20 min
sample preparation (for 10 heated samples) approx. 30 min
test implementation (incubation time) 30 min

Limit of detection: Allergen extraction buffer 0.12 mg/kg (ppm) casein
RIDA[®] Extractor 2 0.71 mg/kg (ppm) casein

Limit of quantification: Allergen extraction buffer 0.5 mg/kg (ppm) casein
RIDA[®] Extractor 2 2.5 mg/kg (ppm) casein

Specificity: The antibodies specifically detect caseins of cow's milk and caseins of sheep's, goat's and buffalo's milk.

There is no cross reactivity to β -lactoglobulin.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis.

The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDA[®] Extractor 2 (R4613)

RIDASCREEN[®] β -Lactoglobulin (R4901)

RIDASCREEN[®]FAST β -Lactoglobulin (R4902)

RIDASCREEN[®]FAST Milk (R4652)

bioavid Lateral Flow Milch / Milk (BL613-10 / -25)

1. Intended use

RIDASCREEN[®]FAST Casein is a sandwich enzyme immunoassay for the quantitative analysis of casein in food like ice cream, wine, chocolate, beverages, infant formula, bakery goods, sausages, cake and bread mix.

2. General

Cow`s milk contains 3.2 % proteins which consist of 10 % β -lactoglobulin (leading protein of whey) and 80 % caseins. The most important allergen for children is β -lactoglobulin while the caseins become to be dominant later in adults. Casein (lat. Caseus = cheese) is a rough flaked curdling protein, which forms micells in the milk and precipitates under acidic conditions.

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011** milk and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to casein. By adding standards and samples to the wells, casein present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of casein takes place by adding Substrate/Chromogen solution. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made

photometrically at 450 nm. The absorbance is proportional to the casein concentration of the sample. The result is expressed in mg/kg casein.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses).

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
1 x Extractor 2	Blue	Concentrate	2x	30 ml
Allergen extraction buffer	Green	Concentrate	10x	100 ml
Additive 1	Blue			2 g
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	0.5 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	1.5 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	4.5 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	13.5 mg/kg	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Concentrate	11x	0.7 ml
Conjugate buffer	Black	Ready to use		7 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

*) **The dilution factor 20** which results after sample preparation for ice cream, wine, chocolate, beverages, infant formula (see 9.1.) has already been considered for the standard concentrations. Therefore, the casein concentrations of the samples can directly be read from the standard curve.

The sample preparation with Extractor 2 (see 9.2. and 9.3) is used, it results in a **dilution factor of 100**. Therefore the concentrations read from the standard curve have to be multiplied additionally by factor 5.

5. Reagents required but not provided

5.1. Equipment:

- scale
- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- water bath
- fluted filter
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- graduated pipettes
- graduated cylinder
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- optionally a further microtiter plate (e.g. low binding from Greiner bio-one Cat.-No. 655101)
- optionally 8-channel pipette for 100µl

5.2. Reagents:

- distilled or deionized water
- 1 M sodium hydroxide (NaOH)
- 1 M hydrochloric acid (HCl)
- for the extraction of walnut, sunflower-, pumpkinkernels and pine nuts: bovine serum albumin (BSA, Serva, fraction V, protease free, Art. No. 11926.01)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The Extractor 2 is harmful to health. It contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

This kit may contain further hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The Substrate/Chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiry date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained Substrate/Chromogen solution prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5
- odorless Extractor 2 after repeated opening of the bottle

9. Preparation of samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of casein and to avoid contamination.

The **Allergen Extraction buffer** is provided as a **10fold concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. twelve weeks.

For the preparation of the **Allergen Extraction buffer containing Additive 1 (A-AEP)** weigh 1.35 g of Additive 1 in a glass beaker and add 15 ml 1 M NaOH. Stir until the Additive 1 is solved. Then fill 700 ml diluted Allergen Extraction buffer (see above) in a measuring cylinder. Add the 15 ml Additive 1 solution by stirring constantly, transfer eventually residues of the Additive 1 solution into the measuring cylinder by rinsing with diluted Allergen Extraction buffer. Adjust the Additive 1 containing Allergen Extraction (A-AEP) buffer to pH 9 with 1 M HCl and fill up to 750 ml with diluted Allergen Extraction buffer.

750 ml A-AEP buffer is sufficient for 45 samples. The buffer can be used for approx. 3 weeks at room temperature 20 - 25 °C (68 - 77 °F) (do **not** store in the refrigerator). Use clean bottles as particles of dust can initiate crystallization. Discard the buffer if crystals are present.

Extractor 2 is provided as 2fold concentrate and has to be diluted 1:2 (1+1) with distilled water (e.g. 30 ml Extractor 2 + 30 ml dist. water). The complete diluted Extractor 2 is sufficient for 15 samples and can be used for approx. 3 month at room temperature (20 - 25 °C). Additional Extractor 2 can be ordered at R-Biopharm with the product number R4613.

All food samples can be extracted according to method 9.2. Food such as ice cream, wine, chocolate, beverages and infant formula can also be extracted according to 9.1. to reduce the limit of quantification.

9.1. Sample extraction with Allergen Extraction buffer (AEP) for food like ice cream, wine, chocolate, beverages, infant formula

Food, which contains walnut, should be extracted after addition of 0.5 g BSA.

- ground well 5 g of the sample (melt chocolate) and mix thoroughly
- weigh 1 g of sample, add 20 ml of diluted Allergen Extraction buffer (the extraction buffer should have already been heated to approx. 60 °C (140 °F))
- or add 19 ml of the diluted Allergen Extraction buffer (the extraction buffer should have already been heated to approx. 60 °C (140 °F)) to 1 ml liquid sample; 1 ml of wine can be extracted with 9 ml diluted Allergen extraction buffer (LOD 0.12 mg/l; LOQ 0.25 mg/l)
- mix intensively and incubate for 10 min at 60 °C (140 °F) by shaking occasionally, let the sample cool down afterwards
- centrifuge: 10 min / 2500 g / if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl of the particle free supernatant or filtrate per well in the assay

9.2. Sample extraction with Extractor 2 and Allergen Extraction buffer containing Additive 1 (A-AEP) for all other samples

- homogenize a representative amount of the sample (5 - 50 g)

solid samples:

- take 1 g and add 4 ml prepared Extractor 2 (see 9.), mix vigorously, close the vial and cook it for 10 min at 100 °C in a water bath
- let the sample cool down shortly
- add 16 ml heated (60 °C) A-AEP (see 9.) to the cooked sample

liquid samples:

- take 1 ml and add 4 ml prepared Extractor 2 (see 9.), mix vigorously, close the vial and cook it for 10 min at 100 °C in a water bath
- let the sample cool down shortly
- add 15 ml heated (60 °C) A-AEP to the cooked sample

further prepare solid and liquid samples as follows:

- mix vigorously (shaker)
- extract for 10 min at 60 °C in a water bath
- cool down (e.g. ice water), centrifuge for 10 min / 2500 g and/or filter (alternative: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed for 10 min in a microcentrifuge)
- dilute the particle free supernatant or the filtrate 1:5 (1+4) with diluted Allergen Extraction buffer, without Additive 1 (see 9.) (samples such as sausage should be diluted with distilled water), e.g. 100 µl + 400 µl (1:100 final)
(**Remark:** use this diluted supernatant immediately in the assay)
- use 100 µl per well in the assay

9.3. Sample extraction with Extractor 2 and Allergen Extraction buffer containing Additive 1 (A-AEP) for sunflower seed-, walnut-, pumpkin seed- and pine nut-containing matrices

- take 1 g and add 0.5 g BSA, 4 ml prepared Extractor 2 (see 9.) and 16 ml heated (60 °C) A-AEP (see 9.), mix vigorously, close the vial and cook it for 10 min at 100 °C in a water bath
- cool down (e.g. ice water), centrifuge for 10 min / 2500 g and/or filter (alternative: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed for 10 min in a microcentrifuge)
- dilute the particle free supernatant or the filtrate 1:5 (1+4) with diluted Allergen Extraction buffer, without Additive 1 (see 9.), e.g. 100 µl + 400 µl (1:100 final)
(**Remark:** use this diluted supernatant immediately in the assay)
- use 100 µl per well in the assay

The sample extracts of 9.1. can be stored at 2 - 8 °C (35 - 46 °F) for 1 day. Extracts of 9.2. and 9.3. can be stored at 2 - 8 °C (35 - 46 °F) for 3 days. Not used extracts from 9.2 and 9.3 can be stored at -20°C (-4°F) for several months.

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as a concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in conjugate buffer (e.g. 200 µl concentrate + 2 ml conjugate buffer, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10fold concentrate. Before dilution of the concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. The concentrate has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific)) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette. It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate wells in duplicate and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.

4. Add 100 µl of the diluted conjugate (see 10.1.) to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 -77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl wash buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of the reddish Substrate/Chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.

Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win.NET (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. Please note the Compliance Criteria, if the curve shape differs from the Certificate.

In comparison with the certificate, higher values of the absorbance ($A_{450 \text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or casein contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450 \text{ nm}}$) > standard 5.

Please note:

For ice cream, wine, chocolate, beverages, infant formula (see 9.1.):

When working in accordance with the regulation stated, the dilution factor is 20. The casein concentration can be read directly from the standard curve (see 4. *) - the sample dilution factor of 20 is already taken into account for the standard concentrations.

For sample dilutions of e.g. 1:10 the dilution factor must be modified to 0.5 for the calculation of the casein concentration.

For extraction with Extractor 2 (see 9.2. and 9.3):

When working in accordance with the regulation stated, the dilution factor is 100. The sample dilution factor of 20 is already taken into account for the standard concentrations and therefore the concentrations read from the standard curve have to be multiplied additionally by factor 5.

For sample dilutions of more than 1:100 the further dilution factor must be considered for the calculation of the casein concentration.

Limitations for the analysis of raw fish:

Raw fish strongly binds casein, therefore the recovery may only be 10 %. This is also true for cooked fish if casein is added before cooking.

In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report.

Recommendation:

In order to ensure a high analytical performance we recommend:

- to analyze each sample material ~~shou~~ld in duplicates
- to use also allergen-free and allergen- containing (spiked) samples as test controls
- in case of extremely acid or alkaline samples, the pH should be adjusted to a neutral pH
- to carry out spiking experiments for an accurate and correct procedure
- For details using the ChemWell[®] or GEMINI automation please contact sales@r-biopharm.de.

Further information and application notes are available on request from your local distributor or R-Biopharm AG.

Further application notes:

- Extraction for sausages using RIDA[®] Extractor 2 (Art. No. R4613)
- Extraction of pine nut, sunflower kernels and pumpkin kernels (addition of BSA)
- Sample preparation for wine (dilution 1:10)
- Sample extraction for heavily liquid absorbing matrices using RIDA[®] Extractor 2 (Art. No. R4613)
- Sample preparation for food samples with the RIDA[®] Extraction Solution (colorless) (Art. No. R7098), that sample preparation results in similar values compared to Extractor 2. Sample extracts with Extractor 2 can also be used both in the RIDASCREEN[®]FAST β -Lactoglobulin (R4902) and in the RIDASCREEN[®]FAST Milk (R4652) assays.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG