RIDASCREEN[®] FAST Lysozym

Art. No. R6452

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RIDASCREEN[®]FAST Lysozym

Brief information

RIDASCREEN[®]FAST Lysozym (Art. No.: R6452) is a sandwich enzyme immunoassay for the quantitative determination of lysozyme (hen's egg protein) in food like wine, cheese and sausage.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation (for 10 samples) approx. 15 min
test implementation (incubation time) 30 min

Limit of detection: Wine0.006 mg/kg (ppm) lysozyme
Cheese, sausage0.016 mg/kg (ppm) lysozyme

Limit of quantification: Wine0.05 mg/kg (ppm) lysozyme
Cheese, sausage0.25 mg/kg (ppm) lysozyme

Specificity: The antibodies specifically detect lysozyme of hen's egg.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. corn bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Produktangebot

bioavid Lateral Ei / Egg (BL608-10)

1. Intended use

RIDASCREEN®FAST Lysozym (Art. No.: R6452) is a sandwich enzyme immunoassay for the quantitative analysis of lysozyme (hen's egg protein) in food like wine, cheese and sausage.

2. General

The allergen can be present as an ingredient or as a contamination in raw and cooked products. According to the **regulation (EU) No. 1169/2011**, egg and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Lysozyme is due to its bacteriostatic effect used for the preservation in wine production as well as in cheese and sausage products (E1105). It is one of the allergenic egg white proteins and is contained to 3.5 % in egg white protein. In addition to ovalbumin, ovomucoid and ovotransferrin lysozyme can lead to allergic reactions.

3. Test principle

The wells of the microtiter strips are coated with specific antibodies to Lysozyme. By adding standards and samples to the wells, any lysozyme present will bind to the specific antibodies. In a washing step components not bound are removed. Then antibody conjugated to peroxidase is added. This antibody conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in a washing step. The detection of lysozyme takes place by adding Substrate/Chromogen. The enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the Lysozyme concentration of the sample. The result is expressed in mg/kg lysozyme.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Allergen Extraction buffer	Green	Concentrate	10x	100 ml
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 ml
Standard 2*	Transparent	Ready to use	0.05 mg/kg	1.3 ml
Standard 3*	Transparent	Ready to use	0.1 mg/kg	1.3 ml
Standard 4*	Transparent	Ready to use	0.2 mg/kg	1.3 ml
Standard 5*	Transparent	Ready to use	0.4 mg/kg	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Ready to use		6 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 ml
Stop solution	Yellow	Ready to use		14 ml

- *) **The dilution factor 20** which results after sample preparation of wine has already been considered for the standard concentrations. Therefore, the lysozyme concentrations of wine samples can directly be read from the standard curve.

In the case of cheese and sausage the **dilution factor 100** results. Therefore, the concentration read from the standard curve has to be multiplied additionally by factor 5.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal vials
- water bath (60 °C / 140 °F)
- fluted paper filter
- variable 20 - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents:

- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

–any bluish coloration of the red stained substrate/chromogen prior to test implementation

–a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of lysozyme and to avoid contamination.

The **Allergen Extraction buffer** is available as a **10x concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate **1:10 (1+9) with distilled water** before use (e.g. 100 ml of buffer concentrate + 900 ml distilled water). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for four weeks.

For the extraction of cheese and sausage samples, a high-salt allergen extraction buffer is needed separately. For this 10 g NaCl is added to 100 ml Allergen Extraction buffer. The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for four weeks.

A representative amount of sample (5 - 50 g) should be homogenized.

9.1. Cheese and sausage

- add 1 g of a representative, homogeneous sample to 20 ml high-salt allergen extraction buffer (preheated to 60 °C / 140 °F) and mix
- incubate samples for 10 min at 60 °C (140 °F) in a water bath
- cool samples to room temperature (20 - 25 °C / 68 - 77 °F) (e.g. in an ice bath)
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39.2 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute supernatant 1:5 (1+4) with finally diluted Allergen Extraction buffer (no salt added)
(**Note:** use this diluted supernatant directly or within 2 hours in the test)
- use 100 µl per well in the test

If further dilutions are necessary, they should be performed with a mixture of Allergen Extraction buffer / high-salt allergen extraction buffer in proportion 5:1 (4 ml + 1 ml).

9.2. Wine

- add 1 ml of a representative, homogeneous sample to 19 ml of Allergen Extraction buffer (preheated to 60 °C / 140 °F) and mix
- incubate samples for 10 min at 60 °C (140 °F) in a water bath
- cool samples to room temperature (20 - 25 °C / 68 - 77 °F) (e.g. in an ice bath)
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39.2 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl per well in the test

To obtain a **lower limit of quantification** carry out a 1:10 dilution instead of a 1:20 dilution. To 1 ml of wine sample add 9 ml of Allergen Extraction Buffer. In the RIDA[®]SOFT Win a dilution factor of 0.5 (instead of 1) has to be taken into account.

When using wine with high **polyphenol content** (e.g. red wine) add 0.25 % casein (e.g. SIGMA, C5890) or 0.25 % gelatin (e.g. SIGMA, G-6144) to the diluted Allergen Extraction buffer. Add e.g. 0.25 g of casein or gelatin to 100 ml final diluted Allergen Extraction buffer and stir for 10 minutes. Add 19 ml of diluted allergen extraction buffer with additive to 1 ml of wine and extract for 10 min at 60°C.

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days or at -20 °C (-4 °F) for several months.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **Wash Buffer** is provided as a 10fold concentrate. Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml distilled water). The diluted wash buuffer is stable at 20 - 25 °C (68 - 77 °F) for four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).

3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat three more times.
4. Add 100 µl of the conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 -77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat three more times.
6. Add 100 µl of the reddish Substrate/Chromogen to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or lysozyme contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5.

Please note:

When working in accordance with the regulation stated for wine, the dilution factor is 20. The allergen concentration can be read directly from the standard curve - the sample dilution factor of 20 is already taken into account for the standard concentrations (see 4. *).

For sample dilutions of more than 1:20 (sample preparation of cheese and sausage) the further dilution factor must be considered for the calculation of the

lysozyme concentration. In case of cheese and sausages the additional factor is 5.

In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food, proteins may be altered or fragmented, this may have an impact on the recovery.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report.

Recommendation:

In order to ensure a high analytical performance we recommend:

- to analyze each sample material in duplicates
- to use also allergen-free and allergen- containing (spiked) samples as test controls
- in case of extremely acid or alkaline samples, the pH should be adjusted to a neutral pH
- to carry out spiking experiments for an accurate and correct procedure
- For details using the ChemWell® or GEMINI automation please contact R-Biopharm AG

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG