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Description

The SureFast® Parasitic Water Panel 4plex is a real-time PCR for the direct, qualitative detection and differentiation of *Giardia intestinalis*, *Entamoeba histolytica* and *Cryptosporidium* spp.. Each reaction contains an internal amplification control. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of four fluorescence emissions at 510 nm, 580 nm, 610 nm and 660 nm (FAM, VIC, ROX and Cy5) at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Roche LightCycler® 480 II/Roche cobas z 480 Analyzer² and Applied Biosystems 7500, BioRad CFX 96.

Limit of Detection

The SureFast® Parasitic Water Panel 4plex real-time PCR has a limit of detection of ≤ 5 DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

DNA-preparation

For DNA-preparation the use of SureFast® PREP Aqua is recommended.

Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- real-time PCR instrument, equipped with four detection channels (510 nm, 580 nm, 610 nm, 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

² For the use of the Roche LightCycler® 480 instrument or cobas z 480 Analyzer a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

2. Setup

	Blockcycler	Rotorcycler/ LightCycler® 480 II
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of Extension Phase Detection system <i>Giardia intestinalis</i> : various devices FAM-channel, Quencher: BHQ LC480 II 465 nm, 510 nm PCR inhibition control: various devices VIC/HEX-channel, Quencher: BHQ LC480 II 533 nm, 580 nm Detection system <i>Entamoeba histolytica</i> : various devices ROX-channel, Quencher: BHQ LC480 II 533 nm, 610 nm Detection system <i>Cryptosporidium</i> spp.: various devices Cy5-channel, Quencher: BHQ LC480 II 618 nm, 660 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

3. Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells.
- Pipette 5 µl of the Positive Control into the designated tubes/wells.
- Centrifuge all tubes/wells shortly at low speed.
- Place tubes/wells into the real-time PCR instrument and start the run according to the setup.

Interpretation of results

The *Giardia intestinalis* DNA is detected in the FAM-channel, the *Entamoeba histolytica* DNA is detected in the ROX-channel, the *Cryptosporidium* spp. DNA is detected in the Cy5-channel and the internal amplification control is detected in the VIC/HEX-channel.

A sample is stated **positive** for the respective parameter (*Giardia intestinalis*/*Entamoeba histolytica*/*Cryptosporidium* spp.), if the sample DNA shows amplification in the respective channel (FAM/ROX/Cy5).

A sample is stated **negative** for the respective parameter (*Giardia intestinalis*/*Entamoeba histolytica*/*Cryptosporidium* spp.), if the sample DNA shows no amplification in the respective channel (FAM/ROX/Cy5) and the internal amplification control (VIC/HEX) of the sample is **positive** (see also table below).

result in the respective channel				
FAM-channel <i>Giardia intestinalis</i>	ROX-channel <i>Entamoeba histolytica</i>	Cy5-channel <i>Cryptosporidium</i> spp.	VIC/HEX-channel amplification control	result
positive	positive/negative	positive/negative	positive/ negative	<i>Giardia intestinalis</i> DNA detected
positive/negative	positive	positive/negative	positive/ negative	<i>Entamoeba histolytica</i> DNA detected
positive/negative	positive/negative	positive	positive/ negative	<i>Cryptosporidium</i> spp. DNA detected
negative	negative	negative	negative	invalid*

* If the sample DNA and the internal amplification control are **negative** the sample contains PCR-inhibiting substances or the sample preparation was not successful. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

Product Information

- Validation Report

Technical Support

For further questions please send an e-mail to info@congen.de.