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**Description**

The test detects *Escherichia coli* DNA. Each reaction contains an internal amplification control (PLUS). The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of two fluorescence emissions at 510 nm and 580 nm (FAM and VIC/HEX) at the same time. The technical validation of instruments was performed on Agilent Mx3005P, BioRad CFX 96, Roche LightCycler® 480 II, Roche LightCycler® 480 Cobas, Roche LightCycler® 2.0<sup>2</sup>, Applied Biosystems 7500, LTF MyGo Pro and Qiagen RotorGene Q.

**Limit of Detection**

The SureFast® Escherichia coli PLUS real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

**Cross reactivity**

Cross reactivity was observed with DNA extracts from *Shigella* spp..

**DNA-preparation**

For DNA-preparation the use of SureFast® PREP Bacteria or SureFast® PREP Aqua is recommended. To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference  $>3$ ).

**Kit components and storage**

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at  $-20^{\circ}\text{C}$  and protected from light.

**Additionally required equipment and materials**

- real-time PCR instrument, equipped with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

**Protocol****1. Preparation of the master-mix**

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at  $-20^{\circ}\text{C}$  and not be mixed by vortexing.

<sup>2</sup> note: For the use of the Roche LightCycler® 2.0 a Color Compensation is necessary. The SureCC Color Compensation Kit II (Art. No. F4010) must be used for the color compensation of such devices.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

## 2. Setup

	Blockcycler	Rotorcycler, LightCycler® 480
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplarisch)	Detection: End of Extension Phase  <u>Detection System:</u> Various devices FAM-channel, Quencher: BHQ Rotor-Gene Q Green LightCycler 2.0 530 nm - none LightCycler 480 II 465 nm - 510 nm  <u>Internal Amplification Control :</u> Various devices VIC/HEX-channel, Quencher: BHQ Rotor-Gene Q Yellow LightCycler 2.0 560 nm - none LightCycler 480 II 533 nm - 580 nm  Passive Reference: None	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/download">http://www.congen.de/en/company/download</a>		

## 3. Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tubes/wells or capillaries of the negative control (it is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the real-time PCR instrument and start the run according to the setup.

**Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows an amplification in the detection system (FAM).

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system and the internal amplification control (inhibition control) of the sample is **positive** (VIC/HEX). If the sample DNA and the internal amplification control are **negative** the sample contains inhibiting substances. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved.

**Product Information**

- Validation Report

**Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).