RIDASCREEN[®] FAST DON

Art. No. R5901 (96 wells)

Art. No. R5902 (48 wells)

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RIDASCREEN®FAST DON

Brief information

RIDASCREEN®FAST DON (Art. No. R5901, 96 wells / R5902, 48 wells) is a competitive enzyme immunoassay for the quantitative determination of deoxynivalenol in cereals, malt and feed. Wheat, barley, malted barley, oats and corn are approved by AOAC Research Institute according to the Performance Tested Method Program and FGIS (Federal Grain Inspection Services) program of the Grain Inspection, Packers and Stockyards Administration of the United States Department of Agriculture (USDA/GIPSA).

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 or 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: extraction and filtration

Time requirement: sample preparation (for 10 samples)approx. 10 min
test implementation (incubation time)8 min

Limit of detection: < 0.2 mg/kg (ppm)

Limit of quantification: 0.2 mg/kg (ppm) / oats: 0.36 mg/kg (ppm)

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDASCREEN® FAST DON SC (R5905)

RIDASCREEN®DON (R5906)

RIDA®QUICK DON (R5904)

TRILOGY® Liquid Standard DON (TSL-317)

TRILOGY® Dried Standard DON (TS-310, TS-317)

1. Intended use

RIDASCREEN®FAST DON is a competitive enzyme immunoassay for the quantitative analysis of DON in cereals, malt and feed.

2. General

Deoxynivalenol belongs to the trichothecene group of mycotoxins and is formed by fungi of the genus *Fusarium*. Deoxynivalenol often occurs in plant products particularly in cereals. Of the trichothecene mycotoxins deoxynivalenol, 3-acetyl- and 15-acetyl-deoxynivalenol are the most frequently occurring toxins in Europe and Northern America. The toxin concentrations found in wheat, corn or rice are often in the ppm range. Due to their high cytotoxic and immunosuppressive properties these toxins pose a risk to human and animal health.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-deoxynivalenol antibodies. Deoxynivalenol standards or sample solutions, deoxynivalenol enzyme conjugate and anti-deoxynivalenol antibodies are added. Free deoxynivalenol and deoxynivalenol enzyme conjugate compete for the deoxynivalenol antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-deoxynivalenol antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the deoxynivalenol concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for as many as 91 (R5901) or 43 (R5902) analyses (plus 5 standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate M	-	Ready to use		96 wells (R5901) 48 wells (R5902)
Standard 1*	white	Ready to use	0 mg/L	1.3 ml
Standard 2*	white	Ready to use	0.222 mg/L	1.3 ml
Standard 3*	white	Ready to use	0.666 mg/L	1.3 ml
Standard 4*	white	Ready to use	2 mg/L	1.3 ml
Standard 5*	white	Ready to use	6 mg/L	1.3 ml
Wash buffer salt Tween		Dissolve the salt		
Conjugate	red	Ready to use		6 ml (R5901) 3 ml (R5902)
Antibody	black	Ready to use		6 ml (R5901) 3 ml (R5902)
Substrate/Chromogen Red Chromogen Pro	brown	Ready to use		10 ml
Stop solution	yellow	Ready to use		14 ml

*) The dilution factor 20 for the sample has already been considered. Therefore, the deoxynivalenol concentrations of samples can be read directly from the standard curve.

5. Materials required but not provided

5.1. Equipment

- microtiter plate spectrophotometer (450 nm)
- graduated cylinder: (plastic or glass) 100 ml, 1 L
- glassware for preparing sample extract: filter funnel and 50 ml flask
- grinder (mill)
- Optional: Ultra-Turrax or equivalent
- optional: shaker
- filter paper: Whatman No. 1 or equivalent
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents

- distilled or deionized water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain deoxynivalenol. Particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and deoxynivalenol solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). **Do not freeze any test kit components.**

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

The test kit can be regularly used at least up to the expiry date (indicated on the kit package), if stored correctly.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected from light.

Bring all reagents and samples to room temperature (20 - 25 °C / 68 - 77 °F) before use.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- weigh 5 g of ground sample into a suitable container and add 100 ml of distilled water *)
- blend the sample by ultra-turrax (or equivalent) for two minutes or shake vigorously for three minutes (manually or with shaker)
- filter the extract through Whatman No. 1 filter (or equivalent)
- use 50 µl of the filtrate per well in the test

*) sample size may be increased if required, but the volume of water must be adapted accordingly, e.g.: 25 g in 500 ml of distilled water or 50 g in 1000 ml of distilled water

USDA/GIPSA extraction method

- weigh 50 g of ground sample into a suitable container and add 250 ml of distilled water
- blend the sample by ultra-turrax (or equivalent) for two minutes or shake vigorously for three minutes (manually or with shaker)
- filter the extract through Whatman No. 1 filter (or equivalent)
- dilute the filtered sample extract 1:4 (1+3) with distilled water (e.g. 1 ml of the extract + 3 ml of distilled water)
- use 50 µl of the diluted filtrate per well in the test

10. Test implementation

10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. The specific reaction starts with the addition of the specific antibody. Not more than three strips should be applied to the test, when a single step pipette is used. More strips (up to 6) could be applied, when a multistep pipette is used.
3. Return all reagents to 2 - 8 °C (35 - 46 °F) immediately after use.

The deoxynivalenol standards are provided ready to use. The dilution factor 20 for the sample has been considered when labeling. Therefore, the deoxynivalenol concentration of samples can be read directly from the standard curve.

As **wash buffer** a PBS tween buffer is needed. Please use the buffer salt contained in the kit (see 4.). Dilute the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the envelope in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

The 10fold concentrate expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C / 68 - 77 °F).

10.2. Test procedure

Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in any ELISA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the ELISA test procedure.

Avoid direct sunlight during all incubations. Therefore covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 µl of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 µl of conjugate to each well.

4. Add 50 µl of antibody to each well. Mix gently by shaking the plate manually and incubate for 5 min (+/- 1) at room temperature (20 - 25 °C; 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells (250 µl per well) with wash buffer (see 10.1.). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
6. Add 100 µl of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 3 min (+/- 0.5) at room temperature (20 - 25 °C; 68 - 77 °F) in the dark.
7. Add 100 µl of stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win.net (Art. No. Z9996), is available to evaluate the RIDASCREEN[®] enzyme immunoassays.

For single determinations we recommend logit/log evaluation and for double or multiple determinations cubic spline should be used.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the deoxynivalenol concentration [mg/kg].

The deoxynivalenol concentration in mg/kg corresponding to the extinction of each sample can be read from the calibration curve.

12. Limit of detection and limit of quantification

The limit of detection (LOD) values were investigated with repeated measurements (n = 10) of blank samples of various commodities.

LOD = mean concentration + 2 x standard deviation of the measurement
 The corresponding limit of quantification (LOQ) values were calculated as
 LOQ = mean conc. + 10 x standard deviation of the measurement.

values (ppb), n = 10	wheat	barley	malted barley	oats	corn
mean	26.3	4.6	30.7	78.0	27.6
s.d.	14.1	3.8	15.9	28.2	11.5
LOD	54	12	62	134	51
LOQ	167	42	190	360	143

All LOD values were found well below 0.2 ppm. LOQ can be set as to 0.2 ppm (standard 2) except for oats (0.36 ppm).

13. Specificity

1. time required for performing the test (incubation time): 8 minutes
2. measuring range: 0.2 - 6.0 ppm
3. approved matrices: wheat, barley, malted barley, oats, corn
4. additionally validated matrices: wheat flour, wheat midds, wheat bran, sorghum, soy flakes, soy meal, mixed animal feed
5. specificity: the test can not differentiate between DON and 3-acetyl-DON (cross reactivity: 213 %), and has a negligibly low or no cross reactivity to other related substances such as Nivalenol, 15-acetyl-DON, Triacetyl-DON, Triacetyl-Nivalenol, Tetraacetyl-DON, and Fusarenon X

Accuracy and precision data with spiked wheat samples (n = 90, each level)

fortification (ppm)	0.5	1.0	2.5	5.0
mean (ppm)	0.53	1.0	2.4	4.1
s.d. (ppm)	0.1	0.22	0.24	0.44

Comparative testing of incurred wheat samples with HPLC (reference method)

ELISA	mean	0.84	4.02	2.30
(n = 30)	s.d.	0.16	0.33	0.15
<hr/>				
HPLC	mean	0.6	3.88	2.28
(n = 5)	s.d	0.12	0.50	0.16

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG