SureFood® QUANTARD Allergen 40

Art. No. S3301

Version 1.6

Kit components and storage

1x QUANTARD Allergen 40 (2.0 g)

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**Description**

SureFood® QUANTARD Allergen 40 is a laboratory reference material developed for quantification in the range of 1 - 400 mg allergenic substance / kg food.

Specification

The laboratory reference material SureFood® QUANTARD Allergen 40 is based on **corn flour** material. It contains all potentially allergenic food ingredients according to regulation (EU) No 1169/2011 (except for **molluscs, lactose and sulphur dioxide**) in a concentration of **40 mg/kg**.

Content of allergenic substances:

- | | | |
|-----------------------------|---|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| • Gluten containing cereals | : | Wheat flour type 405 |
| • Crustaceans | : | Black tiger shrimps (<i>Penaeus monodon</i>) |
| • Egg | : | Egg powder (hen) |
| • Fish | : | Salmon |
| • Peanut | : | Peanuts (<i>Arachis hypogaea</i>) |
| • Soy | : | Soybeans (<i>Glycine max</i>) |
| • Milk | : | Milk powder (cow) |
| • Nuts | : | Hazelnuts (<i>Corylus avellana</i>), Almonds (<i>Prunus dulcis</i>), Walnuts (<i>Juglans regia</i>), Cashewnuts (<i>Anacardium occidentale</i>), Pecan nuts (<i>Carya illinoensis</i>), Brazil nuts (<i>Bertholletia excelsa</i>), Pistachio (<i>Pistacia vera</i>), Macadamia nuts (<i>Macadamia ternifolia</i>) |
| • Celery | : | Celery seeds (<i>Apium graveolens var. secalinum</i>) |
| • Mustard | : | Mustard seeds, yellow (<i>Sinapis alba</i>) |
| • Sesame | : | Sesame seeds (<i>Sesamum indicum</i>) |
| • Lupin | : | Lupin seeds (<i>Lupinus hartwegii</i>) |

Application

The laboratory reference material SureFood® QUANTARD Allergen 40 is validated in combination with the DNA extraction method SureFood® PREP Advanced, protocol 1, and the SureFood® ALLERGEN QUANT real-time PCR detection methods.

Using the SureFood® PREP Advanced kit it is recommended to perform the extraction according to protocol 1 with 100 mg laboratory reference material. Using an external reference for quantification, the sample and the laboratory reference material must be extracted in the same way. Please follow the protocol of the kit for the lysis conditions, the DNA binding and the purification of DNA.

The laboratory reference material has to be extracted in duplicate. The extractions shall be performed in parallel to the sample.

Further information about the quantification of allergenic substances is given in the manual of the SureFood® ALLERGEN QUANT real-time PCR products.

In case of using alternative methods the sample weight and extraction method in combination with the detection system have to be validated.

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R-Biopharm AG