**SureFood® ALLERGEN Crustaceans**

Art. No. S3612

100 rxn

**User Manual**

**March 2018**

 Inhalt /  Content

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Art. No. S3612

March 2018

6	Qualitative Analysis .....	14
6.1	Protocol .....	14
6.1.1	Preparation of the master-mix .....	14
6.1.2	Preparation of the real-time PCR-mix.....	14
6.2	Interpretation of results .....	15
7	Quantitative Analysis .....	16
7.1	Preparation of the master-mix .....	16
7.2	Preparation of the standard DNA dilutions .....	16
7.3	Preparation of the real-time PCR-mix.....	16
7.4	Interpretation of results .....	17
8	Further information .....	18
8.1	Product Information .....	18
8.2	Technical Support .....	18
8.3	Distribution and ordering .....	18

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## 5 General information

### 5.1 Description

The test detects crustaceans DNA according to directive (EC) 1169/2011 qualitatively and / or quantitatively. Each reaction contains an internal amplification control. If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). Also spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

For the quantitative determination the use of the laboratory reference material SureFood® QUANTARD Allergen 40 containing 40 mg crustaceans / kg food sample and a detection system including a standard dilution is required. The real-time PCR assay can be used with established real-time PCR instruments capable of detection of at least two fluorescence dyes. Internal validation was performed on Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, AriaDx and Agilent Mx3005P.

#### Note:

In general this procedure is applicable for all homogenous food samples. The quantitative determination is inapplicable for swabs, sponges and inhomogeneous fluids/liquids. The method is validated for quantitative determination of 1 mg to 400 mg allergenic substance / kg food sample using the reference material SureFood® QUANTARD Allergen 40 containing 40 mg allergenic substance / kg food sample.

The SureFood® ALLERGEN Crustaceans PCR has a limit of detection of  $\leq 0.4$  mg / kg. The limit of quantification is 1 mg / kg using SureFood® PREP Advanced, protocol 1. The limit of detection and limit of quantification were determined using the SureFood® QUANTARD Allergen 40 (matrix: corn flour).

The limit of detection of the complete quantitative determination (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

### 5.2 DNA-preparation

For DNA-preparation the use of SureFood® PREP Advanced, protocol 1 is recommended.

### 5.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
4	Dilution Buffer	1 x 2000 µl	White
5	Standard DNA	1 x 80 µl	Dark Blue

**Store all reagents at -20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple use on the same day.**

# SureFood® ALLERGEN Crustaceans (100 React.)

Art. No. S3612

March 2018

## 5.4 Additionally required equipment and materials

- SureFood® QUANTARD Allergen 40 Art. No.: S3301
- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

## 5.5 Setup

	Blockcycler / Bio Molecular Systems MIC	Rotorcyclcr
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

## 5.6 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Agilent Mx3005P</b>	Crustaceans	FAM	+	
	IAC	HEX	+	
<b>Agilent AriaDx</b>	Crustaceans	FAM	+	
	IAC	HEX	+	
<b>Applied Biosystems 7500</b>	XXXX	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
<b>Bio-Rad CFX96</b>	XXXX	FAM	+	
	IAC	VIC	+	
<b>BioMolecular Systems MIC</b>	XXXX	green	+	
	IAC	yellow	+	

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Qiagen Rotor-Gene Q</b>	XXXX	green	+	
	IAC	yellow	+	
<b>Roche LightCycler® 480 II</b>	XXXX	465-510	+	
	IAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	XXXX	465-510	+	
	IAC	540-580	+	

## 6 Qualitative Analysis

### 6.1 Protocol

#### 6.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative crustaceans detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1 positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

**Note: The Taq Polymerase can be stored at +2 to +8 °C for multiple use on the same day.**

#### 6.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.

- Pipette 5 µl of Positive Control into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the real-time PCR instrument and start the run according to the setup.

## **6.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows amplification in the detection system or the obtained Cp value is  $\leq 35$ .

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is  $> 35$ . The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value  $\leq 2$  compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific crustaceans PCR assay.

### **In general:**

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.

## 7 Quantitative Analysis

### 7.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control. Recommended control reactions for the specific PCR assay: negative control, extraction control, standard curve as positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the quantitative crustaceans detection:

4 reactions for the standard curve

2 reactions for controls (1x no-template control, 1x extraction control)

2 reactions for reference material (SureFood® QUANTARD Allergen 40 DNA)

For each sample: at least 1 reaction for each sample DNA

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

**Note: The Taq Polymerase can be stored at +2 to +8 °C for multiple use on the same day.**

### 7.2 Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different concentrations. Prepare a dilution series of 4 steps. Prepare 4 reaction tubes (labelled S1 to S4) and add 45 µl Dilution Buffer each. The following procedure is recommended:

Standard	dilutions	copy number per µl	final copy number per reaction*
S1	45 µl Dilution Buffer + 5 µl Standard DNA	10.000 copies	50.000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	1.000 copies	5.000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	100 copies	500 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	10 copies	50 copies

**\*Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

If the diluted DNA standards (S1 to S4) are not immediately used, store them at -20°C. The dilution series are stable up to two months at -20°C. Before use allow the reagents to thaw, mix them and centrifuge before use.

### 7.3 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.

- Pipette 5 µl of the SureFood® QUANTARD 40 Allergen DNA and the standard dilutions into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the real-time PCR instrument and start the run according to the setup.

#### 7.4 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

A sample is stated **positive and can be quantified**, if the sample DNA shows an amplification in the detection system and the obtained Cp value is  $\leq 35$ .

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is  $> 35$ . The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value  $\leq 2$  compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific crustaceans PCR assay.

The calculation of mg allergenic substance/kg food sample can be done for samples showing an amplification curve in the detection system with Cp values  $< 35$ . Mark the standards, controls and samples and make the evaluation according to the analysis program recommended by the real-time PCR instrument manufacturer. The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient  $R^2 > 0.98$ . In case of different values for the standard curve, it should not be used for calculation.

#### In general:

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.

For a final calculation of mg allergenic substance / kg food sample following formula is introduced:

$$X \left[ \frac{mg}{kg} \right] = \frac{10 \left( \frac{Cp_{sample}}{s} \right)}{10 \left( \frac{Cp_{QUANTARD}}{s} \right)} * c_{QUANTARD}$$

s = slope  
c = concentration  
Cp<sub>QUANTARD</sub> = Ø of double assay

Example:

$$X \left[ \frac{mg}{kg} \right] = \frac{10 \left( \frac{28.5}{-3.251} \right)}{10 \left( \frac{29.2}{-3.251} \right)} * 40 = 65.7 \frac{mg}{kg}$$

Cp<sub>sample</sub> : 28.5  
Cp<sub>QUANTARD</sub>: Ø 29.2  
s: - 3.251  
c: 40 mg/kg

For this example an amount of **65.7 mg crustaceans / kg food sample** is calculated.

## 8 Further information

### 8.1 Product Information

- **Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices**  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report
- Product Information

### 8.2 Technical Support

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

### 8.3 Distribution and ordering

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)