**CONGEN**

# **SureFood<sup>®</sup> Oat**

Art. No. S7004

100 rxn

## **User Manual**

**Juli 2018**

July 2018

** Inhalt /  Content**

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## 4 General Information

### 4.1 Description

The test detects oat DNA (*Avena sativa*). Each reaction contains an internal amplification control. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of two fluorescence emissions at 510 nm and 580 nm (FAM and VIC/HEX) at the same time. The technical validation of instruments was performed on Applied Biosystems 7500, Agilent Mx3005P, Agilent AriaDx, Bio-Rad CFX96, Roche LightCycler· 480 II, Roche cobas· z 480 Analyzer, Qiagen Rotor-Gene Q and LTF MyGo Pro.

### 4.2 Limit of Detection

The SureFood® Oat real-time PCR has a limit of detection of  $\leq 500$  DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 4.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

### 4.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at  $-20^{\circ}\text{C}$  and protected from light. The Taq Polymerase can be stored at  $+2$  to  $+8^{\circ}\text{C}$  for multiple use on the same day.**

### 4.5 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Agilent Mx3005P</b>	Oat	FAM	+	
	IAC	HEX	+	
<b>Agilent AriaDx</b>	Oat	FAM	+	
	IAC	HEX	+	
<b>Applied Biosystems 7500</b>	Oat	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
<b>Bio-Rad CFX96</b>	Oat	FAM	+	
	IAC	VIC	+	
<b>BioMolecular Systems MIC</b>	Oat	green	+	
	IAC	yellow	+	

July 2018

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>LTF MyGo Pro</b>	Oat	FAM	+	
	IAC	VIC	+	
<b>Qiagen Rotor-Gene Q</b>	Oat	green	+	
	IAC	yellow	+	
<b>Roche LightCycler® 480 II</b>	Oat	465-510	+	
	IAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	Oat	465-510	+	
	IAC	540-580	+	

#### 4.6 Additionally required equipment and materials

- real-time PCR instrument, equipped with two detection channels (510 nm and 580 nm)
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

#### 4.7 Setup

	<b>Blockcycler/ Bio Molecular Systems MIC</b>	<b>Rotorcycler</b>
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	35	35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

## 5 Qualitative Analysis

### 5.1 Protocol

#### 5.1.1 Preparation of the master-mix

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

#### 5.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries and close them.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the PCR instrument and start the run according to the setup.

#### 5.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows an amplification signal in the oat system. A sample is stated **negative**, if the sample DNA shows no amplification in the detection system and the internal amplification control (inhibition control) of the sample is **positive** (VIC/HEX). If the sample DNA and the internal amplification control are **negative** the sample contains PCR inhibiting substances. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

## 6 Further Information

### 6.1 Product Information

- Validation Report

### 6.2 Technical Support

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

### 6.3 Distribution and ordering

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)