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
## 1 General Information

### 1.1 Description

The SureFast® Salmonella Species/Enteritidis/Typhimurium 4plex (R&D Version) is a real-time PCR for the direct, qualitative detection and differentiation of *Salmonella* spp., *Salmonella* Enteritidis and *Salmonella* Typhimurium.

Each reaction contains an internal amplification control. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of four fluorescence emissions at 522 nm, 553 nm, 610 nm and 670 nm (FAM, VIC, ROX and Cy5) at the same time.

### 1.2 Limit of Detection

The SureFast® Salmonella Species/Enteritidis/Typhimurium 4plex (R&D Version) real-time PCR has a limit of detection of  $\leq 5$  DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 1.3 Cross Reactivity

Cross reactivity was observed with DNA extracts from *Salmonella* Amsterdam, *Salmonella* Blegdam, *Salmonella* Moscow and *Salmonella* Nitra in the detection system for *Salmonella* Enteritidis.

### 1.4 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria or SureFast® Speed PREP is recommended. To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing (bacterial growth at  $C_p$  difference  $>3$ ).

### 1.5 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at  $-20^{\circ}\text{C}$  and protected from light.

### 1.6 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Bacteria or SureFast® Speed PREP)
- real-time PCR instrument with four detection channels (522 nm, 553 nm, 610 nm and 670 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps) pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

### 1.7 Setup

	<b>Blockcycler/MIC</b>	<b>Rotorcycler/ LightCycler® 480 II</b>
Initial Denaturation (HOLD) Cycles Denaturation Annealing/Extension (CYCLE)	1 min, 95°C 45 15 sec, 95°C 30 sec, 60°C	1 min, 95°C 45 10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of Extension Phase  Detection System <i>Salmonella</i> spp.: Various devices <b>FAM</b> -channel, Quencher: BHQ LightCycler 480 II 465 nm - 510 nm  Internal Amplification Control: Various devices <b>VIC/HEX</b> -channel, Quencher: BHQ LightCycler 480 II 533 nm - 580 nm  Detection System <i>Salmonella</i> Enteritidis: Various devices <b>ROX</b> -channel, Quencher: BHQ LightCycler 480 II 533 nm - 610 nm  Detection System <i>Salmonella</i> Typhimurium: Various devices <b>Cy5</b> -channel, Quencher: BHQ LightCycler 480 II 618 nm - 660 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/downloads">http://www.congen.de/en/company/downloads</a>		

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

<b>Components for master-mix</b>	<b>Amount per reaction</b>	<b>10 reactions (with 10% excess)</b>
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

## 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

*Salmonella* spp. is detected in the FAM-channel, *Salmonella* Enteritidis is detected in the ROX-channel and *Salmonella* Typhimurium is detected in the Cy5-channel (see table). The internal amplification control will be detected in the VIC/HEX-channel.

A sample is stated **positive**, if the sample DNA shows an amplification in the respective channel. A sample is stated **negative**, if the sample DNA shows no amplification in the respective channel and the internal amplification control (VIC/HEX-channel) of the sample is **positive**.

If the sample DNA and the internal amplification control are **negative** the sample contains PCR-inhibiting substances or the sample preparation was not successful. Under these circumstances an evaluation of the results is not possible. DNA isolation and purification for the sample need to be improved.

result in the respective channel				VIC/HEX channel amplification control	result
FAM channel <i>Salmonella</i> spp.	ROX channel <i>Salmonella</i> Enteritidis	Cy5 channel <i>Salmonella</i> Typhimurium			
<b>positive</b>	negative	negative	<b>positive</b>	<i>Salmonella</i> spp. DNA detected	
negative	<b>positive</b>	negative	<b>positive</b>	<i>Salmonella</i> Enteritidis DNA detected	
negative	negative	<b>positive</b>	<b>positive</b>	<i>Salmonella</i> Typhimurium DNA detected	
negative	negative	negative	negative	invalid	

## 3 Further Information

### 3.1 Product Information

- Validation Report

### 3.2 Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de)