---

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# QuickGEN PCR Kit

## Beer yeast and bacteria differentiation

### 1. Intended use

Differentiation of beer spoilage bacteria and yeast in beer and beer mixing drinks.

The following bacteria and yeasts are differentiated:

*Lactobacillus*, *Pediococcus*, Acetic acid bacteria, Enterobacteriaceae, Bottom-fermented yeast, Top-fermented yeast and wild yeast.

### 2. Test principle

The TaqMan<sup>®</sup> real-time PCR is based on hot-start PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ROX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition Control (ROX/DQ) is amplified together in one tube with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **Tubes 4, 5, 7 and 8 contain lyticase.**

### 3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (24 samples):

1 x Premix	white cap
12 x Dye Strips (freeze-dried, tubes 4, 5, 7 und 8 contain lyticase)	tube strips
1 x ddH <sub>2</sub> O (for PCR negative controls)	colourless cap

#### **4. Storage conditions**

Do **not** freeze the lyophilized Dye Strips.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

#### **5. Materials required, but not provided**

##### **5.1. Instruments**

Real-time PCR machine for low profile tubes

Centrifuge for tube strips

Pipettes

“Vortex”

##### **5.2. Reagents and plastic ware**

sterile ddH<sub>2</sub>O

sterile filter tips

**optional:** Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

## 6. PCR

### 6.1. PCR-Setup

Before every use thoroughly mix Premix and centrifuge briefly.

PCR-reaction setup:

<b>PCR-Components</b>	<b>amount (µL)</b>
Premix	15.0
Sample-DNA	5.0*
<b>Total volume</b>	<b>20.0</b>

\* if using the Simplex Easy DNA- or QuickGEN Yeast Sample Preparation Kit:  
add 2.5 µL DNA and 2.5 µL sterile ddH<sub>2</sub>O

Pipetting scheme:

Tube 1*	Premix 15.0 µL	ddH <sub>2</sub> O 5.0 µL, <b>Attention:</b> pipette no sample-DNA (negative control)
Tube 2*	Premix 15.0 µL	ddH <sub>2</sub> O 5.0 µL, <b>Attention:</b> pipette no sample-DNA (positive control)
Tube 3*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 4*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 5*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 6*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL
Tube 7*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL
Tube 8*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL

\* please notice the orientation of the tubes: they have to be opened to the right side

## 6.2 PCR-Program

**Note:** For the use of UNG the thermal cycler program has to be changed according to manufacturers` instructions.

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35 x
Cycling Annealing/ Elongation	20 sec	62 °C	

\*for preenriched samples program 30 cycles

### For MyGo Pro User:

#### Software 3.3.2

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35 x
Cycling Annealing/ Elongation	20 sec	62 °C	

\*for preenriched samples program 30 cycles

<b>Advanced Settings</b>	Integration time (s)	0.5
	Acquisitions per Cycle	10

#### Software 3.4.8

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35 x
Cycling Annealing/ Elongation	30 sec	62 °C	

\*for preenriched samples program 30 cycles

<b>Advanced Settings</b>	Integration time (s)	0.5
--------------------------	----------------------	-----

## 7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

The negative control (NTC, tube 1) has to be negative in the FAM-, and HEX-channel and positive in the ROX-channel (Inhibition Control). Is the Inhibition Control in the ROX-channel negative and/or the results in the FAM-/ HEX-channel positive, PCR has to be repeated.

The Inhibition Controls in the samples (tubes 4/7) have to be positive in the ROX-channel. Is the Inhibition Control in the ROX-channel negative inhibitors or high amount of DNA are in the sample reaction.

The Positive Control (tube 2) has to be positive in the FAM-channel, if it is negative PCR, has to be repeated.

### \* Wild yeast 1 detects:

<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>	<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Debaromyces hansenii</i>
<i>Hanseniaspora guilliermondii</i>	<i>Hanseniaspora osmophila</i>	<i>Hanseniaspora uvarum</i>	<i>Issotchenkia orientalis</i>	<i>Kazachstania exigua</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrina</i>	<i>Wickerhamomyces anomalus</i>	<i>Pichia fermentans</i>	<i>Pichia membranaefaciens</i>
<i>S. cerevisiae var. diastaticus</i>	<i>Saccharomyces ludwigii</i>	<i>Torulaspota delbrückii</i>		

### \* Wild yeast 2 detects:

<i>Candida glabrata</i>	<i>Candida albicans</i>	<i>Candida kefyr</i>	<i>Candida intermedia</i>	<i>Candida parapsilosis</i>
<i>Candida sake</i>	<i>Candida tropicalis</i>	<i>Naumovozyma dairenensis</i>	<i>Pichia guilliermondii</i>	<i>Zygosaccharomyces bailii</i>
<i>Zygosaccharomyces rouxii</i>				

## PCR-Analysis:

please notice the orientation of the tubes: they have to be opened to the right side

Tube	Sample	FAM	HEX	ROX
1	NTC	-	-	Inhibition Control
2	PTC	Positive Control	-	-
3	1	Enterobacteriaceae	<i>Lactobacillus / Pediococcus</i>	<i>Pediococcus</i>
4	1	Wild yeast 1*	Bottom fermented yeast	Inhibition Control
5	1	Wild yeast 2*	Top fermented yeast	Acetic acid bacteria
6	2	Enterobacteriaceae	<i>Lactobacillus / Pediococcus</i>	<i>Pediococcus</i>
7	2	Wild yeast 1*	Bottom fermented yeast	Inhibition Control
8	2	Wild yeast 2*	Top fermented yeast	Acetic acid bacteria

NTC: negative Control    PTC: positive Control

### Example for interpretation:

If sample 1 shows a positive FAM-signal in tube 4, wild yeast 1 is identified. If the sample contains mixtures, you have positive signals in different channels.

#### Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.