# **RIDA® QUICK Aflatoxin RQS**

**Art. No. R5208**

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# RIDA<sup>®</sup>QUICK Aflatoxin RQS

## Brief information

RIDA<sup>®</sup>QUICK Aflatoxin RQS (Art. No. R5208) is an immunochromatographic test for quantitative analysis of aflatoxin in corn (maize).

The evaluation is carried out with the RIDA<sup>®</sup>SMART APP software (Art. No. ZRSAM1000) or the RIDA<sup>®</sup>QUICK SCAN reader (Art No. ZG5005).

Preparation of samples: homogenize, extract and dilute

Time requirement: Sample preparation (for 10 samples) approx. 10 min  
Test procedure (incubation period)..... 3 min

Limit of detection: < 2 µg/kg (ppb)

Measurement range: 2 - 75 µg/kg and 50 - 300 µg/kg

Specificity: The RIDA<sup>®</sup>QUICK Aflatoxin RQS test detects aflatoxin in corn.

The specificity of the RIDA<sup>®</sup>QUICK Aflatoxin RQS test was established by determining the cross-reactivity to the corresponding substances in the buffer system. The specificity in the sample can differ from the values determined in the buffer system due to the matrix effects. Before cross-reactive substances are analyzed, their detection limit and recovery rate in the relevant matrix must be determined by the user. The test cannot distinguish between analytes and cross-reactive substances.

## Related products

Trilogy<sup>®</sup> quality control materials (naturally contaminated)

Trilogy<sup>®</sup> and R-Biopharm Rhône Mycotoxin Standards

## 1. Intended use

RIDA®QUICK Aflatoxin RQS is a quantitative immunochromatographic test in a test strip format to determine aflatoxin in corn.

## 2. General information

Aflatoxins are secondary metabolic products of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*. These types of fungi most frequently occur in humid, tropical regions. Agricultural products become contaminated in the cultivating countries. Aflatoxins are among the strongest, naturally occurring carcinogenic substances.

## 3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-aflatoxin antibody detects aflatoxin in the sample. During incubation of the test strip, a band pattern forms that is used to determine the concentration of aflatoxin.

The appearance and the intensity of the test line depends on the aflatoxin concentration of the sample. It increases as the aflatoxin concentration increases.

It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the aflatoxin concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software or the RIDA®QUICK SCAN reader.

## 4. Reagents provided

The reagents in the kit are sufficient for 20 determinations. Every test kit contains:

Component	Lid color	Format	Contents
20 x Test strip	-	Ready to use	One analysis per test strip (individually packaged)
1 x Mobile solvent	White	Ready to use	4.5 ml
1 x Dilution buffer	Transparent	Ready to use	25 ml
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for use with RIDA®SMART APP

## **5. Necessary additional reagents - required accessories**

### **5.1. Equipment**

- Grinder
- Weighing scale
- Graduated cylinder
- Pipettes (e.g. R-Biopharm FP 100, Art. No. Z0007)
- Shaker (Rock-It 360 from Trilogy<sup>®</sup> Analytical Laboratory Inc. or similar)
- Centrifuge ( $\geq 2000$  g)
- Latest version of the RIDA<sup>®</sup>SMART APP software (Art. Nr. ZRSAM1000), activated and installed on a compatible smartphone. Alternatively, RIDA<sup>®</sup>QUICK SCAN (Art. Nr. ZG5005) with software version V1.0.1.7 2.0 or higher.

### **5.2. Reagents**

- 70 % Methanol (e.g., 70 ml methanol (100%) / 30 ml dist. water)
- 50 % Ethanol (e.g., 50 ml ethanol (denatured) / 50 ml dist. water)

## **6. Warnings and precautions for the users**

This test should be carried out only by trained personnel. The instructions for performing the test must be strictly complied with.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. This applies particularly if the test strip packaging has been opened.

This kit may contain harmful substances. Please refer to the component safety information in the material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.de](http://www.r-biopharm.de).

## **7. Storage instructions**

Store the test at 2 - 8 °C. Do not freeze components of the test kit.

After the expiration date (see expiration on outer label of test kit), the quality guarantee is no longer valid.

The exchange of individual reagents between kits with different lot numbers is not permitted.

## 8. Sample preparation

**All test kit components used must be bought to room temperature (20 - 25 °C / 68 - 77 °F) prior to use in the test; the samples must be prepared at room temperature. Do not allow methanol or ethanol to evaporate.**

The sample should be stored in a cool place, protected from light.

Before extraction, grind and carefully mix a representative sample (a sample taken in accordance with official sampling regulations).

### 8.1. Measurement range 2 - 75 µg/kg

- Weigh 10 g of ground and homogenized sample and add 20 ml methanol (70 %) or 20 ml ethanol (50 %)
- Seal the vials, shake the sample for 5 minutes (on shaker)
- Centrifuge (1 min, 2000 g) to obtain a particle-free supernatant
- Add 100 µl of particle-free supernatant to 200 µl of mobile solvent and mix carefully (= Dilution I)
- Use 100 µl of Dilution I in the test

#### Note:

**The weight of the sample can be increased if necessary. In this case, the methanol or ethanol volume must be appropriately adjusted, e.g., 20 g sample in 40 ml methanol (70 %) or 40 ml ethanol (50 %).**

### 8.2. Measurement range 50 - 300 µg/kg

For an extended measurement range of 50 - 300 µg/kg, further dilute Dilution I (see sample preparation 8.1.).

- Dilute 100 µl Dilution I with 900 µl dilution buffer and mix carefully (= Dilution II)
- Use 100 µl of Dilution II in the test

## 9. Test procedure

All components used must be bought to room temperature (20 - 25 °C/ 68 - 77 °F) prior to use in the test; the procedure has to be performed at room temperature.

- Apply 100 µl of prepared sample solution (Dilution I or Dilution II, see 8. Sample preparation) to the application area of the test strip
- After exactly 3 min of incubation, evaluate the test strip with the RIDA®SMART APP or the RIDA®QUICK SCAN

## 10. Evaluation

The left line in the reaction area is the test line (see Fig. 1). Its appearance and its intensity depends on the aflatoxin concentration of the sample.

The right-hand line in the reaction area is a control line and must appear after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. Please notify us if the right line is still missing after repeated tests.

**Please note:** If the sample is highly contaminated with aflatoxin, the control line will be only faintly visible. In this case the test result might be reported as „invalid result“ (RIDA®SMART APP) or „invalid“ (RIDA®QUICK SCAN reader). In this instance, dilute the sample as described in 8.2 and repeat the test with a new test strip.

Test strips must always be evaluated after 3 minutes in order to achieve comparable results.


Fig. 1: RIDA®QUICK Aflatoxin RQS test strip

## 10.1. Evaluation with RIDA®SMART APP

First, carefully read the user guide for the RIDA®SMART APP. The evaluation must be performed using these specifications.

Brief instructions on how to use the software application are available online at [www.r-biopharm.de](http://www.r-biopharm.de).

The RIDA®SMART APP cover enclosed in the test kit is required for test strip evaluation. The RIDA®SMART APP cover contains lot-specific information and must be used only with the intended test-kit lot. More information can be found on the analysis certificate enclosed with the test kit.

Please make sure the correct RIDA®SMART APP evaluation cover is present prior to starting the test procedure and test strip evaluation.

RIDA®SMART APP software applications depending on the extraction solution:

Extraction method	RIDA®SMART APP application
70 % Methanol	Methanol 2 - 75 ppb
	Methanol 50 - 300 ppb
50 % Ethanol	Ethanol 2 - 75 ppb
	Ethanol 50 - 300 ppb

## 10.2. Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed.

**Please note:** The Quality Assurance Certificate contains parameters for several standard curves. Depending on the extraction method used, select the corresponding parameters.

To enter the parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the Quality Assurance Certificate. To enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN methods:

Measuring range	RIDA®QUICK SCAN method
2 - 75 µg/kg	Afla ECO Low
50 - 300 µg/kg	Afla ECO High

## 11. Sensitivity

Using the RIDA®QUICK Aflatoxin RQS test in combination with the RIDA®SMART APP or the RIDA®QUICK SCAN, it is possible to detect aflatoxin contaminations between 2 - 300 µg/kg.

**For further product information please contact your local distributor or R-Biopharm AG ([sales@r-biopharm.de](mailto:sales@r-biopharm.de)).**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### **R-Biopharm AG**

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-Mail: [info@r-biopharm.de](mailto:info@r-biopharm.de)

[www.r-biopharm.com](http://www.r-biopharm.com)

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are registered trademarks of R-Biopharm AG