# **RIDA®QUICK Aflatoxin RQS ECO**

**Art. No. R5209**

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# RIDA®QUICK Aflatoxin RQS ECO

## Brief information

RIDA®QUICK Aflatoxin RQS ECO (Art. No. R5209) is an immunochromatographic test for the quantitative analysis of aflatoxin in corn (maize). The test uses an aqueous extraction method.

The evaluation is carried out with the RIDA®SMART APP software (Art. No. ZRSAM1000) or the RIDA®QUICK SCAN reader (RQS, Art. No. ZG5005).

Preparation of samples: homogenize and extract

Time requirement:            Sample preparation (for 10 samples) .... approx. 10 min  
   Test procedure (incubation period) ..... 5 min

Limit of detection:            < 2 µg/kg (ppb)

Measurement range:        2 - 75 µg/kg and 50 - 300 µg/kg

Specificity:                    The RIDA®QUICK Aflatoxin RQS ECO test detects aflatoxin in corn (maize).

The specificity of the RIDA®QUICK Aflatoxin RQS ECO test was established by determining the cross-reactivity to the corresponding substances in the buffer system. The specificity in the sample can differ from the values determined in the buffer system due to the matrix effects. Before cross-reactive substances are analyzed, their detection limit and recovery rate in the relevant matrix must be determined by the user. The test cannot distinguish between analytes and cross-reactive substances.

## Related products

Trilogy® quality control materials (naturally contaminated)

Trilogy® and R-Biopharm Rhône Mycotoxin Standards

## 1. Intended use

RIDA®QUICK Aflatoxin RQS ECO is a quantitative immunochromatographic test in a test strip format to determine aflatoxin in corn. The test uses an aqueous extraction method.

## 2. General information

Aflatoxins are secondary metabolic products of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*. These types of fungi most frequently occur in humid, tropical regions. Agricultural products become contaminated in the cultivating countries. Aflatoxins are among the strongest, naturally occurring carcinogenic substances.

## 3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-aflatoxin antibody detects aflatoxin in the sample. During incubation of the test strip, a band pattern forms that is used to determine the concentration of aflatoxin.

The intensity of the test line depends on the aflatoxin concentration of the sample. It increases as the aflatoxin concentration increases.

It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the aflatoxin concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software or the RIDA®QUICK SCAN reader.

## 4. Reagents provided

Each kit contains sufficient materials for 20 determinations. Each test kit contains:

Component	Cap color	Format	Content
20 x Test strip	-	Ready to use	One analysis per test strip (individually packaged)
1 x Extraction solution	Transparent	Concentrate 10x	45 ml
1 x Dilution buffer	Transparent	Ready to use	20 ml
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for use with RIDA®SMART APP

## **5. Necessary additional reagents - required accessories**

### **5.1. Equipment**

- Grinder
- Weighing scale
- Graduated cylinder
- 1 L bottle
- Pipettes
- Shaker (Rock-It 360 from Trilogy<sup>®</sup> Analytical Laboratory Inc. or similar)
- Centrifuge ( $\geq 2000$  g)
- Latest version of the RIDA<sup>®</sup>SMART APP software (Art. Nr. ZRSAM1000) activated and installed on a compatible smartphone. Alternatively RIDA<sup>®</sup>QUICK SCAN (Art. Nr. ZG5005) with software version V1.0.1.7 2.0 or higher.

### **5.2. Reagents**

- Distilled or deionized water

## **6. Warnings and precautions for the users**

This test should be carried out only by trained personnel. The instructions for performing the test must be strictly complied with.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. This applies particularly if the test strip packaging has been opened.

This kit may contain harmful substances. Please refer to the component safety information in the material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.de](http://www.r-biopharm.de).

## 7. Storage instructions

Store the test at 2 - 8 °C. Do not freeze components of the test kit.

After the expiration date (see expiration on outer label of test kit), the quality guarantee is no longer valid.

The exchange of individual reagents between kits with different lot numbers is not permitted.

## 8. Sample preparation

To obtain the ready-to-use extraction buffer, transfer 1 part of the extraction solution (concentrate 10x) into 9 parts of distilled or deionized water and homogenize for 5 min. The ready-to-use extraction buffer can be stored for approx. 8 - 10 weeks at 2 - 8 °C. If not all test strips are used within 8 - 10 weeks it is recommended to prepare small volumes of the ready-to-use extraction buffer.

**All test kit components used must be brought to room temperature (20 - 25 °C / 68 - 77 °F) prior to use in the test; the samples must be prepared at room temperature.**

The sample should be stored in a cool place, protected from light.

Before extraction, grind and carefully mix a representative sample (a sample taken in accordance with official sampling regulations).

### 8.1. Measurement range 2 - 75 µg/kg

- Weigh 10 g of ground and homogenized sample and add 20 ml of ready-to-use extraction buffer
- Shake for 5 min (on shaker)
- Centrifuge (3 min at 2000 g) until a particle-free solution can be obtained
- Use 100 µl of the particle-free solution in the test

## 8.2. Measurement range 50 - 300 µg/kg

For an extended measurement range of 50 - 300 µg/kg, further dilute the particle-free solution obtained from sample preparation 8.1.

- Dilute 100 µl of the particle-free solution with 700 µl dilution buffer and mix carefully (= Dilution I)
- Use 100 µl of Dilution I in the test

## 9. Test procedure

**All components used must be brought to room temperature (20 - 25 °C/ 68 – 77 °F) prior to use in the test; the samples must be prepared at room temperature.**

- Apply 100 µl of prepared sample solution (particle-free solution or Dilution I, see 8. Sample preparation) to the application area of the test strip
- After exactly 5 min of incubation, evaluate the test strip with the RIDA®SMART APP or the RIDA®QUICK SCAN

## 10. Evaluation

The left line in the reaction area is the test line. Its intensity depends on the aflatoxin concentration of the sample.

The right-hand line in the reaction area is a control line and must appear after each test procedure (see Fig. 1). If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. Please notify us if the control line is still missing after repeated tests.

**Please note:** If the sample was prepared according to 8.1. and is highly contaminated with aflatoxin, the control line will be only faintly visible. In this case the test result might be reported as „invalid result“ (RIDA®SMART APP) or „invalid“ (RIDA®QUICK SCAN reader). In this instance, dilute the sample as described in 8.2 and repeat the test with a new test strip.

Test strips must always be evaluated after 5 minutes in order to achieve comparable results.


Fig.1: RIDA®QUICK Aflatoxin RQS ECO test strip

### 10.1. Evaluation with RIDA®SMART APP

First, carefully read the user guide for the RIDA®SMART APP. The evaluation must be performed using these specifications.

Brief instructions on how to use the software application are available online at [www.r-biopharm.de](http://www.r-biopharm.de).

The RIDA®SMART APP cover enclosed in the test kit is required for test strip evaluation. The RIDA®SMART APP cover contains lot-specific information and must be used only with the intended test-kit lot. More information can be found on the Quality Assurance Certificate enclosed with the test kit.

Please make sure the correct RIDA®SMART APP evaluation cover is present prior to starting the test procedure and test strip evaluation.

RIDA®SMART APP software applications:

Detection range	RIDA®SMART APP application
2 - 75 µg/kg	2 - 75 ppb
50 - 300 µg/kg	50 - 300 ppb

### 10.2. Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed.

**Please note:** The Quality Assurance Certificate contains several 2D barcodes. Depending on the measurement range used, select the corresponding barcode for evaluation using the RIDA®QUICK SCAN.

To enter the barcode parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the Quality Assurance Certificate; to enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN methods:

Measuring range	RIDA®QUICK SCAN method
2 - 75 µg/kg	Afla ECO Low
50 - 300 µg/kg	Afla ECO High

## 11. Sensitivity

Using the RIDA®QUICK Aflatoxin RQS ECO test in combination with the RIDA®SMART APP software or the RIDA®QUICK SCAN reader, it is possible to detect aflatoxin contamination between  $\geq 2 - 300 \mu\text{g/kg}$ .

**For further product information and applications please contact your local distributor or R-Biopharm AG ([sales@r-biopharm.de](mailto:sales@r-biopharm.de)).**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-Mail: [info@r-biopharm.de](mailto:info@r-biopharm.de)

[www.r-biopharm.com](http://www.r-biopharm.com)

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are registered trademarks of R-Biopharm AG