
CONGEN

# SureFast<sup>®</sup> Salmonella PLUS

Art. No. F5111  
100 rxn

## User Manual

November 2019

 **Inhalt**

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 **Content**

1	General Information .....	8
1.1	Description .....	8
1.2	Limit of Detection .....	8
1.3	DNA-preparation .....	8
1.4	Kit components and storage .....	8
1.5	Additionally required equipment and materials .....	8
1.6	Setup .....	9
1.7	Detection channel Set-up .....	9
2	Qualitative Analysis .....	10
2.1	Protocol .....	10
2.1.1	Preparation of the master-mix .....	10
2.1.2	Preparation of the real-time PCR-mix .....	10
2.2	Interpretation of results .....	11
3	Further Information .....	11
3.1	Product Information .....	11
3.2	Technical Support .....	11

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## 1 General Information

### 1.1 Description

The SureFast® Salmonella PLUS is a real-time PCR for the direct, qualitative detection of a specific DNA sequence of *Salmonella* spp.. The SureFast® Salmonella PLUS real-time PCR kit has been validated and certified in combination with the SureFast® PREP Salmonella kit by AOAC (licence number: 041103).

Each reaction contains an internal amplification control (IAC).

The real-time PCR assay can be used with real-time PCR instruments, equipped for detection of two fluorescence emissions in channels FAM and VIC/HEX at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Bio-Rad CFX96, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Roche LightCycler® 2.0, Applied Biosystems 7500, Cepheid SmartCycler®, BMS MIC, LTF MyGo Pro and Qiagen Rotor-Gene Q.

### 1.2 Limit of Detection

The SureFast® Salmonella PLUS real-time PCR has a limit of detection of  $\leq 5$  DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Salmonella or SureFast® Speed PREP is recommended. To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at –20°C and protected from light.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit  
(e.g. SureFast® PREP Salmonella Art. No. F1007 / Speed PREP Art. No. F1054)
- Real-time PCR instrument with two detection channels (510 nm and 580 nm)
- Real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- Pipettes with filter tips
- Unpowdered disposable gloves
- Vortex mixer
- Micro centrifuge with a rotor for the reaction tubes

**1.6 Setup**

	<b>Blockcycler &amp; BMS MIC</b>	<b>Rotorcycler &amp; LightCycler® 480 II &amp; LTF MyGo Pro</b>
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.7 Detection channel set-up**

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	<i>Salmonella</i> spp.	FAM	+	
	IAC	HEX	+	
<b>Applied Biosystems 7500</b>	<i>Salmonella</i> spp.	FAM	None	Use the passive reference dye ROX setting "none".
	IAC	VIC	None	
<b>Bio-Rad CFX96</b>	<i>Salmonella</i> spp.	FAM	+	
	IAC	VIC/HEX	+	
<b>BMS MIC</b>	<i>Salmonella</i> spp.	green	+	
	IAC	yellow	+	
<b>LTF MyGo Pro</b>	<i>Salmonella</i> spp.	FAM	+	
	IAC	VIC	+	
<b>Qiagen Rotor-Gene Q</b>	<i>Salmonella</i> spp.	green	+	
	IAC	yellow	+	
<b>Roche LightCycler® 2.0</b>	<i>Salmonella</i> spp.	530	None	The SureCC Color Compensation Kit II (Art. No. F4010) is required. Set the Seek Temperature to 58°C.
	IAC	560	None	
<b>Roche LightCycler® 480 II</b>	<i>Salmonella</i> spp.	465-510	+	
	IAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	<i>Salmonella</i> spp.	465-510	+	
	IAC	540-580	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

#### Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µl of Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/ wells or capillaries into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

*Salmonella* spp. DNA is detected in the FAM-channel. In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for *Salmonella* spp., if the sample DNA shows amplification in the FAM-channel.

A sample is stated **negative** for *Salmonella* spp., if the sample DNA shows no amplification in the FAM-channel and the internal control (VIC/HEX-channel) of the sample is **positive**.

If the sample DNA and the internal amplification control are **negative** the sample contains PCR-inhibiting substances or the sample preparation was not successful. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved.

## 3 Further Information

### 3.1 Product Information

Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage. Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads)

- Validation Report upon request

### 3.2 Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).