CONGEN

# SureFast<sup>®</sup> Salmonella ONE

Art. No. F5211

100 rxn

## User Manual

MICROVAL 

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## 1 General Information

### 1.1 Description

The SureFast® Salmonella ONE kit can be applied for the fast and simple isolation and detection of *Salmonella* spp. in different types of food and is intended to be used by all kinds of food testing laboratories.

The method consists of the following steps: 1) cultural enrichment, 2) DNA extraction, 3) specific real-time PCR detection and 4) data analysis and interpretation. SureFast® Salmonella ONE has been validated and certified by MicroVal (License-No. 2014LR43) and AOAC (License-No. 081803).

Each reaction contains an internal amplification control. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of two fluorescence emissions in the FAM and VIC/HEX channel at the same time.

Internal validation was performed on Agilent Mx3005P, BioRad CFX 96, Roche LightCycler® 480 II, Roche cobas z 480 Analyzer, Roche LightCycler® 2.0<sup>2</sup>, Applied Biosystems® 7500, Qiagen RotorGene Q, Cepheid SmartCycler®, Bio Molecular Systems MIC and LTF MyGo Pro.

External validation by MicroVal and AOAC was performed on Roche LightCycler® 480 II, Applied Biosystems® 7500, BioRad CFX 96, Bio Molecular Systems MIC and Agilent AriaMx.

### 1.2 Limit of Detection

The SureFast® Salmonella ONE real-time PCR has a limit of detection of  $\leq 5$  DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
L	Lysis Buffer	2 x 25 ml	Clear

Store all reagents at -20 °C and protected from light. The Taq Polymerase can be stored at +2 to +8 °C for multiple use on the same day. The Lysis Buffer should be stored at -20 to +8 °C. Reagents should be used until the expiry date indicated on the label.

### 1.4 Additionally required equipment, materials and media

- Real-time PCR instrument, equipped with two detection channels (FAM and VIC/HEX channel)
- Real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- Pipettes with filter tips
- 2 mL tubes
- Disposable gloves
- Vortexmixer
- Centrifuge with a rotor for the reaction tubes
- Thermomixer/heating block (up to 95 °C)
- Non-selective enrichment broth (e. g. buffered peptone water)

<sup>2</sup> note: For the use of the Roche LightCycler® 2.0 a Color Compensation is necessary. The SureCC Color Compensation Kit II (Art. No. F4010) must be used for the color compensation of such devices.

### 1.5 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

## 2 Protocol

### 2.1 Enrichment

Sample preparation of food and feed matrices should apply to ISO 6887 and ISO 6579 standards. For enrichment, a volume of 25 g or 25 ml sample material is homogenized with 225 ml of a non-selective enrichment broth (e. g. buffered peptone water). If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio. Enrichment is performed for 18 h ± 2 h at 37 °C ± 1°C. Alternatively, other suitable, validated enrichment procedures can be used.

To assess the process of bacterial growth, it is recommended to compare samples at the beginning and at the end of the culturing (bacterial growth at Cp difference >3).

### 2.2 DNA preparation

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 500 µl of the upper third from the enrichment culture into a 2.0 mL tube (not provided with the kit).

Add 500 µl of Lysis Buffer (**Code L**) to the 500 µl of the enrichment culture. Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

**Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction.**

The lysate is ready-to-use for PCR. If the lysate is not used immediately in the PCR or is intended for a longer storage, transfer 100 µl of the lysate in a new tube (not provided with the kit) and store at -20 °C.

### 2.3 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control. The master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix each master-mix well and centrifuge shortly before use.

**Note:** The Taq Polymerase may be in a frozen and/or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

**2.4 Setup**

	<b>Blockcycler/MIC SmartCycler®</b>	<b>Rotorcycler/LightCycler® 480 II</b>
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup (exemplary)	Detection: End of Extension Phase  Detection System: <i>Salmonella</i> spp. Various devices <b>FAM</b> -Channel, Quencher: BHQ Rotor-Gene Q Green LightCycler 2.0 530 nm - none LightCycler 480 II 465 nm - 510 nm  Internal Amplification Control: Various devices <b>VIC/HEX</b> -Channel, Quencher: BHQ Rotor-Gene Q Yellow LightCycler 2.0 560 nm - none LightCycler 480 II 533 nm - 580 nm  Passive Reference: None	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: <a href="http://www.congen.de/en/company/downloads">http://www.congen.de/en/company/downloads</a>		

**2.5 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control consist of the master-mix).
- Pipette 5 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µl of the Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/wells or capillaries into the real-time PCR instrument and start the run according to the setup.

**2.6 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive**, if the sample DNA shows an amplification in the detection system (FAM).

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system and the internal amplification control (inhibition control) of the sample is **positive** (VIC/HEX) with a shift in Cp-Value ≤ 2 compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and

November 2018

purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Salmonella* PCR assay.

**Note:** During the MicroVal study increased inhibition levels were observed for fresh herbs, cocoa, spices and aromatics and egg products. In these cases, please proceed as recommended above. If you require further assistance please contact [info@congen.de](mailto:info@congen.de).

### 2.7 Cultural Confirmation

Positive PCR results should be confirmed by cultural detection methods (e.g. ISO 6579).

## 3 Product Information

- Validation Report
- Certificate for MicroVal (License-No. 2014LR43) and AOAC (License-No. 081803)

## 4 Validation

SureFast® Salmonella ONE was validated per the AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces and the ISO 16140-2:2016 [Microbiology of the food chain -- Method validation -- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method]. The reference method for the matrix study was ISO 6579:2002 (Microbiology of food and animal feeding stuffs -- Horizontal method for the detection of *Salmonella* spp.).

The kit is an extension of the SureFast® Salmonella PLUS method and has been AOAC RI validated for the following additional food matrices: raw ground beef, frozen poultry meat, raw milk, fresh spinach, pasteurized liquid whole egg, pet food pellets.

The MicroVal method comparison study according to ISO 16140-2:2016 was carried out using 25 g portions of sample material. Sample preparations were done according to ISO 6887 parts. In the sensitivity study a total of 391 samples from six different food categories were investigated. SureFast® Salmonella ONE has been MicroVal validated for the following food categories: poultry meat, meat products, dairy products, vegetables (excluding sprouts), egg products, feed.

## 5 Technical Support

For further questions please contact:

CONGEN Biotechnologie GmbH  
Robert-Roessle-Str. 10  
13125 Berlin, Germany

Phone: +49 (0)30 9489 3500  
Fax: +49 (0)30 9489 3510  
E-Mail: [info@congen.de](mailto:info@congen.de)