**User Manual**


# **SureFast<sup>®</sup> PREP Bacteria**

## **Art. No. F1021**

**100 extractions**

Efficient DNA preparation of bacteria

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**Table of contents**

Description ..... 7  
Principle ..... 7  
Kit components (per box) and storage ..... 7  
Additionally required equipment and materials ..... 7  
Preparations ..... 8  
    General ..... 8  
    Before each preparation ..... 8  
Extraction protocol ..... 8  
Product Information ..... 9  
Technical Support ..... 9  
Safety information.....10

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1x Lysis Buffer (20 ml)	<b>(Code L)</b>	1x Receiver Tubes 2,0 ml, klar (50 x)	<b>(Code R)</b>
1x Binding Buffer (10 ml)	<b>(Code B)</b>	1x Receiver Tubes 1,5 mL, klar (50 x)	<b>(Code T)</b>
1x Wash Buffer (60 ml)*	<b>(Code W)</b>	1x Spin Filter, klar (50 x)	<b>(Code S)</b>
1x Elution Buffer (10 ml)	<b>(Code E)</b>		

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**Description**

The kit is intended to be used for the isolation of bacteria - DNA from food (enrichments, avulsions and swabs).

In samples with a very low microbial activity (e.g heat pasteurized egg white powder) a contamination with not viable or not culturable bacteria cells might be observed. In case of other conditions the samples show a fast metabolization of these cells. To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing

**Principle**

1. Preparation of the basic material
2. Lysis at 99°C
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

**Kit components (per box) and storage**

1x Lysis Buffer (20 ml)	<b>(Code L)</b>	1x Receiver Tubes 2.0 ml, clear (50x)	<b>(Code R)</b>
1x Binding Buffer (10 ml)	<b>(Code B)</b>	1x Receiver Tubes 1.5 ml, clear (50x)	<b>(Code T)</b>
1x Wash Buffer (60 ml)*	<b>(Code W)</b>	1x Spin Filter, clear (50x)	<b>(Code S)</b>
1x Elution Buffer (10 ml)	<b>(Code E)</b>		

\* After adding ethanol (purity  $\geq$  96 %; not supplied with the kit).

All reagents of the kit should be stored dry and at room temperature (14-25°C).

**Additionally required equipment and materials**

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 99°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol (purity  $\geq$  96 %) for preparation of Pre-Wash Buffer and Wash Buffer
- disposal bags or waste bin
- autoclave
- biological safety cabinet

**Preparations****General**

- Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

**Before each preparation**

- Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer (**Code E**) under calculation of a reserve volume into a reaction tube (not supplied with the kit) and equilibrate to 60°C (The Elution Buffer is necessary for step 7).

**Extraction protocol****1. Preparation of the basic material**

Transfer 1.0 ml of the enrichment or avulsions under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).

Centrifuge for 5 min at 12.000 rpm.

Discard the fluid and inactivate per sterilization.

For swabs, place them into a 1.5 ml reaction tube (not provided with the kit) and cut the shaft, so that you can close the cap of the extraction tube.

**2. Lysis of the basic material**

Add 400 µl of Lysis Buffer (**Code L**), mix it briefly and close the reaction tube with caps.

Incubate under continuously shaking for 10 minutes at 99°C.

**3. Setting of optimal binding conditions**

Centrifuge the sample lysate for 1 min at 12.000 rpm.

Transfer ca. 300 µl of the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit).

**Note:** Depending on the nature of the sample a smaller amount of liquid may be obtained after centrifugation.

In this case take the actual volume for the following steps or vary the volume of the Lysis Buffer before lysis.

**4. Binding of the nucleic acids on a Spin Filter**

Add 200 µl Binding Buffer (**Code B**) to the filtrate and mix the sample.

Place a Spin Filter (**Code S**) into a 2.0 ml Receiver Tube (**Code R**).

Transfer the solution directly onto the filter and incubate at room temperature for 1 min.

Centrifuge for 1 min at 12.000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

**5. Purification of the bound nucleic acids**

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12.000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add again 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12.000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

**6. Drying of the Spin Filter**

Remove the residual Ethanol by centrifugation for 2 min at 12.000 rpm.

**7. Elution of nucleic acids from the Spin Filter**

Place the Spin Filter into a 1.5 ml Receiver Tube (**Code T**) and add 100 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min at room temperature and centrifuge for 1 min at 10.000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.

**Product Information**

- Material Safety Data Sheet
- Flow chart  
(Download: <http://www.congen.de/en/company/download>)

**Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).

## Safety information

### Binding Buffer

Danger

### Lysis Buffer

Warning

H225 H319 H336  
P210 P233 P305+P351+P338

H319 H412  
P280 P305+P351+P338

- H225:** Highly flammable liquid and vapour.  
**H319:** Causes serious eye irritation.  
**H336:** May cause drowsiness or dizziness.  
**H412:** Harmful to aquatic life with long lasting effects.  
**P210:** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
**P233:** Keep container tightly closed.  
**P280:** Wear protective gloves/protective clothing/eye protection/face protection.  
**P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information we offer a Material Safety Data Sheet. Please contact your distributor or send an e-mail to [info@congen.de](mailto:info@congen.de).