Compact Dry YMR medium for yeast and mould / Medium für Hefe und Schimmelpilze / milieu pour levure et moisissures / medio para hongos y moho / medio per fermenti e muffe / meio para levedura e fungos

100 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9801

40 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9802

English

Deutsch

Français

Compact Dry YMR is a ready to use, chromogenic plate for detection of yeasts and moulds

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Sample pretreatment

Viable count in water or liquid foodstuff

Drop 1 ml of specimen (dilute if necessary) on the middle of the Compact Dry plate.

Viable count in solid foodstuff

Add buffer solution (e.g. 0.9% saline solution, PBS, buffered peptone water, etc.) to the sample and homogenize by stomacher®. Drop 1 ml of specimen (dilute if necessary) on the middle of the dry sheet of the Compact Dry plate.

Viable count in swab test specimen

Use a steril wet swab to wipe the surface which should be analyzed and put it back into the device with wiping solution. After shaking the device (approx. 10 sec) drop 1 ml of wiping solution (dilute if necessary) on the middle of the Compact Dry plate. Use the Compact Dry Swab (Art. No. ZCS1002953, content 40 pcs) as single application (1 ml) or Promedia ST-25 (Art. No. Z0302, content 10 pcs) as multiple application (10 ml) for this procedure.

Test instructions

1. Open the cap and drop 1 ml of specimen on the middle of the Compact Dry plate.
2. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
3. Put the cap again on the plate and write the information needed on the memorandum section.
4. Turn over the capped plate and put in the incubator.
5. After incubation count the number of colored colonies underneath the plate. White paper placed under the plate helps to count the colonies.

Incubation time 48 - 72 hours
Incubation temperature 25 - 30 °C

Interpretation of the results

Blue colored colonies indicate yeasts.
Molds form cottony colonies with characteristic colors.

Storage and shelf life

Keep at room temperature (+ 1 to + 30 °C).
Total shelf life 12 months after manufacturing.

Notes

- Some yeasts do not form blue colonies.
- Antibiotics in the medium inhibit the growth of bacteria.
- On Compact Dry plates up to 300 cfu per plate can be detected. High concentrations on plates (> 300 cfu) will cause the entire growth area to become green / blue. In this case dilute the sample.
- After use please follow the current disposal regulations.
- The growth area is 20 cm². The back of the plate has a grid carved of 1 cm x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates are produced at an ISO 9001 / ISO 13485 certified site.

- AOAC submitted
- MicroVal approval No. 2016LR61
- NordVal approval No. 050

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