
**CONGEN**

# SureFast® PREP Bacteria

Art. No. F1021  
100 extractions

## User Manual

Efficient DNA preparation of bacteria

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**Content**

1	General Information .....	8
1.1	Description .....	8
1.2	Kit components and storage .....	8
1.3	Additionally required equipment and materials .....	8
2	Protocol .....	9
2.1	Principle .....	9
2.2	Preparations .....	9
2.3	Extraction protocol .....	9
	1. Preparation of the basic material .....	9
	2. Lysis of the basic material .....	9
	3. Setting of optimal binding conditions .....	10
	4. Binding of the nucleic acids on a Spin Filter .....	10
	5. Purification of the bound nucleic acids .....	10
	6. Drying of the Spin Filter .....	10
	7. Elution of nucleic acids from the Spin Filter .....	10
3	Further Information .....	10
3.1	Product Information .....	10
3.2	Technical Support .....	10
4	Safety Information .....	11

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Kit Code	Reagenz/Material	Menge (je Box)
<b>L</b>	Lysis Buffer	1 x 25 ml
<b>B</b>	Binding Buffer*	1 x 4 ml
<b>W</b>	Wash Buffer**	1 x 18 ml
<b>E</b>	Elution Buffer	1 x 10 ml
<b>S</b>	Spin Filter	1 x 50 Filter
<b>R</b>	Receiver Tubes 2,0 ml	1 x 50 Tubes
<b>T</b>	Receiver Tubes 1,5 ml	1 x 50 Tubes

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## 1 General Information

### 1.1 Description

The kit is intended to be used for the isolation of bacteria - DNA from food (enrichments, avulsions and swabs).

The SureFast® PREP Bacteria kit has been validated and certified in combination with the SureFast® Salmonella PLUS real-time PCR kit (Art. No. F5111) by AOAC (licence number: 041103).

In samples with a very low microbial activity (e.g. heat pasteurized egg white powder) the detection of viable or not culturable bacteria cells is possible. In case of other conditions the samples show a fast metabolization of these cells.

To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing.

### 1.2 Kit components and storage

Kit Code	Reagent /Material	Amount (per box)
<b>L</b>	Lysis Buffer	1 x 25 ml
<b>B</b>	Binding Buffer*	1 x 4 ml
<b>W</b>	Wash Buffer**	1 x 18 ml
<b>E</b>	Elution Buffer	1 x 10 ml
<b>S</b>	Spin Filter	1 x 50 Filter
<b>R</b>	Receiver Tubes 2.0 ml	1 x 50 Tubes
<b>T</b>	Receiver Tubes 1.5 ml	1 x 50 Tubes

\* Add isopropanol (not supplied with the kit, see 1.3 Additionally required equipment and materials)

\*\* Add ethanol (not supplied with the kit, see 1.3 Additionally required equipment and materials)

All reagents of the kit should be stored dry and at room temperature (14-25°C).

### 1.3 Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 99°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol for preparation of Wash Buffer (puriss, purity ≥96%)
- isopropanol for preparation of Binding Buffer (e.g. Carl Roth 2-Propanol Rotipurán® >99.7%; Applichem 2-Propanol for microbiology; Sigma-Aldrich 2-Propanol)
- disposal bags or waste bin
- autoclave
- biological safety cabinet

## 2 Protocol

### 2.1 Principle

1. Preparation of the basic material
2. Lysis of the basic material
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

### 2.2 Preparations

#### General

Add 11 ml isopropanol to the Binding Buffer (**Code B**) and mix thoroughly.

Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

#### Before each preparation

Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer under calculation of a reserve volume into a reaction tube (not provided with the kit).

Equilibrate the Elution Buffer to 60°C. The Elution Buffer is necessary for step 7.

### 2.3 Extraction protocol

#### 1. Preparation of the basic material

Transfer 1.0 ml of the enrichment or avulsions under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).

Centrifuge for 5 min at 12,000 rpm.

Discard the fluid and inactivate per sterilization.

For swabs, place them into a 1.5 ml reaction tube (not provided with the kit) and cut the shaft, so that you can close the cap of the extraction tube.

#### 2. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**) to the reaction tube and mix it briefly and close the reaction tube with caps.

Incubate on a heating block under continuously shaking for 10 min at 99°C.

**3. Setting of optimal binding conditions**

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Transfer ca. 300 µl of the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit).

**Note:** Depending on the nature of the sample a smaller amount of liquid may be obtained after centrifugation. In this case take the actual volume for the following steps or increase the volume of the Lysis Buffer before lysis.

**4. Binding of the nucleic acids on a Spin Filter**

Add 200 µl Binding Buffer (**Code B**) to the supernatant and mix.

Place a Spin Filter (**Code S**) into a new 2.0 ml Receiver Tube (**Code R**). Transfer the solution onto the Spin Filter.

Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm.

After centrifugation discard the filtrate and place the Spin Filter back into the Receiver Tube.

**5. Purification of the bound nucleic acids**

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

**6. Drying of the Spin Filter**

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

**7. Elution of nucleic acids from the Spin Filter**

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 100 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min at room temperature.

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

## 3 Further Information

**3.1 Product Information**

- Validation data upon request
- Material Safety Data Sheet (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))

**3.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).

## 4 Safety Information

### Lysis Buffer

#### Warning

H319-H412-P264-P273  
P280-P305+P351+P338-  
P337+P313-P501

<b>H319:</b>	Causes serious eye irritation.
<b>H412:</b>	Harmful to aquatic life with long lasting effects.
<b>P264:</b>	Wash hands, forearms and face thoroughly after handling.
<b>P273:</b>	Avoid release to the environment.
<b>P280:</b>	Wear protective gloves/protective clothing/eye protection/face protection.
<b>P305+P351+P338:</b>	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
<b>P337+P313:</b>	If eye irritation persists: Get medical advice/attention.
<b>P501:</b>	Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

For further information we offer a Material Safety Data Sheet, see at [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/).  
Alternatively please contact your distributor or send an e-mail to [info@congen.de](mailto:info@congen.de).