
**CONGEN**

# **SureFast® PREP Bacteria**

Art. No. F1021  
100 extractions

## **User Manual**

Efficient DNA preparation of bacteria

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1.1 Description .....7

1.2 Kit components and storage .....7

1.3 Additionally required equipment and materials .....7

2 Protocol .....8

2.1 Principle .....8

2.2 Preparations .....8

2.3 Extraction protocol .....8

    1. Preparation of the basic material .....8

    2. Lysis of the basic material .....8

    3. Pre-filtration and setting of optimal binding conditions .....9

    4. Binding of the nucleic acids on a Spin Filter .....9

    5. Purification of the bound nucleic acids .....9

    6. Drying of the Spin Filter .....9

    7. Elution of nucleic acids from the Spin Filter .....9

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Kit Code	Reagenz/Material	Menge (je Box)
L	Lysis Buffer	1 x 25 ml
B	Binding Buffer*	1 x 4 ml
W	Wash Buffer**	1 x 18 ml
E	Elution Buffer	1 x 10 ml
S	Spin Filter	1 x 50 Filter
R	Receiver Tubes 2,0 ml	1 x 50 Tubes
T	Receiver Tubes 1,5 ml	1 x 50 Tubes

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## 1 General Information

### 1.1 Description

The kit is intended to be used for the isolation of bacteria - DNA from food (enrichments, avulsions and swabs).

In samples with a very low microbial activity (e.g. heat pasteurized egg white powder) the detection of not viable or not culturable bacteria cells is possible. In case of other conditions the samples show a fast metabolization of these cells.

To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing.

### 1.2 Kit components and storage

Kit Code	Reagent /Material	Amount (per box)
L	Lysis Buffer	1 x 25 ml
B	Binding Buffer*	1 x 4 ml
W	Wash Buffer**	1 x 18 ml
E	Elution Buffer	1 x 10 ml
S	Spin Filter	1 x 50 Filter
R	Receiver Tubes 2,0 ml	1 x 50 Tubes
T	Receiver Tubes 1,5 ml	1 x 50 Tubes

\* Add isopropanol (not supplied with the kit, see 1.3 Additionally required equipment and materials)

\*\*Add ethanol (not supplied with the kit, see 1.3 Additionally required equipment and materials)

All reagents of the kit should be stored dry and at room temperature (14-25°C).

### 1.3 Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 99°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol for preparation of Wash Buffer (puriss, purity  $\geq 96\%$ )
- isopropanol for preparation of Binding Buffer (e.g. Carl Roth 2-Propanol Rotipurán® >99.7 %; Applichem 2-Propanol for microbiology ; Sigma-Aldrich 2-Propanol)
- disposal bags or waste bin
- autoclave
- biological safety cabinet

## 2 Protocol

### 2.1 Principle

1. Preparation of the basic material
2. Lysis of the basic material
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

### 2.2 Preparations

#### General

Add 11 ml isopropanol to the Binding Buffer (**Code B**) and mix thoroughly.

Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

#### Before each preparation

Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer under calculation of a reserve volume into a reaction tube (not provided with the kit).

Equilibrate the Elution Buffer to 60°C. The Elution Buffer is necessary for step 7.

### 2.3 Extraction protocol

#### 1. Preparation of the basic material

Transfer 1.0 ml of the enrichment or avulsions under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).

Centrifuge for 5 min at 12.000 rpm.

Discard the fluid and inactivate per sterilization.

For swabs, place them into a 1.5 ml reaction tube (not provided with the kit) and cut the shaft, so that you can close the cap of the extraction tube.

#### 2. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**) to the reaction tube and mix it briefly and close the reaction tube with caps.

Incubate on a heating block under continuously shaking for 10 minutes at 99°C.

**3. Setting of optimal binding conditions**

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Transfer ca. 300 µl of the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit).

**Note:** Depending on the nature of the sample a smaller amount of liquid may be obtained after centrifugation. In this case take the actual volume for the following steps or vary the volume of the Lysis Buffer before lysis.

**4. Binding of the nucleic acids on a Spin Filter**

Add 200 µl Binding Buffer (**Code B**) to the filtrate and mix.

Place a Spin Filter (**Code S**) into a new 2.0 ml Receiver Tube (**Code R**). Transfer the solution onto the Spin Filter.

Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm.

After centrifugation discard the filtrate and place the Spin Filter back into the Receiver Tube.

**5. Purification of the bound nucleic acids**

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

**6. Drying of the Spin Filter**

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

**7. Elution of nucleic acids from the Spin Filter**

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 100 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min at room temperature.

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

**3 Further Information****3.1 Product Information**

- Validation data upon request
- Material Safety Data Sheet (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))

**3.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).

## 4 Safety information

### Lysis Buffer

#### Warning

H319-H412-P264-P273  
P280-P305+P351+P338-  
P337+P313-P501

- H319:** Causes serious eye irritation.
- H412:** Harmful to aquatic life with long lasting effects.
- P264:** Wash hands, forearms and face thoroughly after handling.
- P273:** Avoid release to the environment.
- P280:** Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
- P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P337+P313:** If eye irritation persists: Get medical advice/attention.
- P501 :** Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation

For further information we offer a Material Safety Data Sheet, see at [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/).  
Alternatively Please contact your distributor or send an e-mail to [info@congen.de](mailto:info@congen.de).