**CONGEN**

# **SureFast<sup>®</sup> PREP Aqua**

Art. No. F1023  
100 extractions

## **User Manual**

Efficient DNA preparation of bacteria DNA from water samples

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**Content**

1	General Information .....	8
1.1	Description .....	8
1.2	Kit components and storage .....	8
1.3	Additionally required equipment and materials .....	8
2	Protocol .....	9
2.1	Principle .....	9
2.2	Preparations .....	9
	General .....	9
	Before each preparation .....	9
2.3	Extraction protocol .....	9
	1. Preparation of the basic material .....	9
	2. Lysis of the basic material .....	9
	3. Setting of optimal binding conditions .....	10
	4. Binding of the nucleic acids on a Spin Filter .....	10
	5. Purification of the bound nucleic acids .....	10
	6. Drying of the Spin Filter .....	10
	7. Elution of nucleic acids from the Spin Filter .....	10
3	Further Information .....	11
3.1	Product Information .....	11
3.2	Technical Support .....	11

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# 1 General Information

## 1.1 Description

This kit is intended to be used for the isolation of bacteria DNA from water samples.

## 1.2 Kit components and storage

Kit Code	Reagent / Material	Amount (per box)
<b>F</b>	Membrane Filter (white with grid; pore size 0.45 µm) ***	50 x
<b>L</b>	Lysis Buffer	1 x 50 ml
<b>B</b>	Binding Buffer*	1 x 4 ml
<b>W</b>	Wash Buffer**	1 x 18 ml
<b>E</b>	Elution Buffer	1 x 10 ml
<b>S</b>	Spin Filter	1 x 50 Filter
<b>R</b>	Receiver Tubes 2.0 ml	1 x 50 Tubes
<b>T</b>	Receiver Tubes 1.5 ml	1 x 50 Tubes

\* Adding isopropanol (not supplied with the kit, see 1.3 additionally required equipment and materials)

\*\* Adding ethanol (not supplied with the kit, see 1.3 additionally required equipment and materials)

\*\*\* The membrane included in the kit is made of mixed cellulose esters. As an alternative for the membrane filter included in the kit, polycarbonate membrane filters can be used for the detection of Legionella spp. according to ISO/TS 12869.

All reagents of the kit should be stored dry and at room temperature (14-25°C).

## 1.3 Additionally required equipment and materials

- membrane filtering system
- reaction tubes free from DNA and DNase 1.5 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- sterile tweezers and spatula
- Vortex mixer
- Thermomixer/ heating block (up to 95°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol for preparation of Wash Buffer (puriss., purity ≥ 96%)
- isopropanol for preparation of Binding Buffer (e.g. Carl Roth 2-Propanol Rotipur® >99.7; Applichem 2-Propanol for microbiologie; Sigma-Aldrich 2-Propanol)
- disposal bags or waste bin
- autoclave
- biological safety cabinet

## 2 Protocol

### 2.1 Principle

1. Membrane filtration
2. Lysis of bacterial cells
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

### 2.2 Preparations

#### General

Add 11 ml isopropanol to the Binding Buffer **(Code B)** and mix thoroughly.

Add 42 ml ethanol to the Wash Buffer **(Code W)** and mix thoroughly.

#### Before each preparation

Preheating the Elution Buffer **(Code E)** - Transfer the needed amount of Elution Buffer **(Code E)** under calculation of a reserve volume into a reaction tube (not provided with the kit).

Equilibrate the Elution Buffer to 60°C. The Elution Buffer is necessary for step 7.

### 2.3 Extraction protocol

#### 1. Preparation of the basic material

Filtrate 500 ml (0.1 – 1.0 l depending on the degree of contamination) of the water sample using the Membrane Filter **(Code F)**. Place the filter (white with grid) with the grid showing upwards onto the filter holder.

**Note:** If the SureFast® PREP Aqua kit is used together with the SureFast® Legionella Screen PLUS real time PCR kit for the detection of *Legionella* spp. in combination with the German drinking water ordinance 2023 limit of 100 CFU/100 ml (Positive Control), 100 ml of the sample have to be filtered.

Fold the filter with a sterile tweezer (because the filter needs to be folded more than once, a sterile spatula to fix the filter might be helpful) and transfer it into a 1.5 ml reaction tube (not supplied with the kit).

#### 2. Lysis of the basic material

Add 1 ml of Lysis Buffer **(Code L)** to the reaction tube and mix it briefly to ensure that the whole filter is completely surrounded by Lysis Buffer. Close the reaction tube with caps.

Incubate on a heating block under continuously shaking for 10 minutes at 95°C.

### 3. Setting of optimal binding conditions

Transfer 130 µl of the lysate into a new 1.5 ml reaction tube (not supplied with the kit).

Add 200 µl Binding Buffer (**Code B**) to the filtrate and mix.

### 4. Binding of the nucleic acids on a Spin Filter

Place a Spin Filter (**Code S**) into a new 2.0 ml Receiver Tube (**Code R**).

Transfer the solution onto the Spin Filter. Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm.

After centrifugation discard the filtrate and place the Spin Filter back into the Receiver Tube.

### 5. Purification of the bound nucleic acids

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

### 6. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

### 7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a 1.5 ml Receiver Tube (**Code T**) and add 50 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min at room temperature.

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

### **3 Further Information**

#### **3.1 Product Information**

- Validation data upon request
- Material Safety Data Sheet (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))

#### **3.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).