CONGEN

SureFast® Speed PREP

Art. No. F1054
100 extractions

User Manual

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1 General Information

1.1 Description

The kit is intended for the rapid and efficient extraction of bacteria- and parasite-DNA from tissue sample and enriched food.

To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

1.2 Kit components and storage

Kit Code	Reagent /Material	Amount (per box)
L	Lysis Buffer	2 x 25 ml

The Lysis Buffer should be stored at -28 to +8 °C.

1.3 Additionally required equipment and materials

- reaction tubes free from DNA and DNase 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 95°C)
- micro centrifuge (up to 12,000 rpm)
- disposal bags or waste bin
- autoclave
- biological safety cabinet

1.4 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures shall be considered as potentially infectious and must be handled and disposed of in accordance with applicable regulations for potentially infectious materials. All reagents and materials that come into contact with such samples shall be treated as biohazardous waste and disposed of accordingly.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

2 Protocol

2.1 Principle

1. Preparation of the basic material
2. Lysis of the basic material

2.2 Extraction protocol

2.2.1 Preparation of the basic material

Extraction of Enrichment Cultures

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 500 µl of the upper third from the enrichment under sterile conditions into a 2 ml reaction tube (not provided with the kit).

Extraction of Tissues

Transfer 50 mg tissue sample into a 2.0 ml reaction tube (not provided with the kit).

2.2.2 Lysis of the basic material

Transfer 500 µl of the upper third from the Lysis Buffer (**Code L**) into a 2.0 ml tube (not provided with the kit) and add 500 µl of the upper third from the enrichment or 50 mg of the tissue sample to the Lysis Buffer (**Code L**).

Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction.

The lysate is ready-to-use for the PCR. If the lysate is not used immediately in the PCR or is intended for a longer storage, transfer 100 µL of the lysate in a new tube (not provided with the kit) and store at -20 °C.

3 Further Information

3.1 Product Information

- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation data upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.