
CONGEN

SureFast® Mag PREP Food

Art. No. F1060
96 extractions

User Manual

Efficient DNA preparation from food and feed

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Content

1	General Information	11
1.1	Description	11
1.2	Principle	11
1.3	Kit components and storage.....	11
1.4	Auto Plate design and contents	11
1.5	Additionally, required equipment and materials	12
1.6	Precautions for users.....	12
2	Protocol.....	13
2.1	Preparations	13
2.1.1	General	13
2.1.2	Before each preparation	13
2.2	Preparation of the basic material	13
2.2.1	Extraction of animal and plant DNA from food, feed and tissues	13
2.2.2	Extraction of bacterial DNA from bacterial culture enrichments	14
2.3	Extraction with TANBead® Maelstrom 8 Autostage	14
2.4	Program for TANBead® Maelstrom 8 Autostage	15
2.5	Extraction with TANBead® Maelstrom 4800	16
2.6	Program for TANBead® Maelstrom 4800.....	16
3	Further Information.....	17
3.1	Product Information	17
3.2	Technical Support	17
4	Safety Information	17

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
Reihe	Reagenz	Menge
1+7	Lysis Buffer	700 µl
2+8	Washing Buffer 1	800 µl
3+9	Magnetic Beads	800 µl
4+10	Washing Buffer 2	800 µl
5+11	Washing Buffer 2	800 µl
6+12	Elution Buffer	130 µl

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1 General Information

1.1 Description

The kit is intended to be used for isolation of animal and plant DNA (desoxyribonucleic acid) from raw materials, low- and high-processed food and feed as well as bacterial DNA from bacterial culture enrichments of food with the TANBead® Maelstrom 8 Autostage or Maelstrom 4800.

1.2 Principle

1. Preparation of starting material and Lysis at 65°C
2. Extraction with TANBead® Maelstrom 8 Autostage or Maelstrom 4800

1.3 Kit components and storage

Kit Code	Reagent	Amount
I	Incubation Buffer	2 x 60 ml
E	Elution Buffer	1 x 20 ml
K	Proteinase K	1 x 1 ml
A	Auto Plates	6 x
Sp	Spin Tips	1 x 96 Tips

All components of the extraction kit should be stored at room temperature (+15 to +35 °C). We recommend to store Proteinase K at +2 to +8 °C after opening.

1.4 Auto Plate design and contents


Figure 2: Auto Plate scheme (16 samples)

sample 1-8 columns A - H 1
sample 9-16 columns A - H 7

Column	Reagent	Volume
1+7	Lysis Buffer	700 µl
2+8	Washing Buffer 1	800 µl
3+9	Magnetic Beads	800 µl
4+10	Washing Buffer 2	800 µl
5+11	Washing Buffer 2	800 µl
6+12	Elution Buffer	130 µl

1.5 Additionally, required equipment and materials

- TANBead® Maelstrom 8 Autostage / Maelstrom 4800 and manual
- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes DNA and DNase-free 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- powder-free disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/heating block (up to 100°C)
- micro centrifuge (up to 12.000 rpm)
- plate centrifuge
- disposal bags or waste bin

1.6 Precautions for users

- This test must only be performed by qualified and authorized laboratory personnel.
- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Strictly follow the working instructions.
- When handling samples, wear suitable protective equipment (disposable gloves, lab coat if necessary). After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled and disposed according to the national safety regulations.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

For more details see Material Safety Data Sheet, at www.congen.de/en/eifu/

2 Protocol

2.1 Preparations

2.1.1 General

- Avoid vigorous shaking of Auto Plate to prevent foam formation.
- Keep Auto Plates and reagents closed. Vaporization results in changes of pH or effects the extraction efficiency.
- Reagents are transparent and colorless except for columns 3+9. These columns contain the Magnetic Beads which are colored black. Changes in color indicates contamination, a fresh Auto Plate should be used.
- Reagents include guanidinium-salt, do not use bleach-containing detergent.

2.1.2 Before each preparation

- In case the room temperature is below 20 °C, the Auto Plate should be pre-heated 5 to 10 minutes at 42-60 °C.
- Spin down the Auto Plate to remove the liquid from the cover.
- Integrity of Auto Plate should be checked before use and Spin Tips have to be in place. Ensure that Spin Tips are firmly fixed on the Maelstrom.

2.2 Preparation of the basic material

2.2.1 Extraction of animal and plant DNA from raw materials, food and feed

50 mg* of a representative, homogenized sample should be weighed into a 2 ml reaction tube (not included in kit).

Add 600 µl Incubation Buffer (**Code I**) and 10 µl Proteinase K (**Code K**) to the sample. Then mix well using the vortex-mixer.

Incubation at 65 °C for 30 min with continuous shaking on a thermo-mixer/heating block.

***Note:** In case of highly processed and watery samples use 100 mg of the sample. For strongly swelling samples it can be necessary to add additional volume of Incubation Buffer during the lysis.

Centrifugation of the lysate for 1 min at 12.000 rpm. The supernatant is processed in the next step

2.3 Extraction with TANBead® Maelstrom 8 Autostage or 2.5 Extraction with TANBead® Maelstrom 4800.

2.2.2 Extraction of bacterial DNA from bacterial culture enrichments

Note: To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of culturing.

2.2.2.1 Protocol 1: Recommended extraction of gram-positive bacteria

Transfer 1.0 ml of bacterial enrichment under sterile conditions into 2 ml reaction tube (not included in kit).

Centrifugation of enrichment for 5 min at 12.000 rpm.

Carefully discard the supernatant and deactivate e.g. by autoclaving.

Add 600 µl Incubation Buffer (**Code I**) and 10 µl Proteinase K (**Code K**) to the pellet. Then mix well using the vortex-mixer.

Incubation at 65 °C for 30 min and following at 99 °C for 10 min under continuously shaking.

Centrifugation of the lysate for 1 min at 12.000 rpm. The supernatant is processed in the next step 2.3 Extraction with TANBead® Maelstrom 8 Autostage or 2.5 Extraction with TANBead® Maelstrom 4800.

2.2.2.2 Protocol 2: Recommended extraction of gram-negative bacteria

Transfer 1.0 ml of bacterial enrichment under sterile conditions into 2 ml reaction tube (not included in kit).

Centrifugation of enrichment for 5 min at 12.000 rpm.

Carefully discard the supernatant and deactivate e.g. by autoclaving.

Add 600 µl Incubation Buffer (**Code I**) to the pellet. Then mix well using the vortex-mixer.

Incubation at 99 °C for 10 min under continuously shaking.

Centrifugation of the lysate for 1 min at 12.000 rpm. The supernatant is processed in the next step 2.3 Extraction with TANBead Maelstrom 8 Autostage or 2.5 Extraction with TANBead® Maelstrom 4800.

Note: It is also possible to extract DNA from gram-negative bacteria using protocol 1.

2.3 Extraction with TANBead® Maelstrom 8 Autostage

1. Switch on the Maelstrom 8 Channel Handler refer to Maelstrom 8 manual.
2. Switch on the TANBead® Maelstrom 8 Autostage refer to Maelstrom 8 manual.
3. Carefully remove the aluminium foil cover of the Auto Plate.
4. Transfer maximal 600 µl of the supernatant into Auto Plate column A-H 1 or 7.
5. Place the Auto Plate in the instrument. The slanted edge should be placed in the left lower corner.
6. **Attach the Spin Tips to the Maelstrom 8 Channel Handler and place the Handler in the Autostage. Ensure that Spin Tips are firmly fixed on the Maelstrom 8 Channel Handler.**
7. Select or adjust the program „6T2-1/7“. The settings are listed below.
8. As soon as the program ends, the Auto Plate can be removed from the instrument.
9. The purified nucleic acid has to be transferred into a fresh reaction tube.
10. The Spin Tips and Auto Plate may now be discarded.

2.4 Program for TANBead® Maelstrom 8 Autostage

For DNA extraction please use program 6T2-1/7.

Note: For changing program you can connect the Maelstrom 8 Channel Handler to the computer with the enclosed micro-USB cable.

Step	Position	Action	RPM	Time [s]	Temperature [°C]	Temperature Control
1	3/9	Mix	3000	60	55	Yes
2	3/9	Collect	0	30	55	Yes
3	2/8	Mix	3000	60	55	Yes
4	1/7	Mix	3000	600	55	Yes
5	2/8	Collect	0	30	55	Yes
6	1/7	Mix	3000	600	55	Yes
7	1/7	Collect	0	30	55	Yes
8	2/8	Mix	3000	300	45	Yes
9	2/8	Collect	0	30	45	Yes
10	3/9	Mix	3000	300	45	Yes
11	3/9	Collect	0	30	45	Yes
12	4/10	Mix	3000	300	45	Yes
13	4/10	Collect	0	30	45	Yes
14	5/11	Mix	3000	300	45	Yes
15	5/11	Collect	0	30	45	Yes
16	5/11	Vapor	0	300	45	Yes
17	6/12	Mix	2700	600	45	Yes
18	6/12	Collect	0	60	45	Yes
19	5/11	Mix	3000	30	0	No

Note:

- The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.
- The Elution Buffer (**Code E**) can be used to dilute the extracted DNA.
- If the eluate contains magnetic beads, spin down the beads for 1 min at 10.000 rpm and transfer the supernatant into a new tube (not provided with the kit).
- SureFood® PREP Advanced (S1053) or SureFood® PREP Basic (S1052) is recommended if magnetic beads get lost during the extraction.

2.5 Extraction with TANBead® Maelstrom 4800

1. Switch on the Maelstrom 4800 refer to Maelstrom 4800 User Manual.
2. Enter the user code to get to the main menu.
3. **Attach the Spin Tips to the device. Ensure that Spin Tips are firmly fixed on the Maelstrom 4800.**
4. Carefully remove the aluminium foil cover of the Auto Plate.
5. Transfer maximal 600 µl of the supernatant into Auto Plate column A-H 1 or 7.
6. Place the Auto Plate in the instrument. The slanted edge should be placed in the left lower corner.
7. Select or adjust the program „6T2-4800“. The settings are listed below.
8. As soon as the program ends, the Auto Plate can be removed from the instrument.
9. The purified nucleic acid has to be transferred into a fresh reaction tube.
10. The Spin Tips and Auto Plate (after extraction of 16 samples) may now be discarded.

2.6 Program for TANBead® Maelstrom 4800

For DNA extraction please use program 6T2-4800.

Note: There is an edit function for programs on the device (refer to Maelstrom 4800 User Manual).

Step	Well	Temperature [°C]	Mixing [min]	Mix speed [rpm]	Collect time [min]	Vapor time [min]	Pause
1	3		1	3000	0.5	0	Off
2	2		1	3000	0	0	Off
3	1	55	10	3000	0	0	Off
4	2		0	500	0.5	0	Off
5	1	Off	10	3000	0.5	0	Off
6	2		5	3000	0.5	0	Off
7	3		5	3000	0.5	0	Off
8	4		5	3000	0.5	0	Off
9	5		5	3000	0.5	5	Off
10	6	Off	10	2700	1	0	Off
11	3		0.5	3000	0	0	Off

Note:

- The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.
- The Elution Buffer (**Code E**) can be used to dilute the extracted DNA.
- If the eluate contains magnetic beads, spin down the beads for 1 min at 10.000 rpm and transfer the supernatant into a new tube (not provided with the kit).
- SureFood® PREP Advanced (S1053) or SureFood® PREP Basic (S1052) is recommended if magnetic beads get lost during the extraction.

3 Further Information

3.1 Product Information

- Product-related documents (Download: www.congen.de/en/eifu/)
- Programming Guide
- TANBead Maelstrom 8 or Maelstrom 4800 manual
- Validation Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.

4 Safety Information

Lysis Buffer

Danger
H302+312+332-314-318-412
P273-280

Washing Buffer 1

Danger
H302+312+332-314-318-412
P273-280

Proteinase K

Danger
H334-317
P101-102-261-280-284-304+340-342+311

H302:	Harmful if swallowed.
H312:	Harmful in contact with skin.
H314:	Causes severe skin burns and eye damage.
H317:	May cause an allergic skin reaction.
H318:	Causes serious eye damage.
H332:	Harmful if inhaled.
H334:	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H412:	Harmful to aquatic life with long lasting effects.
P101:	If medical advice is needed, have product container or label at hand.
P102:	Keep out of reach of children.
P261:	Avoid breathing dust/fume/gas/mist/vapors/spray.
P273:	Avoid release to the environment.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P284:	[In case of inadequate ventilation] wear respiratory protection.
P304+P340:	f Inhaled: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P342+311:	If experiencing respiratory symptoms: Call a POISON CENTER/doctor/physician.

For further information we offer a Material Safety Data Sheet, see at www.congen.de/en/eifu/.

Alternatively please contact your distributor or send an e-mail to info@congen.de.