
CONGEN

SureFast® Mag PREP Pathogens

Art. No. F1062
96 extractions

User Manual

Efficient preparation of viral or bacterial DNA/RNA

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Content

1	General Information	11
1.1	Description	11
1.2	Principle	11
1.3	Kit components and storage.....	11
1.4	Auto Plate design and contents	11
1.5	Additionally required equipment and materials	12
1.6	Precautions for users	12
2	Protocol.....	13
2.1	Preparations	13
2.1.1	General	13
2.1.2	Before each preparation	13
2.2	Preparation of the basic material	13
2.2.1	Extraction of viral DNA/RNA from rinsing liquids	13
2.2.2	Extraction of viral DNA/RNA from swabs and membrane filters	13
2.2.3	Extraction of bacterial DNA from enrichments of food and feed	13
2.3	Extraction with TANBead Maelstrom 8 Autostage	14
2.4	Program for TANBead Maelstrom 8 Autostage	15
2.5	Extraction with TANBead Maelstrom 4800	15
2.6	Program for TANBead Maelstrom 4800	16
3	Further Information.....	16
3.1	Product Information	16
3.2	Technical Support	16
4	Safety Information	17

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Proben 1-8 in Reihe A-H 1

Proben 9-16 in Reihe A-H 7

Reihe	Reagenz	Menge
1+7	Lysis Buffer	600 µl
2+8	Washing Buffer 1	800 µl
3+9	Washing Buffer 2	800 µl
4+10	Washing Buffer 2	800 µl
5+11	Magnetic Beads	800 µl
6+12	Elution Buffer	80 µl

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1 General Information

1.1 Description

The kit is intended to be used for the isolation of DNA and RNA viruses nucleic acids from food (for example wash up fluids from fruits, salads etc.), filters from water samples and swabs (for example from work surfaces etc.) as well as bacterial DNA from bacterial culture enrichments with the TANBead Maelstrom 8 Autostage or TANBead Maelstrom 4800.

1.2 Principle

1. Preparation of starting material
2. Lysis and extraction with TANBead Maelstrom 8 Autostage or TANBead Maelstrom 4800

1.3 Kit components and storage

Kit Code	Reagent	Amount (per box)
E	Elution Buffer	1 x 1.5 ml
K	Proteinase K	1 x 1 ml
A	Auto Plates	6 x
Sp	Spin Tips	1 x 96 Tips

All components of the extraction kit should be stored at room temperature (+15 to +35 °C). We recommend to store Proteinase K at +2 to +8 °C after opening.

1.4 Auto Plate design and contents


Figure 1: Auto Plate scheme (16 samples)

sample 1-8 columns A - H 1

sample 9-16 columns A - H 7

Column	Reagent	Volume
1+7	Lysis Buffer	600 µl
2+8	Washing Buffer 1	800 µl
3+9	Washing Buffer 2	800 µl
4+10	Washing Buffer 2	800 µl
5+11	Magnetic Beads	800 µl
6+12	Elution Buffer	80 µl

1.5 Additionally required equipment and materials

- TANBead Maelstrom 8 Autostage /Maelstrom 4800 and manual
- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes NA and NAse-free 2.0 ml
- NA and NAse-free water
- waterproof pen and tags for labeling the reaction tubes
- powder-free disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/heating block
- micro centrifuge (up to 12,000 rpm)
- micro pipette, filter tips (10 µl, 200 µl, 1000 µl)
- disposal bags or waste bin
- for nucleic extraction from swabs / filters, Tris/HCL buffer

1.6 Precautions for users

- This test must only be performed by qualified and authorized laboratory personnel.
- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Strictly follow the working instructions.
- When handling samples, wear suitable protective equipment (disposable gloves, lab coat if necessary). After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled and disposed according to the national safety regulations.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

For more details see Material Safety Data Sheet, at www.congen.de/en/eifu/

2 Protocol

2.1 Preparations

2.1.1 General

Avoid vigorous shaking of Auto Plate to prevent foam formation.

Keep Auto Plates and reagents closed. Vaporisation results in changes of pH or effects the extraction efficiency.

Reagents are transparent and colorless. Changes in color indicates contamination, a fresh Auto Plate should be used.

Reagents include guanidinium-salt, do not use bleach-containing detergence.

2.1.2 Before each preparation

In case the room temperature is below 20 °C, the Auto Plate should be pre-heated 5 to 10 minutes at 42- 60 °C.

Spin down the Auto Plate to remove the liquid from the cover.

Integrity of Auto Plate should be checked before use and Spin Tips have to be in place. Make sure that the Spin Tips are firmly fixed on the device.

2.2 Preparation of the basic material

2.2.1 Extraction of viral DNA/RNA from rinsing liquids

Use directly 300 µl of liquid sample for the extraction. For samples which have a smaller volume please fill up to a total volume of 300 µl with RNase free PCR grade H₂O.

2.2.2 Extraction of viral DNA/RNA from swabs and membrane filters

Add 400 µl Tris/HCL buffer to the swab/membrane filter. Subsequently shake for 10 minutes in a thermomix at 800 rpm (room temperature). Use 300 µl of the sample in the extraction.

2.2.3 Extraction of bacterial DNA from enrichments of food and feed

Note: To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of culturing.

Use 300 µl from the top 1 to 3 cm of the enrichment culture in the extraction.

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

In order to increase DNA yield, it is possible to pellet 1.0 ml of bacterial enrichment* under sterile conditions after centrifugation at 12.000 rpm for 5 min. Discard the supernatant and resuspend the pellet in 300 µl PCR water. Use 300 µl of the sample in the extraction.

***Note:** To avoid contamination by strongly positive samples, the resuspended pellet can first be heated at 95°C for 10min. Then continue with extraction.

2.3 Extraction with TANBead Maelstrom 8 Autostage

1. Switch on the Maelstrom 8 Channel Handler refer to Maelstrom 8 manual.
2. Switch on the TanBead Maelstrom 8 Autostage refer.
3. Carefully remove the aluminum foil cover of the Auto Plate (Code A).
4. Transfer maximally 300 µl of the sample together with 10 µl Proteinase K (Code K) into Auto Plate (Code A) column A-H 1 or column 1 or 7.
5. Place the Auto Plate in the instrument. The slanted edge should be placed in the left lower corner.
6. **Attach the Spin Tips to the Maelstrom 8 Channel Handler and place the Handler in the Autostage. Ensure that the Spin Tips are firmly fixed on the Maelstrom 8 Channel Handler.**
7. Select or adjust the program „665-1/7“. The settings are listed on the next page.
8. As soon as the program ends, the Auto Plate can be removed from the instrument.
9. The purified nucleic acid has to be transferred into a fresh reaction tube.
10. The Spin Tips and Auto Plate (after extraction of 16 samples) may now be discarded.

Note: The eluted DNA/RNA is ready-to-use for the PCR. The DNA/RNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.

The Elution Buffer (**Code E**) can be used to dilute the extracted DNA/RNA.

2.4 Program for TANBead Maelstrom 8 Autostage

For DNA/RNA extraction please use program 665-1/7.

Note: For changing program you can connect the Maelstrom 8 Channel Handler to the computer with the enclosed micro-USB cable.

Step	Position	Action	RPM	Time [s]	Temperature [°C]	Temperature Control
1	5/11	Collect	0	30	60	Yes
2	1/7	Mix	3000	480	60	Yes
3	1/7	Collect	0	30	60	Yes
4	2/8	Mix	3000	60	45	Yes
5	2/8	Collect	0	30	45	Yes
6	3/9	Mix	3000	60	45	Yes
7	3/9	Collect	0	30	45	Yes
8	4/10	Mix	3000	60	45	Yes
9	4/10	Collect	0	30	45	Yes
10	4/10	Vapor	0	300	45	Yes
11	6/12	Mix	3000	300	45	Yes
12	6/12	Collect	0	30	45	Yes
13	5/11	Mix	3000	60	0	No

2.5 Extraction with TANBead Maelstrom 4800

1. Switch on the Maelstrom 4800 refer to Maelstrom 4800 User Manual.
2. Enter the user code to get to the main menu.
3. **Attach the Spin Tips to the device. Ensure that Spin Tips are firmly fixed on the Maelstrom 4800.**
4. Carefully remove the aluminum foil cover of the Auto Plate.
5. Transfer maximally 300 µl of the sample together with 10 µl Proteinase K (Code K) into Auto Plate column A-H 1 or 7.
6. Place the Auto Plate in the instrument. The slanted edge should be placed in the left lower corner.
7. Select or adjust the program „665-1/7“. The settings are listed under 2.6 Program for TANBead Maelstrom 4800.
8. As soon as the program ends, the Auto Plate can be removed from the instrument.
9. The purified nucleic acid has to be transferred into a fresh reaction tube.
10. The Spin Tips and Auto Plate (after extraction of 16 samples) may now be discarded.

Note: The eluted DNA/RNA is ready-to-use for the PCR. The DNA/RNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.

The Elution Buffer (**Code E**) can be used to dilute the extracted DNA/RNA.

2.6 Program for TANBead Maelstrom 4800

For DNA/RNA extraction please use program 665-1/7

Note: For changing program you can connect the Maelstrom 4800 Channel Handler to the computer with the enclosed micro-USB cable.

Step	Position	Temperature [°C]	Mixing [min]	Mix [RPM]	Collect time [min]	Vapor time [min]	Pause
1	5		0.5	2500	0.5	0	Off
2	1	55	12	1500	0.5	0	Off
3	2		2	2500	0.5	0	Off
4	3		1	1500	0.5	0	Off
5	4		1	1500	0.5	10	Off
6	6		5	2500	1.5	0	Off
7	3		0.5	2500	0	0	Off

3 Further Information

3.1 Product Information

- Product-related documents (Download: www.congen.de/en/eifu/)
- Programming Guide
- TANBead Maelstrom 8 or Maelstrom 4800 manual
- Validation Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.

4 Safety Information

Lysis Buffer

Danger

H302+312+332-314-318-412
P273-280

Washing Buffer 1

Danger

H302+312+332+314-318-412
P273-280

Proteinase K

Danger

H334-317
P101-102-261-280-284-304+340-342+311

H302:	Harmful if swallowed.
H312:	Harmful in contact with skin.
H314:	Causes severe skin burns and eye damage.
H317:	May cause an allergic skin reaction.
H318:	Causes serious eye damage.
H332:	Harmful if inhaled.
H334:	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H412:	Harmful to aquatic life with long-lasting effects.
P101:	If medical advice is needed, have product container or label at hand.
P102:	Keep out of reach of children.
P261:	Avoid breathing dust/fumes/gas/mist/vapors/spray.
P273:	Avoid release to the environment.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P284:	[In case of inadequate ventilation] wear respiratory protection.
P304+340:	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P342+311:	If experiencing respiratory symptoms: Call a POISON CENTER/doctor/physician.

For further information we offer a Material Safety Data Sheet, see at www.congen.de/en/eifu/.

Alternatively please contact your distributor or send an e-mail to info@congen.de.