
CONGEN

SureFast® Salmonella PLUS

Art. No. F5111
100 rxn

User Manual

January 2023

 **Inhalt**

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Content

1	General Information	8
1.1	Description	8
1.2	Limit of Detection	8
1.3	DNA-preparation	9
1.4	Kit components and storage	9
1.5	Additionally required equipment and materials	9
1.6	Micro centrifuge with a rotor for the reaction tubes Precautions for users	9
1.7	Setup.....	10
1.8	Detection channel set-up	10
2	Qualitative Analysis	11
2.1	Protocol	11
2.1.1	Preparation of the master-mix	11
2.1.2	Preparation of the real-time PCR-mix	11
2.2	Interpretation of results	12
3	Further Information	12
3.1	Product Information	12
3.2	Technical Support	12

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1 General Information

1.1 Description

The SureFast® Salmonella PLUS is a real-time PCR for the direct, qualitative detection of a specific DNA sequence of *Salmonella* spp.. The SureFast® Salmonella PLUS real-time PCR kit has been validated and certified in combination with the SureFast® PREP Salmonella kit by AOAC (licence number: 041103).

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed.

The real-time PCR assay can be used with real-time PCR instruments, equipped for detection of two fluorescence emissions in channels FAM and VIC/HEX at the same time.

The technical validation of instruments was performed on Agilent AriaMx, Roche LightCycler® 480 II, Applied Biosystems® 7500, Bio-Rad CFX96 und Bio Molecular Systems MIC.

The internal technical verification of instruments (not part of the PTM validation) was performed on Roche LightCycler® 2.0, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER, Cepheid SmartCycler, LTF MyGo Pro and Agilent Mx3005P.

1.2 Limit of Detection

The SureFast® Salmonella PLUS real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Salmonella or SureFast® Speed PREP is recommended. To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µL	Yellow
2	Taq Polymerase	1 x 11 µL	Red
3	Positive Control	1 x 200 µL	Light Blue

Store all reagents at –20°C and protected from light.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFast® PREP Salmonella Art. No. F1007 / Speed PREP Art. No. F1054)
- Real-time PCR instrument with two detection channels (510 nm and 580 nm)
- Real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- Pipettes with filter tips
- Unpowdered disposable gloves
- Vortex mixer

1.6 Micro centrifuge with a rotor for the reaction tubes Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II & LTF MyGo Pro
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	<i>Salmonella</i> spp.	FAM	+	
	IAC	HEX	+	
Applied Biosystems 7500	<i>Salmonella</i> spp.	FAM	None	Use the passive reference dye ROX setting "none".
	IAC	VIC	None	
Bio-Rad CFX96/Dx	<i>Salmonella</i> spp.	FAM	+	
	IAC	VIC/HEX	+	
BMS MIC	<i>Salmonella</i> spp.	green	+	
	IAC	yellow	+	
LTF MyGo Pro	<i>Salmonella</i> spp.	FAM	+	
	IAC	VIC	+	
Qiagen Rotor-Gene Q	<i>Salmonella</i> spp.	green	+	Note: Please use only 0.1 mL reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
R-Biopharm RIDA®CYCLER	<i>Salmonella</i> spp.	green	+	
	IAC	yellow	+	
Roche LightCycler® 2.0	<i>Salmonella</i> spp.	530	None	The SureCC Color Compensation Kit II (Art. No. F4010) is required. Set the Seek Temperature to 58°C.
	IAC	560	None	
Roche LightCycler® 480 II	<i>Salmonella</i> spp.	465-510	+	
	IAC	533-580	+	
Roche cobas® z 480 Analyzer	<i>Salmonella</i> spp.	465-510	+	
	IAC	540-580	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative *Salmonella* spp. detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative *Salmonella* spp. detection in enrichments:

5 reactions for controls (1x no-template control, 1x extraction control, 1x positive control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µL	218.9 µL
Taq Polymerase	0.1 µL	1.1 µL
Total volume	20 µL	220 µL

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µL of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µL of sample DNA into the designated tubes/wells or capillaries.
- Pipette 5 µL of Positive Control into the designated tubes/wells or capillaries.
- Centrifuge all tubes/wells or capillaries shortly at low speed.
- Place tubes/ wells or capillaries into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

Salmonella spp. DNA is detected in the FAM-channel. In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for *Salmonella* spp., if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for *Salmonella* spp., if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value ≤ 2 compared to the negative control. If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances' DNA isolation and purification of the sample need to be improved. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Salmonella* spp. PCR assay.

3 Further Information

3.1 Product Information

Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage. Download: www.congen.de/en/company/downloads

- Validation Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de