
CONGEN

SureFast[®] Salmonella PLUS

Art. No. F5111
100 rxn

User Manual

January 2025

 **Inhalt**

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Content

1	General Information	10
1.1	Description	10
1.2	Limit of Detection	10
1.3	DNA-preparation	11
1.4	Kit components and storage	11
1.5	Additionally required equipment and materials	11
1.6	Precautions for users	11
1.7	Setup.....	12
1.8	Detection channel set-up	12
2	Qualitative Analysis	14
2.1	Protocol	14
2.1.1	Preparation of the master-mix	14
2.1.2	Preparation of the real-time PCR-mix	14
2.2	Interpretation of results	15
3	Validation	15
4	Limitations of the method	16
5	Further Information	16
5.1	Product Information	16
5.2	Technical Support	16

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1 General Information

1.1 Description

The SureFast® Salmonella PLUS is a real-time PCR for the direct, qualitative detection of a specific DNA sequence of *Salmonella* spp.. The SureFast® Salmonella PLUS real-time PCR kit has been validated and certified in combination with the SureFast® PREP Bacteria kit by AOAC (licence number: 041103).

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed.

The real-time PCR assay can be used with real-time PCR instruments, equipped for detection of two fluorescence emissions in channels FAM and VIC/HEX at the same time.

The technical validation of instruments was performed on Agilent AriaMx, Roche LightCycler® 480 II, Applied Biosystems® 7500, Bio-Rad CFX96 und Bio Molecular Systems MIC.

The internal technical verification of instruments (not part of the PTM validation) was performed on Roche LightCycler® 2.0, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER, Cepheid SmartCycler, LTF MyGo Pro and Agilent Mx3005P.

1.2 Limit of Detection

The SureFast® Salmonella PLUS real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria (Art. No. F1021) or SureFast® Speed PREP (Art. No. F1054) is recommended. To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µL	Yellow
2	Taq Polymerase	1 x 11 µL	Red
3	Positive Control	1 x 200 µL	Light Blue

Store all reagents at -28 to -16°C and protected from light.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFast® PREP Bacteria Art. No. F1021 / Speed PREP Art. No. F1054)
- Real-time PCR instrument with two detection channels (510 nm and 580 nm)
- Real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- Pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- Micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II & LTF MyGo Pro
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	<i>Salmonella</i> spp.	FAM	+	
	IAC	HEX	+	
Applied Biosystems 7500	<i>Salmonella</i> spp.	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
Bio-Rad CFX96 / Dx / Opus	<i>Salmonella</i> spp.	FAM	+	
	IAC	VIC/HEX	+	
Bio Molecular Systems MIC	<i>Salmonella</i> spp.	green	+	
	IAC	yellow	+	
LTF MyGo Pro	<i>Salmonella</i> spp.	FAM	+	
	IAC	VIC	+	
Qiagen Rotor-Gene Q	<i>Salmonella</i> spp.	green	+	Note: Please use only 0.1 mL reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	

Real-time PCR device	Detection	Detection channel	Quencher	Note
R-Biopharm RIDA®CYCLER	<i>Salmonella</i> spp.	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter
	IAC	yellow	+	
Roche LightCycler® 480 II	<i>Salmonella</i> spp.	465-510	+	
	IAC	533-580	+	
Roche cobas® z 480 Analyzer	<i>Salmonella</i> spp.	465-510	+	
	IAC	540-580	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive Control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control.

The reaction mix contains an internal amplification control (inhibition control) per reaction.

Reactions needed for the qualitative *Salmonella* spp. detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative *Salmonella* spp. detection in enrichments:

5 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

The tube of the Taq Polymerase should be kept at -28 to -16°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µL	218.9 µL
Taq Polymerase	0.1 µL	1.1 µL
Total volume	20 µL	220 µL

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µL of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µL of sample DNA into the designated tubes/wells.
Close the tubes.
- Pipette 5 µL of Positive Control into the designated tubes/wells.
Close the tubes.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Salmonella spp. DNA is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** for *Salmonella* spp., if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for *Salmonella* spp., if the sample DNA shows no amplification in the FAM-channel and if the IAC (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value ≤ 2 compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Salmonella* spp. PCR assay.

3 Validation

SureFast® Salmonella PLUS was validated per the AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces and the ISO 16140-2:2016 [Microbiology of the food chain -- Method validation -- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method]. The reference method for the matrix study was ISO 6579:2002 (Microbiology of food and animal feeding stuffs -- Horizontal method for the detection of *Salmonella* spp.).

The kit has been AOAC RI validated for the following food matrices: Salami, pork minced meat, bacon, fresh ground chicken, fresh chicken carcass, frozen marinated chicken fillets, semi-skim milk powder, raw goat milk cheese, chocolate ice cream, salad with mayonnaise, paella, cream-based pastry, dehydrated poultry proteins, wheat-based flour, raw ground beef, frozen poultry meat, raw milk, fresh spinach, pasteurized liquid whole egg and pet food pellets.

4 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- A positive test result does not necessarily indicate the presence of viable organism. It is indicative for the presence of the target DNA (*Salmonella* spp. DNA).

5 Further Information

5.1 Product Information

- Detailed information on the setup of several real-time PCR devices (Download: www.congen.de/en/downloads).
- Validation Report upon request

5.2 Technical Support

For further questions please send an e-mail to info@congen.de