
**CONGEN**

**SureFast®**  
**Emetic Bacillus cereus PLUS**

Art. No. F5127  
100 rxn

**User Manual**

**May 2026**

 **Inhalt**

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## **Content**

1	General Information .....	10
1.1	Description .....	10
1.2	Limit of Detection .....	10
1.3	DNA-preparation .....	11
1.4	Kit components and storage .....	11
1.5	Additionally required equipment and materials .....	11
1.6	Precautions for users .....	11
1.7	Setup.....	12
1.8	Detection channel Set-up .....	12
2	Qualitative Analysis .....	14
2.1	Protocol .....	14
2.1.1	Preparation of the master-mix .....	14
2.1.2	Preparation of the real-time PCR-mix .....	15
2.2	Interpretation of results .....	15
3	Limitations of the method .....	16
4	Further Information .....	16
4.1	Product Information .....	16
4.2	Technical Support .....	16

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## **1 General Information**

### **1.1 Description**

The SureFast® Emetic *Bacillus cereus* PLUS is a real-time PCR for the direct, qualitative detection of a specific cereulide synthetase DNA sequence of emetic *Bacillus cereus* in enriched food.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The internal technical verification of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Applied Biosystems 7500, Bio-Rad CFX96, R-Biopharm RIDA®CYCLER, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Roche LightCycler® PRO and Qiagen Rotor-Gene Q.

### **1.2 Limit of Detection**

The SureFast® Emetic *Bacillus cereus* PLUS real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria (Art. No. F1021) or Speed PREP (Art. No. F1054) is recommended.

To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

For cultivation the method of the DIN EN ISO 7932:2020-11 is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit  
(e.g. SureFast® PREP Bacteria Art. No. F1021 / Speed PREP Art. No. F1054)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- Incubator
- Enrichment broth
- micro centrifuge with a rotor for the reaction tubes

### 1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled and disposed according to the national safety regulations.
- Do not use the kit after the expiration date.
- All reagents and materials used must be disposed properly after use. Please refer to the relevant national regulation for disposal.

**1.7 Setup**

	<b>Blockcycler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler &amp; LightCycler® 480 II</b>
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.8 Detection channel Set-up**

For additional details see cycler operating instructions for real-time PCR device.

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	Emetic <i>Bacillus cereus</i>	FAM	+	
	IAC	HEX	+	
<b>Agilent AriaDx /Mx</b>	Emetic <i>Bacillus cereus</i>	FAM	+	
	IAC	HEX	+	
<b>Applied Biosystems 7500</b>	Emetic <i>Bacillus cereus</i>	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
<b>Bio-Rad CFX96/Dx/ Opus</b>	Emetic <i>Bacillus cereus</i>	FAM	+	Baseline Settings: <ul style="list-style-type: none"> <li>• Baseline subtracted curve fit</li> <li>• Apply fluorescence drift correction</li> </ul>
	IAC	VIC/HEX	+	
<b>R-Biopharm RIDA®CYCLER</b>	Emetic <i>Bacillus cereus</i>	green	+	<b>Ignore cycles before</b> , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	

# SureFast® Emetic *Bacillus cereus* PLUS (100 rxn)

Art. No. F5127

May 2026

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Qiagen Rotor-Gene Q</b>	Emetic <i>Bacillus cereus</i>	green	+	<b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
<b>Roche LightCycler® 480 II</b>	Emetic <i>Bacillus cereus</i>	465-510	+	
	IAC	533-580	+	
<b>Roche LightCycler® PRO</b>	Emetic <i>Bacillus cereus</i>	465-510	+	
	IAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	Emetic <i>Bacillus cereus</i>	465-510	+	
	IAC	540-580	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive Control.

For the analysis of enrichments additional controls are recommended: zero control and medium control.

The reaction mix contains an internal amplification control (inhibition control) per reaction.

#### **Reactions needed for the qualitative Emetic *Bacillus cereus* detection:**

3 reactions for controls\* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

#### **Reactions needed for the qualitative Emetic *Bacillus cereus* detection in enrichments:**

5 reactions for controls\* (1x negative control, 1x extraction control, 1x Positive Control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### **Example for the calculation and preparation of 10 reactions:**

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### **\* Description of the controls**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control
- Zero control: master-mix and sample prior enrichment
- Medium control: only medium – no sample

**2.1.2 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

**2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Emetic *Bacillus cereus* DNA is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** for Emetic *Bacillus cereus*, if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for Emetic *Bacillus cereus*, if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC/HEX-channel) of the sample is positive.

If the sample DNA shows no amplification or sigmoidal amplification curve in the VIC/HEX channel, this indicates that inhibitors are present in the sample DNA that inhibit the PCR. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to repeat the extraction process. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific PCR assay.

**The following table shows the specification ranges of the kit controls**

	<b>Specification range</b>	<b>Internal amplification control (IAC)</b>
negative control	negative	$28 \leq C_p \leq 36$
Positive Control (FAM – Emetic <i>Bacillus cereus</i> )	$25 \leq C_p \leq 33$	not relevant

### **3 Limitations of the method**

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.
- A positive test result does not necessary indicate the presence of viable organism. It is indicative for the presence of the target DNA (*Emetic Bacillus cereus* DNA).

### **4 Further Information**

#### **4.1 Product Information**

- Detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/downloads](http://www.congen.de/en/downloads))
- Product-related documents (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))
- Validation Report upon request

#### **4.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).