
**CONGEN**

**SureFast®**  
**Clostridium estertheticum PLUS**

Art. No. F5160  
100 rxn

**User Manual**

**April 2024**

 **Inhalt**

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## **Content**

1	General Information .....	8
1.1	Description .....	8
1.2	Limit of Detection .....	8
1.3	DNA-preparation .....	9
1.4	Kit components and storage .....	9
1.5	Additionally required equipment and materials .....	9
1.6	Precautions for users .....	9
1.7	Setup.....	10
1.8	Detection channel Set-up .....	10
2	Qualitative Analysis .....	11
2.1	Protocol .....	11
2.1.1	Preparation of the master-mix .....	11
2.1.2	Preparation of the real-time PCR-mix .....	11
2.2	Interpretation of results .....	12
3	Limitations of the method .....	12
4	Further Information .....	12
4.1	Product Information .....	12
4.2	Technical Support .....	12

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## 1 General Information

### 1.1 Description

The SureFast® Clostridium estertheticum PLUS is a real-time PCR for the direct, qualitative detection of a specific *Clostridium estertheticum* DNA sequence in food.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Bio-Rad CFX96, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Applied Biosystems 7500, Qiagen Rotor-Gene Q, LTF MyGo Pro and R-Biopharm RIDA®CYCLER.

### 1.2 Limit of Detection

The SureFast® Clostridium estertheticum DNA PLUS real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria (Art. No. F1021) is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at -28 to -16 °C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Bacteria Art. No. F1021)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

### 1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

# SureFast® Clostridium estertheticum PLUS (100 rxn)

Art. No. F5160

April 2024

## 1.7 Setup

	<b>Blockcycler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler &amp; LightCycler® 480 II &amp; LTF MyGo Pro</b>
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

## 1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Agilent Mx3005P</b>	<i>Clostridium estertheticum</i>	FAM	+	
	IAC	HEX	+	
<b>Agilent AriaDx/Mx</b>	<i>Clostridium estertheticum</i>	FAM	+	
	IAC	HEX	+	
<b>Applied Biosystems 7500</b>	<i>Clostridium estertheticum</i>	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
<b>BioRad CFX96/Dx/Opus</b>	<i>Clostridium estertheticum</i>	FAM	+	
	IAC	VIC	+	
<b>R-Biopharm RIDA®CYCLER</b>	<i>Clostridium estertheticum</i>	green	+	
	IAC	yellow	+	
<b>LTF MyGo Pro</b>	<i>Clostridium estertheticum</i>	FAM	+	
	IAC	VIC	+	
<b>Qiagen Rotor-Gene Q</b>	<i>Clostridium estertheticum</i>	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
<b>Roche LightCycler® 480 II</b>	<i>Clostridium estertheticum</i>	465-510	+	
	IAC	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	<i>Clostridium estertheticum</i>	465-510	+	
	IAC	540-580	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (Inhibition control) per reaction.

#### **Reactions needed for the qualitative *Clostridium estertheticum* detection:**

3 reactions for controls (1x negative control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

#### **Example for the calculation and preparation of 10 reactions:**

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

*Clostridium estertheticum* DNA is detected in the FAM-channel. In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for *Clostridium estertheticum*, if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for *Clostridium estertheticum*, if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value  $\leq 2$  compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Clostridium estertheticum* PCR assay.

## 3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- A positive test result does not necessarily indicate the presence of viable organism. It is indicative for the presence of the target DNA (*Clostridium estertheticum* DNA).

## 4 Further Information

### 4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request

### 4.2 Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).