
CONGEN

SureFast® Salmonella ONE

Art. No. F5211
100 rxn

User Manual

January 2023

 **Inhalt**

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Content

| | | |
|-------|---|----|
| 1 | General Information | 11 |
| 1.1 | Description | 11 |
| 1.2 | Limit of Detection | 11 |
| 1.3 | Kit components and storage | 12 |
| 1.4 | Additionally required equipment and materials | 12 |
| 1.5 | Precautions for users | 12 |
| 1.6 | Setup..... | 13 |
| 1.7 | Detection channel Set-up | 13 |
| 2 | Qualitative Analysis | 14 |
| 2.1 | Protocol | 14 |
| 2.1.1 | Cultural Enrichment | 14 |
| 2.1.2 | DNA-preparation..... | 14 |
| 2.1.3 | Preparation of the master-mix | 15 |
| 2.1.4 | Preparation of the real-time PCR-mix | 15 |
| 2.2 | Interpretation of results | 16 |
| 2.3 | Cultural Confirmation | 16 |
| 3 | Validation | 16 |
| 4 | Further Information | 17 |
| 4.1 | Product Information..... | 17 |
| 4.2 | Technical Support | 17 |

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1 General Information

1.1 Description

The SureFast® Salmonella ONE can be applied for the fast and simple isolation and detection of *Salmonella* spp. In different types of food and is intended to be used by all kinds of food testing laboratories.

The method consists of the following steps: 1) cultural enrichment, 2) DNA-preparation, 3) real-time PCR detection and 4) data analysis and interpretation. SureFast® Salmonella ONE has been validated and certified by MicroVal (License-No. 2014LR43) and AOAC (License-No. 081803).

Each reaction contains an internal amplification control (IAC). The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time.

The technical validation of instruments was performed on Agilent AriaMx, Roche LightCycler® 480 II, Applied Biosystems® 7500, Bio-Rad CFX96 und Bio Molecular Systems MIC.

The internal technical verification of instruments (not part of the PTM validation) was performed on Roche LightCycler® 2.0, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER, Cepheid SmartCycler, LTF MyGo Pro and Agilent Mx3005P.

1.2 Limit of Detection

The SureFast® Salmonella ONE real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 Kit components and storage

| Kit Code | Reagent | Amount | Lid Color |
|----------|------------------|-------------|------------|
| 1 | Reaction Mix | 2 x 1100 µL | Yellow |
| 2 | Taq Polymerase | 1 x 80 µL | Dark Red |
| 3 | Positive Control | 1 x 200 µL | Light Blue |
| L | Lysis Buffer | 2 x 25 mL | Clear |

Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

The Lysis Buffer should be stored at -20 to +8 °C. Reagents should be used until the expiry date indicated on the label.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.4 Additionally required equipment and materials

- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- 2 mL tubes
- powder-free disposable gloves
- Vortex mixer
- centrifuge with a rotor for the reaction tubes
- Thermomixer/heating block (up to 95 °C)
- Non-selective enrichment broth (e. g. buffered peptone water)

1.5 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

SureFast® Salmonella ONE (100 rxn)

Art. No. F5211

January 2023

1.6 Setup

| | Blockcyler & R-Biopharm RIDA®CYCLER | Rotorcyler & LightCycler® 480 II & LTF MyGo Pro |
|--|--|--|
| Initial Denaturation (HOLD) | 1 min, 95°C | 1 min, 95°C |
| Cycles | 45 | 45 |
| Denaturation | 15 sec, 95°C | 10 sec, 95°C |
| Annealing/Extension (CYCLE) | 30 sec, 60°C | 15 sec, 60°C |
| Temperature Transition Rate/ Ramp Rate | Maximum | Maximum |

1.7 Detection channel Set-up

| Real-time PCR device | Detection | Detection channel | Quencher | Note |
|------------------------------------|------------------------|--------------------------|-----------------|--|
| Agilent Mx3005P | <i>Salmonella</i> spp. | FAM | + | |
| | IAC | HEX | + | |
| Applied Biosystems 7500 | <i>Salmonella</i> spp. | FAM | None | Check the passive reference option ROX is none. |
| | IAC | VIC | None | |
| Bio Molecular Systems MIC | <i>Salmonella</i> spp. | green | + | |
| | IAC | yellow | + | |
| Bio-Rad CFX96 / Dx | <i>Salmonella</i> spp. | FAM | + | |
| | IAC | VIC/HEX | + | |
| R-Biopharm RIDA®CYCLER | <i>Salmonella</i> spp. | green | + | |
| | IAC | yellow | + | |
| LTF MyGo Pro | <i>Salmonella</i> spp. | FAM | + | |
| | IAC | VIC | + | |
| Qiagen Rotor-Gene Q | <i>Salmonella</i> spp. | green | + | Note: Please use only 0.1 mL reaction tubes. The gain settings must be set to 5 (factory default) for all channels. |
| | IAC | yellow | + | |
| Roche LightCycler® 2.0 | <i>Salmonella</i> spp. | 530 | None | The SureCC Color Compensation Kit II (Art. No. F4010) is required. Check the Seek Temperature is 58°C. |
| | IAC | 560 | None | |
| Roche LightCycler® 480 II | <i>Salmonella</i> spp. | 465-510 | + | |
| | IAC | 533-580 | + | |
| Roche cobas® z 480 Analyzer | <i>Salmonella</i> spp. | 465-510 | + | |
| | IAC | 540-580 | + | |

2 Qualitative Analysis

2.1 Protocol

2.1.1 Cultural Enrichment

Sample preparation of food and feed matrices should apply to ISO 6887 and ISO 6579 standards. For enrichment, a volume of 25 g or 25 mL sample material is homogenized with 225 mL of a non-selective enrichment broth (e. g. buffered peptone water). If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio. Enrichment is performed for 18 h ± 2 h at 37 °C ± 1 °C. Alternatively, other suitable, validated enrichment procedures can be used.

To assess the process of bacterial growth, it is recommended to compare samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

2.1.2 DNA-preparation

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 500 µL of the upper third from the enrichment culture into a 2.0 mL tube (not provided with the kit) and add 500 µL of Lysis Buffer (**Code L**) to this sample.

Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction.

The lysate is ready-to-use for PCR. If the lysate is not used immediately in the PCR or is intended for a longer storage, transfer 100 µL of the lysate in a new tube (not provided with the kit) and store at -20 °C.

2.1.3 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative *Salmonella* spp. detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative *Salmonella* spp. detection in enrichments:

5 reactions for controls (1x no-template control, 1x extraction control, 1x positive control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

| Components of the master-mix | Amount per reaction | 10 reactions (with 10% excess) |
|------------------------------|---------------------|--------------------------------|
| Reaction Mix | 19.3 µL | 212.3 µL |
| Taq Polymerase | 0.7 µL | 7.7 µL |
| Total volume | 20 µL | 220 µL |

Mix each master-mix well and centrifuge shortly before use.

2.1.4 Preparation of the real-time PCR-mix

- Pipette 20 µL of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µL of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µL of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Salmonella spp. DNA is detected in the FAM-channel. In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for *Salmonella* spp., if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for *Salmonella* spp., if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value ≤ 2 compared to the negative control. If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Salmonella* spp. PCR assay.

Note: During the MicroVal study increased inhibition levels were observed for fresh herbs, cocoa, spices and aromatics and egg products. In these cases, please proceed as recommended above. If you require further assistance please contact info@congen.de.

2.3 Cultural Confirmation

Positive PCR results should be confirmed by cultural detection methods (e.g. ISO 6579).

3 Validation

SureFast® Salmonella ONE was validated per the AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces and the ISO 16140-2:2016 [Microbiology of the food chain -- Method validation -- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method]. The reference method for the matrix study was ISO 6579:2002 (Microbiology of food and animal feeding stuffs -- Horizontal method for the detection of *Salmonella* spp.).

The kit is an extension of the SureFast® Salmonella PLUS method and has been AOAC RI validated for the following additional food matrices: raw ground beef, frozen poultry meat, raw milk, fresh spinach, pasteurized liquid whole egg, pet food pellets.

The MicroVal method comparison study according to ISO 16140-2:2016 was carried out using 25 g portions of sample material. Sample preparations were done according to ISO 6887 parts. In the sensitivity study a total of 391 samples from six different food categories were investigated. SureFast® Salmonella ONE has been MicroVal validated for the following food categories: poultry meat, meat products, dairy products, vegetables (excluding sprouts), egg products, feed.

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Verification Report upon request

4.2 Technical Support

For further questions please send an e-mail to info@congen.de.