CONGEN


SureFast® Salmonella ONE

Art. No. F5211
100 rxn

User Manual

MICROVAL 
European validation and certification organisation

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Content

1	General Information	12
1.1	Description	12
1.2	Limit of Detection	12
1.3	Kit components and storage	13
1.4	Additionally required equipment and materials	13
1.5	Precautions for users	13
1.6	Setup.....	14
1.7	Detection channel Set-up	14
2	Qualitative Analysis	16
2.1	Protocol	16
2.1.1	Cultural Enrichment	16
2.1.2	DNA-preparation.....	16
2.1.3	Preparation of the master-mix	17
2.1.4	Preparation of the real-time PCR-mix	17
2.2	Interpretation of results	18
2.3	Cultural Confirmation	18
3	Validation	19
4	Limitations of the method	19
5	Further Information	19
5.1	Product Information	19
5.2	Technical Support	19

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1 General Information

1.1 Description

The SureFast® Salmonella ONE can be applied for the fast and simple isolation and detection of *Salmonella* spp. in different types of food and is intended to be used by all kinds of food testing laboratories.

The method consists of the following steps: 1) cultural enrichment, 2) DNA-preparation, 3) real-time PCR detection and 4) interpretation of results. SureFast® Salmonella ONE has been validated and certified by MicroVal (License-No. 2014LR43) and AOAC (License-No. 081803).

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time.

The technical validation of instruments was performed on Agilent AriaMx, Roche LightCycler® 480 Instrument II, Applied Biosystems® 7500, Bio-Rad CFX96 and Bio Molecular Systems MIC.

The internal technical verification of instruments (not part of the PTM validation) was performed on Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER and Agilent Mx3005P.

1.2 Limit of Detection

The SureFast® Salmonella ONE real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
L	Lysis Buffer	2 x 25 ml	Clear

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

The Lysis Buffer should be stored at -28 to +8°C. Reagents should be used until the expiry date indicated on the label.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.4 Additionally required equipment and materials

- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- SureFast® PREP Bacteria (Art. No. F1021), optional
- real-time PCR consumables (plates, tubes, foils, caps)
- incubator (up to 37°C)
- pipettes with filter tips, for volumes from 0.7 µL – 1000 µL
- microcentrifuge tubes - 2 ml tubes with lock cap
- powder-free disposable gloves
- vortex mixer
- microcentrifuge (up to 12,000 rpm) with a rotor for the reaction tubes
- thermomixer/heating block (up 99°C)
- enrichment broth (i.e. BPW)
- suitable equipment for sample comminution and homogenization
- balance and spatula for weighing the samples
- autoclave
- biological safety cabinet

1.5 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.
- Dispose of contents to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulations.

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER ¹	Rotorcycler & LightCycler® 480 Instrument II ²
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P¹	<i>Salmonella</i> spp.	FAM	+	
	IAC	HEX	+	
Applied Biosystems 7500¹	<i>Salmonella</i> spp.	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
Bio-Rad CFX96 / Dx / Opus (incl. Deep Well)¹	<i>Salmonella</i> spp.	FAM	+	Baseline settings: • Baseline subtracted curve fit • Apply fluorescence drift correction
	IAC	VIC/HEX	+	
Bio Molecular Systems MIC¹	<i>Salmonella</i> spp.	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please cyclor operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	
R-Biopharm RIDA®CYCLER¹	<i>Salmonella</i> spp.	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please cyclor operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	
Qiagen Rotor- Gene Q²	<i>Salmonella</i> spp.	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	

SureFast® Salmonella ONE (100 rxn)

Art. No. F5211

March 2026

Real-time PCR device	Detection	Detection channel	Quencher	Note
Roche LightCycler® 480 Instrument II ²	<i>Salmonella</i> spp.	465-510	+	
	IAC	533-580	+	
Roche cobas® z 480 Analyzer ²	<i>Salmonella</i> spp.	465-510	+	
	IAC	540-580	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Cultural Enrichment

Sample preparation of food and feed matrices should apply to ISO 6887 and ISO 6579 standards. For enrichment, a volume of 25 g or 25 ml sample material is homogenized with 225 ml of a non-selective enrichment broth (e. g. buffered peptone water). If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio. Enrichment is performed for $18 \text{ h} \pm 2 \text{ h}$ at $37^\circ\text{C} \pm 1^\circ\text{C}$. Alternatively, other suitable, validated enrichment procedures can be used.

To assess the process of bacterial growth, it is recommended to compare samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

2.1.2 DNA-preparation

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 500 µl of Lysis Buffer (**Code L**) into a 2.0 ml tube (not provided with the kit) and add 500 µl of the upper third from the enrichment to the Lysis Buffer (**Code L**).

Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 min at 95°C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 min.

Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction.

The lysate is ready-to-use for PCR. If it is not used immediately, it can be stored at $2\text{--}8^\circ\text{C}$ for up to 4h. If it is intended for a longer storage, transfer 100 µl of the lysate in a new tube (not provided with the kit) and store at -16 to -28°C .

2.1.3 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive Control.

For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control.

The reaction mix contains an internal amplification control (inhibition control) per reaction.

Reactions needed for the qualitative *Salmonella* spp. detection:

3 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative *Salmonella* spp. detection in enrichments:

5 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control,
1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

*** Description of the controls**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – use only Lysis buffer
- Positive Control: master-mix and within the kit's provided Positive Control
- Zero control: master-mix and sample prior enrichment
- Medium control: only medium – no sample

2.1.4 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Salmonella spp. DNA is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** for *Salmonella* spp., if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for *Salmonella* spp., if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC/HEX-channel) of the sample is **positive**.

If the sample DNA shows no amplification or sigmoidal amplification curve in the VIC/HEX channel, this indicates that inhibitors are present in the sample DNA that inhibit the PCR. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to centrifuge the lysate before the PCR or to repeat the extraction process. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Salmonella* spp. PCR assay.

Note: During the MicroVal study increased inhibition levels were observed for fresh herbs, cocoa, spices and aromatics as well as egg products. In these cases, please proceed as recommended above.

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
negative control	negative	$25 \leq C_p \leq 34$
Positive Control (FAM – <i>Salmonella</i> spp.)	$25 \leq C_p \leq 33$	not relevant

2.3 Cultural Confirmation

Positive PCR results should be confirmed by cultural detection methods (e.g. ISO 6579). Therefore 0,1 ml of the primary enrichment will be transferred to 10 ml in secondary enrichment (e.g. Rappaport-Vassiliadis) and processed according to ISO after incubation.

3 Validation

SureFast® Salmonella ONE was validated per the AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces and the ISO 16140-2:2016 [Microbiology of the food chain -- Method validation -- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method]. The reference method for the matrix study was ISO 6579:2002 (Microbiology of food and animal feeding stuffs -- Horizontal method for the detection of *Salmonella* spp.).

The kit is an extension of the SureFast® Salmonella PLUS method and has been AOAC RI validated for the following additional food matrices: raw ground beef, frozen poultry meat, raw milk, fresh spinach, pasteurized liquid whole egg, pet food pellets.

The MicroVal method comparison study according to ISO 16140-2:2016 was carried out using 25 g portions of sample material. Sample preparations were done according to ISO 6887 parts. In the sensitivity study a total of 391 samples from six different food categories were investigated. SureFast® Salmonella ONE has been MicroVal validated for the following food categories: poultry meat, meat products, dairy products, vegetables (excluding sprouts), egg products and feed.

4 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- A positive test result does not necessary indicate the presence of viable organism. It is indicative for the presence of the target DNA (*Salmonella* DNA).
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.

5 Further Information

5.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- AOAC and MicroVal certificates (Download: www.congen.de/en/quality-standards/)
- Validation Report upon request

5.2 Technical Support

For further questions please send an e-mail to info@congen.de.