
CONGEN

SureFast®
Listeria 3plex ONE

Art. No. F5217
100 rxn

User Manual

July 2023

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	Kit components and storage	10
1.4	Additionally required equipment and materials	10
1.5	Precautions for users	10
1.6	Setup.....	11
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Cultural Enrichment	12
2.1.2	DNA-preparation.....	12
2.1.3	Preparation of the master-mix	13
2.1.4	Preparation of the real-time PCR-mix	13
2.2	Interpretation of results	14
3	Further Information	14
3.1	Product Information	14
3.2	Technical Support.....	14

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1 General Information

1.1 Description

The SureFast® Listeria 3plex ONE can be applied for the fast and simple isolation and differentiation of *Listeria* spp. and *Listeria monocytogenes* in different types of food and is intended to be used by all kinds of food testing laboratories.

The method consists of the following four steps: 1) cultural enrichment, 2) DNA-preparation, 3) specific real-time PCR detection and 4) interpretation of results.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of fluorescence emissions at the channels **VIC/HEX**, **ROX** and **Cy5** at the same time. The internal technical verification of instruments was performed on Roche LightCycler® 480 II, Bio-Rad CFX96 Dx, Bio-Rad CFX96 Opus, R-Biopharm RIDA®CYCLER, Applied Biosystems 7500 Fast Dx and Agilent AriaDx.

1.2 Limit of Detection

The SureFast® Listeria 3plex ONE real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
L	Lysis Buffer	2 x 25 ml	Clear

Store all reagents at -20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

The Lysis Buffer should be stored at -20 to +8 °C. Reagents should be used until the expiry date indicated on the label.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.4 Additionally required equipment and materials

- real-time PCR instrument with **four** detection channels (510 nm, **580 nm, 610 nm** and **660 nm**)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- 2 ml tubes
- powder-free disposable gloves
- Vortex mixer
- centrifuge with a rotor for the reaction tubes
- Thermomixer/heating block (up to 95 °C)
- enrichment broth (e. g. Half Fraser)

1.5 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

SureFast® Listeria 3plex ONE (100 rxn)

Art. No. F5217

July 2023

1.6 Setup

	Blockcyler & R-Biopharm RIDA®CYCLER	Rotorcyler & LightCycler® 480 II
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaDx/Mx	IAC	VIC/HEX	+	
	<i>Listeria</i> spp.	ROX	+	
	<i>Listeria monocytogenes</i>	Cy5	+	
Applied Biosystems 7500 Fast Dx	IAC	VIC/HEX	None	Check the passive reference option ROX is none.
	<i>Listeria</i> spp.	ROX	None	
	<i>Listeria monocytogenes</i>	Cy5	None	
Bio-Rad CFX96/Dx/Opus	IAC	VIC/HEX	+	
	<i>Listeria</i> spp.	ROX	+	
	<i>Listeria monocytogenes</i>	Cy5	+	
R-Biopharm RIDA®CYCLER	IAC	yellow	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cyler operating instructions, section 12.1.2 Cycling analysis parameter.
	<i>Listeria</i> spp.	orange	+	
	<i>Listeria monocytogenes</i>	red	+	
Roche LightCycler® 480 II	IAC	533-580	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required. Evaluate the ROX channel with Fit Point.
	<i>Listeria</i> spp.	533-610	+	
	<i>Listeria monocytogenes</i>	618-660	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Cultural Enrichment

Sample preparation of food and feed matrices should apply to ISO 6887. For enrichment, a volume of 25 g or 25 ml sample material is homogenized with 225 ml of enrichment broth (Half Fraser). If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio. Enrichment is performed at 37 °C for 26 - 28 h. If SureFast® PREP Bacteria (F1021) is used for DNA-extraction enrichment has to be performed at 37 °C for 18 - 20 h. Alternatively, other suitable, validated enrichment procedures can be used.

To assess the process of bacterial growth, it is recommended to compare samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

2.1.2 DNA-preparation

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Pipette 500 µl Lysis Buffer (**Code L**) in a 2.0 ml Tube (not provided with the kit).

Afterwards transfer 200 µl of the upper third from the enrichment culture and add it to the 500 µl Lysis Buffer.

Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction. If the particles do not sediment, the lysate should be centrifuged for 1 minute at 12,000 rpm. Afterwards transfer 100 µl of the supernatant into a new tube (not provided with the kit).

The lysate is ready-to-use for PCR. If the lysate is not used within four hours in the PCR or is intended for a longer storage, transfer 100 µl of the lysate in a new tube (not provided with the kit) and store at -20 °C.

2.1.3 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control and positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative *Listeria spp. and Listeria monocytogenes* detection in enrichments:

5 reactions for controls (1x no-template control, 1x extraction control, 1x positive control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.4 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Listeria spp. DNA is detected in the ROX-channel and *Listeria monocytogenes* DNA is detected in the Cy5-channel. In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value ≤ 2 compared to the negative control. If the sample DNA in the VIC/HEX-channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Listeria* spp. and *Listeria monocytogenes* PCR assay.

Result in the respective channel			Interpretation
ROX channel <i>Listeria</i> spp.	Cy5 channel <i>Listeria</i> <i>monocytogenes</i>	VIC/HEX channel IAC	
positive	negative	positive/negative	<i>Listeria</i> spp.- DNA detected
positive	positive	positive/negative	<i>Listeria monocytogenes</i> -DNA detected
negative	negative	positive	negative, <i>Listeria</i> - DNA is not detected
negative	positive	positive/negative	invalid
negative	negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.