
CONGEN

# SureFast® Listeria 3plex ONE

Art. No. F5217  
100 rxn

## User Manual

September 2025

 **Inhalt**

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## **Content**

|       |                                                     |    |
|-------|-----------------------------------------------------|----|
| 1     | General Information .....                           | 13 |
| 1.1   | Description .....                                   | 13 |
| 1.2   | Limit of Detection.....                             | 13 |
| 1.3   | Kit components and storage.....                     | 14 |
| 1.4   | Additionally required equipment and materials ..... | 14 |
| 1.5   | Precautions for users.....                          | 14 |
| 1.6   | Setup .....                                         | 15 |
| 1.7   | Detection channel Set-up.....                       | 15 |
| 2     | Qualitative Analysis .....                          | 16 |
| 2.1   | Protocol.....                                       | 16 |
| 2.1.1 | Cultural Enrichment .....                           | 16 |
| 2.1.2 | DNA-preparation.....                                | 17 |
| 2.1.3 | Preparation of the master-mix .....                 | 18 |
| 2.1.4 | Preparation of the real-time PCR-mix .....          | 18 |
| 2.2   | Interpretation of results.....                      | 19 |
| 2.3   | Cultural Confirmation.....                          | 19 |
| 3     | Validation .....                                    | 20 |
| 4     | Limitations of the method .....                     | 22 |
| 5     | Further Information.....                            | 22 |
| 5.1   | Product Information .....                           | 22 |
| 5.2   | Technical Support .....                             | 22 |

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*Listeria sensu stricto* species are the primary targets for food safety. This group includes *L. monocytogenes* and the species that are genotypically and phenotypically similar to *L. monocytogenes*. Conversely, the *Listeria sensu lato* are more distantly related to *L. monocytogenes*; many lack the hallmark characteristics attributed to the *Listeria* genus. It is expected that all *Listeria sensu stricto* are detected, while some variation in detection may occur with the *Listeria sensu lato* species.

The method consists of the following four steps: 1) cultural enrichment, 2) DNA preparation, 3) specific real-time PCR detection and 4) interpretation of results. Step 2) DNA preparation can also be performed using SureFast® PREP Bacteria (F1021).

SureFast® Listeria 3plex ONE has been validated and certified both individually and in the context of an extraction using SureFast® PREP Bacteria (F1021) by MicroVal (License-No. 2023LR114) and AOAC (License-No. 062501).

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR-inhibiting substances, the signal of the amplification control can be affected, or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition, spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels **VIC/HEX**, **ROX** and **Cy5** at the same time.

The external technical validations of instruments (AOAC, MicroVal) were performed on Bio-Rad CFX96 Deep Well DX, Bio-Rad CFX96 Opus Deep Well and R-Biopharm RIDA®CYCLER.

The internal technical verification of instruments (not part of the PTM validation) was additionally performed on Bio-Rad CFX96 DX, Roche LightCycler® 480 Instrument II, Qiagen Rotor-Gene Q and Agilent AriaDx.

### 1.2 Limit of Detection

The SureFast® Listeria 3plex ONE real-time PCR has a limit of detection of  $\leq 5$  DNA copies per reaction.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

**1.3 Kit components and storage**

| Kit Code | Reagent          | Amount      | Lid Color  |
|----------|------------------|-------------|------------|
| 1        | Reaction Mix     | 2 x 1050 µl | Yellow     |
| 2        | Taq Polymerase   | 1 x 80 µl   | Dark Red   |
| 3        | Positive Control | 1 x 200 µl  | Light Blue |
| L        | Lysis Buffer     | 2 x 25 ml   | Clear      |

**Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

The Lysis Buffer should be stored at -28 to +8°C. Reagents should be used until the expiration date indicated on the label. Reagents can sustain up to 20 freeze/thaw cycles without influencing the assay performance.

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

**1.4 Additionally required equipment and materials**

- real-time PCR instrument with **four** detection channels (510 nm, **580 nm, 610 nm and 660nm**)
- SureFast® PREP Bacteria (Art. No. F1021), optional
- real-time PCR consumables (plates, tubes, foils, caps)
- incubator (up to 37°C)
- pipettes with filter tips, for volumes from 0.7 µL – 1000 µL
- microcentrifuge tubes - 2 ml tubes with lock cap
- powder-free disposable gloves
- vortex mixer
- microcentrifuge (up to 12,000 rpm) with a rotor for the reaction tubes
- thermomixer/heating block (up 99°C)
- enrichment broth (i.e. Half Fraser)
- suitable equipment for sample comminution and homogenization
- balance and spatula for weighing the samples
- autoclave
- biological safety cabinet

**1.5 Precautions for users**

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious, as well as all reagents and materials being exposed to the samples and must be handled according to the national safety regulations.
- Do not use the kit after the expiration date.
- Due to the significantly increased risk of infection and mortality, direct contact with *L. monocytogenes* samples is strongly discouraged for pregnant women and people with weakened immune system.
- Dispose of contents to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulations.

## 1.6 Setup

|                                        | Blockcycler <sup>1</sup> &<br>R-Biopharm RIDA®CYCLER <sup>1</sup> | Rotorcycler <sup>2</sup> &<br>LightCycler® 480<br>Instrument II <sup>2</sup> |
|----------------------------------------|-------------------------------------------------------------------|------------------------------------------------------------------------------|
| Initial Denaturation (HOLD)            | 1 min, 95°C                                                       | 1 min, 95°C                                                                  |
| Cycles                                 | 45                                                                | 45                                                                           |
| Denaturation                           | 15 sec, 95°C                                                      | 10 sec, 95°C                                                                 |
| Annealing/Extension (CYCLE)            | 30 sec, 60°C                                                      | 15 sec, 60°C                                                                 |
| Temperature Transition Rate/ Ramp Rate | Maximum                                                           | Maximum                                                                      |

## 1.7 Detection channel Set-up

For additional details see cycler operating instructions for real-time PCR device

| Real-time PCR device                                        | Detection                     | Detection channel | Quencher | Note                                                                                                                                                                                                           |
|-------------------------------------------------------------|-------------------------------|-------------------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Agilent AriaDx/Mx</b> <sup>1</sup>                       | IAC                           | VIC/HEX           | +        |                                                                                                                                                                                                                |
|                                                             | <i>Listeria</i> spp.          | ROX               | +        |                                                                                                                                                                                                                |
|                                                             | <i>Listeria monocytogenes</i> | Cy5               | +        |                                                                                                                                                                                                                |
| <b>Qiagen Rotor-Gene Q</b> <sup>2</sup>                     | IAC                           | yellow            | +        | <b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.                                                                                     |
|                                                             | <i>Listeria</i> spp.          | orange            | +        |                                                                                                                                                                                                                |
|                                                             | <i>Listeria monocytogenes</i> | red               | +        |                                                                                                                                                                                                                |
| <b>Bio-Rad CFX96/Dx/Opus (incl. Deep Well)</b> <sup>1</sup> | IAC                           | VIC/HEX           | +        | Baseline Settings:<br><ul style="list-style-type: none"> <li>Baseline subtracted curve fit</li> <li>Apply fluorescence drift correction</li> </ul>                                                             |
|                                                             | <i>Listeria</i> spp.          | ROX               | +        |                                                                                                                                                                                                                |
|                                                             | <i>Listeria monocytogenes</i> | Cy5               | +        |                                                                                                                                                                                                                |
| <b>R-Biopharm RIDA®CYCLER</b> <sup>1</sup>                  | IAC                           | yellow            | +        | <b>Ignore cycles before</b> , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter. |
|                                                             | <i>Listeria</i> spp.          | orange            | +        |                                                                                                                                                                                                                |
|                                                             | <i>Listeria monocytogenes</i> | red               | +        |                                                                                                                                                                                                                |
| <b>Roche LightCycler® 480 Instrument II</b> <sup>2</sup>    | IAC                           | 533-580           | +        | The SureCC Color Compensation Kit I (Art. No. F4009) is required. Evaluate the ROX channel with Fit Point.                                                                                                     |
|                                                             | <i>Listeria</i> spp.          | 533-610           | +        |                                                                                                                                                                                                                |
|                                                             | <i>Listeria monocytogenes</i> | 618-660           | +        |                                                                                                                                                                                                                |

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Cultural Enrichment

Sample preparation of food and feed matrices should apply to ISO 6887. Sampling for surface samples should be according to DIN EN ISO 18593:2018 - Microbiology of the food chain - Horizontal methods for surface sampling.

For enrichment, a volume of 25 g or 25 ml sample material is mixed with 225 ml of enrichment broth (Half Fraser, room temperature), homogenized for 30 seconds (e.g. in a Stomacher homogenizer) and incubated as indicated below.

If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio.

#### Enrichment conditions

| Sample type                                         | Enrichment                                                                                                                                                                                | Incubation                                                                                           |                                                                                                      |
|-----------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
|                                                     |                                                                                                                                                                                           | SureFast® Listeria 3plex ONE Lysis                                                                   | SureFast® PREP Bacteria                                                                              |
| Food and feed matrices, except dairy                | 25g sample + 225 ml Half Fraser Broth                                                                                                                                                     | 37°C ± 1°C<br>for 26-28 h                                                                            | 37°C ± 1°C<br>for 18-20 h                                                                            |
| Environmental samples                               | <ul style="list-style-type: none"><li>• 25ml Process Rinse Water + 225ml Half Fraser Broth</li><li>• swab* + 9 ml Half Fraser Broth</li><li>• sponge* + 90 ml Half Fraser Broth</li></ul> | 37°C ± 1°C<br>for 26-28 h                                                                            | 37°C ± 1°C<br>for 18-20 h                                                                            |
| Heat processed cheeses, raw milk and dairy products | 25g sample + 225 ml Half Fraser Broth<br><br>Following incubation, transfer of 0.1 ml of primary enrichment to 10 mL of Fraser Broth and further incubate                                 | 37°C ± 1°C<br>for 26-28 h<br>(Half Fraser Broth)<br>+<br>37°C ± 1°C<br>for 22-26 h<br>(Fraser Broth) | 37°C ± 1°C<br>for 18-20 h<br>(Half Fraser Broth)<br>+<br>37°C ± 1°C<br>for 22-26 h<br>(Fraser Broth) |

#### Note:

\*Swabs containing 1 ml Hi-Cap neutralizing broth; Sponges containing 10 mL Hi-Cap neutralizing broth.

To assess the process of bacterial growth, it is recommended to compare samples before enrichment (zero control) and at the end of the culturing (bacterial growth at Cp difference > 3 to zero control).

## 2.1.2 DNA-preparation

The DNA preparation can be conducted with the materials included in this kit following the protocol below or with the SureFast® PREP Bacteria (F1021) Kit according to the kit instructions.

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Pipette 500 µl of Lysis Buffer (**Code L**) into a 2.0 ml reaction tube (not provided with the kit).

Afterwards, transfer 200 µl of the upper third from the enrichment culture to the Lysis Buffer (**Code L**).

Vortex the reaction tube briefly.

Incubate in a heating block for 10 minutes at 95°C without shaking.

Remove the reaction tube from the heating block and leave it at room temperature for 1 minute.

**Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction. If the particles do not sediment, the lysate should be centrifuged for 1 minute at 12,000 rpm. Afterwards transfer 100 µl of the supernatant into a new reaction tube (not provided with the kit).**

The lysate is ready-to-use for PCR. If it is not used immediately, it can be stored at 2-8°C for up to 4h. If it is intended for a longer storage, transfer 100 µl of the lysate in a new tube (not provided with the kit) and store at -20°C

**2.1.3 Preparation of the master-mix**

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Required control reactions for the specific PCR assay: negative control, extraction control, and Positive Control.

For the analysis of enrichments, additional controls are needed: zero control and medium control.

The reaction mix contains an internal amplification control (Inhibition Control) per reaction.

**Reactions needed for the qualitative *Listeria* spp. and *Listeria monocytogenes* detection:** (e.g. when using samples without enrichment a.o. picked colonies or swabs - not part of the PTM validation\*):

3 reactions for controls\*\* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

**Reactions needed for the qualitative *Listeria* spp. and *Listeria monocytogenes* detection in enrichments:**

5 reactions for controls\*\* (1x negative control, 1x extraction control, 1x Positive Control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

**Example for the calculation and preparation of 10 reactions:**

| Components of the master-mix | Amount per reaction | 10 reactions (with 10% excess) |
|------------------------------|---------------------|--------------------------------|
| Reaction Mix                 | 19.3 µl             | 212.3 µl                       |
| Taq Polymerase               | 0.7 µl              | 7.7 µl                         |
| <b>Total volume</b>          | <b>20 µl</b>        | <b>220 µl</b>                  |

**Mix each master-mix well and centrifuge shortly before use.**

\* After DNA extraction with SureFast® PREP Bacteria (F1021), according to kit instructions

**\*\* Description of the controls**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – use only Lysis Buffer
- Positive Control: master-mix and within the kit's provided Positive Control
- Zero control: master-mix and sample prior incubate enrichment (analysis is recommended to rule out false positives in samples with late Cps)
- Medium control: only medium - no sample (analysis recommended to rule out contamination)

**2.1.4 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates (depending on the thermocycler used) briefly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation shall be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions must show the correct results.

*Listeria* spp. is detected in the ROX-channel and *Listeria monocytogenes* DNA is detected in the Cy5-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive**.

If the sample DNA shows no amplification or sigmoidal amplification curve in the VIC/HEX channel, this indicates that inhibitors are present in the sample DNA that inhibit the PCR. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to centrifuge the lysate before the PCR or to repeat the extraction process. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Listeria* spp. and *Listeria monocytogenes* PCR assay.

| Result in the respective channel    |                                                        |                          | Interpretation                                                |
|-------------------------------------|--------------------------------------------------------|--------------------------|---------------------------------------------------------------|
| ROX channel<br><i>Listeria</i> spp. | Cy5 channel<br><i>Listeria</i><br><i>monocytogenes</i> | VIC/HEX channel<br>IAC   |                                                               |
| <b>positive</b>                     | negative                                               | <b>positive/negative</b> | <i>Listeria</i> spp. DNA detected                             |
| <b>positive</b>                     | <b>positive</b>                                        | <b>positive/negative</b> | <i>Listeria monocytogenes</i><br>DNA detected                 |
| negative                            | <b>positive</b>                                        | <b>positive</b>          | questionable* for <i>Listeria</i><br><i>monocytogenes</i> DNA |
| negative                            | <b>positive</b>                                        | <b>negative</b>          | invalid                                                       |
| negative                            | negative                                               | <b>positive</b>          | negative, <i>Listeria</i> DNA<br>is not detected              |
| negative                            | negative                                               | negative                 | invalid                                                       |

\*Result for *Listeria monocytogenes* is questionable and analysis for this sample shall be repeated. This can appear with low levels of target at the limit of detection.

## 2.3 Cultural Confirmation

Positive PCR results should be confirmed by cultural detection methods according to the ISO 11290-1:2017 reference standard.

### 3 Validation

#### MicroVal certification (License no. 2023LR114)

SureFast® Listeria 3plex ONE was MicroVal validated (License No. 2023LR114) for a broad range of food and environmental samples per ISO 16140-2:2016 [Microbiology of the food chain -- Method validation -- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method].

The reference method for the matrix study was ISO 11290-1:2017 (Microbiology of the food chain - Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp.).

#### Summary of Validation Claims for MicroVal certification 2023LR114

MicroVal claim for the detection of *Listeria* spp. and / or *Listeria monocytogenes*

| Category                                                       | Enrichment                                                                                                                           | Sample size                | SureFast® Listeria 3plex ONE Lysis                                                     | SureFast® PREP Bacteria                                                                |
|----------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|----------------------------|----------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| Heat-Processed Cheeses/ Raw Milk and Dairy Products            | 1:10 dilution<br>Half Fraser Broth.<br>Following incubation,<br>transfer of 0.1 ml of<br>primary enrichment to<br>10 ml Fraser Broth | 25 g                       | 37 ± 1°C<br>26-28 h (Half Fraser<br>Broth) +<br>37 ± 1°C<br>22- 26 H<br>(Fraser Broth) | 37 ± 1°C<br>18-20 h (Half Fraser<br>Broth) +<br>37 ± 1°C<br>22- 26 H<br>(Fraser Broth) |
| Raw and Ready-to-eat (RTE) Meat Products                       | 1:10 dilution<br>Half Fraser Broth<br><br>(25 g sample enriched<br>with 225 mL of Half<br>Fraser Broth)                              | 25 g                       | 37 ± 1°C<br>26-28 h                                                                    | 37 ± 1°C<br>18-20 h                                                                    |
| Ready- to- eat (RTE) Poultry Products                          |                                                                                                                                      |                            |                                                                                        |                                                                                        |
| Ready- to- eat (RTE), Ready- to-reheat (RTRH) Fishery Products |                                                                                                                                      |                            |                                                                                        |                                                                                        |
| Fresh Produce and Fruits                                       |                                                                                                                                      |                            |                                                                                        |                                                                                        |
| Multi-component Foods or Meal Components                       |                                                                                                                                      |                            |                                                                                        |                                                                                        |
| Environmental Samples                                          | 1:10 dilution<br>Half Fraser Broth<br><br>(1 sterile swab +<br>9 ml Half Fraser<br>1 sponge +<br>90 ml Half Fraser)                  | Sponge,<br>swab or<br>25 g |                                                                                        |                                                                                        |

Certification was performed on Bio-Rad CFX96 Deep Well DX (Maestro software version 2.0), Bio-Rad CFX96 Opus Deep Well (Maestro software version 2.0) and R-Biopharm RIDA®CYCLER (RIDA®CYCLER software version 1.2.1).

## AOAC PTM certification (License no. 062501)

SureFast® Listeria 3plex ONE was validated per the AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces and the ISO 16140-2:2016 [Microbiology of the food chain -- Method validation -- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method].

The reference method for the matrix study was ISO 11290-1:2017 (Microbiology of the food chain - Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp.).

## Summary of Validation Claims for AOAC certification 062501

AOAC claim for the detection of *Listeria* spp. and / or *Listeria monocytogenes*

| Matrix                                         | Enrichment                                                                                                | SureFast® Listeria 3plex ONE Lysis | SureFast® PREP Bacteria       | <i>L. spp</i> | <i>L. mono</i> |
|------------------------------------------------|-----------------------------------------------------------------------------------------------------------|------------------------------------|-------------------------------|---------------|----------------|
| Raw Milk                                       | 1:10 dilution                                                                                             | 37 ± 1°C                           | 37 ± 1°C                      | X             | X              |
| Brie Cheese                                    | Half Fraser Broth.                                                                                        | 26-28 h (Half Fraser Broth) +      | 18-20 h (Half Fraser Broth) + | X             | X              |
| Gorgonzola Cheese                              | Following                                                                                                 | 37 ± 1°C                           | 37 ± 1°C                      | X             | X              |
| Pasteurized Milk                               | incubation, transfer of 0.1 ml of primary                                                                 | 22- 26 h                           | 22- 26 h                      | X             | X              |
| Gouda Cheese                                   | enrichment to 10 ml                                                                                       | (Fraser Broth)                     | (Fraser Broth)                | X             | X              |
| Ice Cream                                      | Fraser Broth                                                                                              |                                    |                               | X             |                |
| Liquid Whole Egg                               | 1:10 dilution<br>Half Fraser Broth<br><br>(25 g / 25 ml sample enriched with 225 ml of Half Fraser Broth) | 37 ± 1°C<br>26-28 h                | 37 ± 1°C<br>18-20 h           | X             | X              |
| Frankfurter Sausage                            |                                                                                                           |                                    |                               | X             | X              |
| Deli Ham                                       |                                                                                                           |                                    |                               | X             | X              |
| Smoked Deli Turkey                             |                                                                                                           |                                    |                               | X             | X              |
| Smoked Salmon                                  |                                                                                                           |                                    |                               | X             |                |
| Frozen cooked Shrimp                           |                                                                                                           |                                    |                               | X             | X              |
| Bagged Salad                                   |                                                                                                           |                                    |                               | X             | X              |
| Ground Soy Vegan Meat Substitute               |                                                                                                           |                                    |                               | X             | X              |
| Canned Enoki Mushrooms                         |                                                                                                           |                                    |                               | X             | X              |
| Guacamole                                      |                                                                                                           |                                    |                               | X             |                |
| Process Rinse Water                            |                                                                                                           |                                    |                               | X             |                |
| Stainless Steel Environmental Sample (Swabs)   | 1:10 dilution<br>Half Fraser Broth                                                                        |                                    |                               | X             | X              |
| Stainless Steel Environmental Sample (Sponges) | (1 sterile swab + 9 ml Half Fraser<br>1 sponge + 90 ml Half Fraser)                                       |                                    |                               | X             | X              |

Certification was performed on Bio-Rad CFX96 Deep Well DX (Maestro software Version 2.0), Bio-Rad CFX96 Opus Deep Well (Maestro software Version 2.0) and R-Biopharm RIDA®CYCLER (RIDA®CYCLER software version 1.2.1).

## **4 Limitations of the method**

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- For samples at the LoD, a positive signal may occur in the *Listeria monocytogenes* channel and not in the *Listeria* spp. channel. Analysis should be repeated.
- The following *Listeria* spp. sensu lato strains, which are primarily known as environmental bacteria and not as pathogens, are not detected: *Listeria grandensis*, *Listeria floridensis*, *Listeria rocourtiae*, *Listeria cornellensis*, *Listeria aquatica*, and *Listeria portnoyi*.
- A positive test result does not necessary indicate the presence of viable organisms. It is indicative for the presence of the target DNA (*Listeria*-DNA).
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.

## **5 Further Information**

### **5.1 Product Information**

- Detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/downloads](http://www.congen.de/en/downloads))
- Product-related documents (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))
- AOAC and MicroVal certificates (Download: [www.congen.de/en/quality-standards/](http://www.congen.de/en/quality-standards/))
- Validation Report upon request

### **5.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).