
CONGEN

**SureFast®
STEC 4plex ONE**

Art. No. F5265
100 rxn

User Manual

August 2023

 **Inhalt**

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Content

1	General Information	10
1.1	Description	10
1.2	Limit of Detection	10
1.3	Kit components and storage	11
1.4	Additionally required equipment and materials	11
1.5	Setup	11
1.6	Detection channel Set-up	12
2	Qualitative Analysis	13
2.1	Protocol	13
2.1.1	Cultural Enrichment	13
2.1.2	DNA-preparation	13
2.1.3	Preparation of the master-mix	14
2.1.4	Preparation of the real-time PCR-mix	14
2.2	Interpretation of results	15
3	Further Information	16
3.1	Product Information	16
3.2	Technical Support	16

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1 General Information

1.1 Description

The SureFast® STEC 4plex ONE can be applied for the fast and simple isolation, detection and differentiation of specific DNA sequences of the *Escherichia coli* virulence factors *stx1* (subtype a-d), *stx2* (subtype a-g), *eae* and the *Escherichia coli* serotype O157 in different types of food.

The method consists of the following three steps: 1) cultural enrichment, 2) DNA-preparation, 3) specific real-time PCR detection and 4) interpretation of results.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX und Cy5 at the same time. The technical verification of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Bio-Rad CFX96, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Applied Biosystems 7500, Qiagen Rotor-Gene Q and Bio Molecular Systems MIC.

1.2 Limit of Detection

The SureFast® STEC 4plex ONE real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
L	Lysis Buffer	2 x 25 ml	Clear

Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

The Lysis Buffer should be stored at -20 to +8 °C. Reagents should be used until the expiry date indicated on the label.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.4 Additionally required equipment and materials

- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- 2 ml tubes
- powder-free disposable gloves
- Vortex mixer
- centrifuge with a rotor for the reaction tubes
- Thermomixer/heating block (up to 95 °C)
- enrichment broth (e. g. modified tryptone soy broth (mTSB) with Novobiocin)

1.5 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.6 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	<i>E. coli stx1/stx2</i>	FAM	+	
	IAC	HEX	+	
	<i>E. coli</i> O157	ROX	+	
	<i>E. coli eae</i>	Cy5	+	
Agilent AriaDx	<i>E. coli stx1/stx2</i>	FAM	+	
	IAC	HEX	+	
	<i>E. coli</i> O157	ROX	+	
	<i>E. coli eae</i>	Cy5	+	
Applied Biosystems 7500	<i>E. coli stx1/stx2</i>	FAM	None	Check the passive reference option ROX is none.
	IAC	HEX	None	
	<i>E. coli</i> O157	ROX	None	
	<i>E. coli eae</i>	Cy5	None	
Bio-Rad CFX96	<i>E. coli stx1/stx2</i>	FAM	+	
	IAC	VIC/HEX	+	
	<i>E. coli</i> O157	ROX	+	
	<i>E. coli eae</i>	Cy5	+	
R-Biopharm RIDA®CYCLER	<i>E. coli stx1/stx2</i>	green	+	
	IAC	yellow	+	
	<i>E. coli</i> O157	orange	+	
	<i>E. coli eae</i>	red	+	
Qiagen Rotor-Gene Q	<i>E. coli stx1/stx2</i>	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
	<i>E. coli</i> O157	orange	+	
	<i>E. coli eae</i>	red	+	
Roche LightCycler® 480 II	<i>E. coli stx1/stx2</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	533-580	+	
	<i>E. coli</i> O157	533-610	+	
	<i>E. coli eae</i>	618-660	+	
Roche cobas® z 480 Analyzer	<i>E. coli stx1/stx2</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	540-580	+	
	<i>E. coli</i> O157	540-610	+	
	<i>E. coli eae</i>	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Cultural Enrichment

Sample preparation of food and feed matrices should apply to ISO 6887 and ISO 13136 standards. For enrichment, a volume of 25 g or 25 ml sample material is homogenized with 225 ml of an enrichment broth (e. g. modified tryptone soy broth (mT5B) with Novobiocin). If test portions deviate from 25 g, the enrichment broth volume should be chosen in a 1:10 ratio. Enrichment is performed for 18 h - 24 h at 37 °C ± 1°C. Alternatively other suitable, validated enrichment procedures can be used.

To assess the process of bacterial growth, it is recommended to compare samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

2.1.2 DNA-preparation

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 200 µl of the upper third from the enrichment culture into a 2.0 ml tube (not provided with the kit) and add 500 µl of Lysis Buffer (**Code L**) to this sample.

Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction. If the particles do not sediment, the lysate should be centrifuged for 1 minute at 12.000 rpm. Afterwards transfer 100 µl of the supernatant into a new tube (not provided with the kit).

The lysate is ready-to-use for PCR. If the lysate is not used within four hours in the PCR or is intended for a longer storage, transfer 100 µl of the lysate in a new tube (not provided with the kit) and store at -20 °C.

2.1.3 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative STEC detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative STEC detection in enrichments:

5 reactions for controls (1x no-template control, 1x extraction control, 1x positive control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.4 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Escherichia coli stx1/stx2 DNA is detected in the FAM-channel, *Escherichia coli* O157 DNA is detected in the ROX-channel and *Escherichia coli eae* DNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value ≤ 2 compared to the negative control. If the sample DNA in the VIC/HEX channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Escherichia coli stx1/stx2*, *Escherichia coli* O157 or *Escherichia coli eae* PCR assay.

Result in the respective channel				Interpretation
FAM channel <i>Escherichia coli stx1/stx2</i>	ROX channel <i>Escherichia coli</i> O157	Cy5 channel <i>Escherichia coli eae</i>	VIC/HEX channel IAC	
positive	negative	negative	positive/negative	<i>Escherichia coli stx1/stx2</i> DNA detected
positive	positive	negative	positive/negative	<i>Escherichia coli</i> O157 DNA detected
negative	negative	positive	positive/negative	<i>Escherichia coli eae</i> DNA detected
negative	negative	negative	positive	Negative, <i>Escherichia coli stx1/stx2</i> , O157, <i>eae</i> – DNA is not detected
negative	negative	negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.