**CONGEN**

**SureFast<sup>®</sup>**  
**Pseudomonas aeruginosa PLUS**

Art. No. F5503  
100 rxn

**User Manual**

**December 2024**

 **Inhalt**

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## **Content**

|       |   |    |
|-------|---|----|
| 1     | General Information .....                           | 9  |
| 1.1   | Description .....                                   | 9  |
| 1.2   | Limit of Detection .....                            | 9  |
| 1.3   | DNA-preparation .....                               | 10 |
| 1.4   | Kit components and storage .....                    | 10 |
| 1.5   | Additionally required equipment and materials ..... | 10 |
| 1.6   | Precautions for users .....                         | 10 |
| 1.7   | Setup.....  | 11 |
| 1.8   | Detection channel Set-up .....                      | 11 |
| 2     | Qualitative Analysis .....                          | 13 |
| 2.1   | Protocol .....                                      | 13 |
| 2.1.1 | Preparation of the master-mix .....                 | 13 |
| 2.1.2 | Preparation of the real-time PCR-mix .....          | 13 |
| 2.2   | Interpretation of results .....                     | 14 |
| 3     | Limitations of the method .....                     | 14 |
| 4     | Further Information .....                           | 14 |
| 4.1   | Product Information .....                           | 14 |
| 4.2   | Technical Support .....                             | 14 |

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## 1 General Information

### 1.1 Description

The SureFast® *Pseudomonas aeruginosa* PLUS is a real-time PCR for the direct, qualitative detection of a specific DNA sequence of *Pseudomonas aeruginosa* in enriched food.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical verification of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, LTF MyGo Pro, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER, Roche LightCycler® 480 II, and Roche cobas® z 480 Analyzer.

### 1.2 Limit of Detection

The SureFast® *Pseudomonas aeruginosa* PLUS real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria (Art. No. F1021) is recommended.

To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

### 1.4 Kit components and storage

| Kit Code | Reagent          | Amount      | Lid Color  |
|----------|------------------|-------------|------------|
| 1        | Reaction Mix     | 2 x 1050 µl | Yellow     |
| 2        | Taq Polymerase   | 1 x 80 µl   | Dark Red   |
| 3        | Positive Control | 1 x 200 µl  | Light Blue |

**Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Bacteria Art. No. F1021)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

### 1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

**1.7 Setup**

|  | <b>Blockcycler &amp;<br/>R-Biopharm RIDA®CYCLER</b> | <b>Rotorcycler &amp;<br/>LightCycler® 480 II &amp;<br/>LTF MyGo Pro</b> |
|--|---|---|
| Initial Denaturation (HOLD)            | 1 min, 95°C   | 1 min, 95°C   |
| Cycles                                 | 45  | 45  |
| Denaturation                           | 15 sec, 95°C  | 10 sec, 95°C  |
| Annealing/Extension (CYCLE)            | 30 sec, 60°C  | 15 sec, 60°C  |
| Temperature Transition Rate/ Ramp Rate | Maximum   | Maximum   |

**1.8 Detection channel Set-up**

| <b>Real-time PCR device</b>    | <b>Detection</b>              | <b>Detection channel</b> | <b>Quencher</b> | <b>Note</b>  |
|--------------------------------|-------------------------------|--------------------------|-----------------|--|
| <b>Agilent Mx3005P</b>         | <i>Pseudomonas aeruginosa</i> | FAM                      | +               |  |
|                                | IAC                           | HEX                      | +               |  |
| <b>Agilent AriaDx /Mx</b>      | <i>Pseudomonas aeruginosa</i> | FAM                      | +               |  |
|                                | IAC                           | HEX                      | +               |  |
| <b>Applied Biosystems 7500</b> | <i>Pseudomonas aeruginosa</i> | FAM                      | None            | Check the passive reference option ROX is none.  |
|                                | IAC                           | VIC                      | None            |  |
| <b>Bio-Rad CFX96/Dx/Opus</b>   | <i>Pseudomonas aeruginosa</i> | FAM                      | +               |  |
|                                | IAC                           | VIC/HEX                  | +               |  |
| <b>R-Biopharm RIDA®CYCLER</b>  | <i>Pseudomonas aeruginosa</i> | green                    | +               | <b>Ignore cycles before</b> , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter. |
|                                | IAC                           | yellow                   | +               |  |
| <b>LTF MyGo Pro</b>            | <i>Pseudomonas aeruginosa</i> | FAM                      | +               |  |
|                                | IAC                           | VIC                      | +               |  |
| <b>Qiagen Rotor-Gene Q</b>     | <i>Pseudomonas aeruginosa</i> | green                    | +               | <b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.   |
|                                | IAC                           | yellow                   | +               |  |

# SureFast® *Pseudomonas aeruginosa* PLUS (100 rxn)

Art. No. F5503

December 2024

| Real-time PCR device        | Detection                     | Detection channel | Quencher | Note |
|-----------------------------|-------------------------------|-------------------|----------|------|
| Roche LightCycler® 480 II   | <i>Pseudomonas aeruginosa</i> | 465-510           | +        |      |
|                             | IAC                           | 533-580           | +        |      |
| Roche cobas® z 480 Analyzer | <i>Pseudomonas aeruginosa</i> | 465-510           | +        |      |
|                             | IAC                           | 540-580           | +        |      |

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, Positive Control.

For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control.

The reaction mix contains an internal amplification control (inhibition control) per reaction.

#### Reactions needed for the qualitative *Pseudomonas aeruginosa* detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

#### Reactions needed for the qualitative *Pseudomonas aeruginosa* detection in enrichments:

5 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### Example for the calculation and preparation of 10 reactions:

| Components of the master-mix | Amount per reaction | 10 reactions (with 10% excess) |
|------------------------------|---------------------|--------------------------------|
| Reaction Mix                 | 19.3 µl             | 212.3 µl                       |
| Taq Polymerase               | 0.7 µl              | 7.7 µl                         |
| <b>Total volume</b>          | <b>20 µl</b>        | <b>220 µl</b>                  |

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

*Pseudomonas aeruginosa* DNA is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** for *Pseudomonas aeruginosa*, if the sample DNA shows amplification in the FAM-channel. High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for *Pseudomonas aeruginosa*, if the sample DNA shows no amplification in the FAM-channel and if the internal control (VIC/HEX-channel) of the sample is positive.

If the sample DNA in the VIC/HEX-Channel shows no amplification or an irregular amplification curve, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Pseudomonas aeruginosa* PCR assay.

## 3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- A positive test result does not necessary indicate the presence of viable organism. It is indicative for the presence of the target DNA (*Pseudomonas aeruginosa* DNA).

## 4 Further Information

### 4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/downloads](http://www.congen.de/en/downloads))
- Product-related documents (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))
- Validation Report upon request

### 4.2 Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).