
**CONGEN**

# **SureFast<sup>®</sup> Fecal Screen 4plex**

Art. No. F 5504  
100 rxn

## **User Manual**

**March 2024**

## Inhalt

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## Content

1	General Information .....	9
1.1	Description .....	9
1.2	Limit of Detection .....	9
1.3	DNA-preparation .....	10
1.4	Kit components and storage .....	10
1.5	Additionally required equipment and materials .....	10
1.6	Precautions for users .....	10
1.7	Setup.....	11
1.8	Detection channel Set-up .....	11
2	Qualitative Analysis .....	12
2.1	Protocol .....	12
2.1.1	Preparation of the master-mix .....	12
2.1.2	Preparation of the real-time PCR-mix .....	12
2.2	Interpretation of results .....	13
3	Limitations of the method .....	14
4	Further Information .....	14
4.1	Product Information .....	14
4.2	Technical Support .....	14

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## 1 General Information

### 1.1 Description

The SureFast® Fecal Screen 4plex is a real-time PCR for the direct, qualitative detection and differentiation of specific DNA sequences of Enterobacteriaceae, Enterococci, *Escherichia coli* in water.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride).

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical verification of instruments was performed on Roche LightCycler® 480 II, Qiagen Rotor-Gene Q, Bio-Rad CFX96Dx, R-Biopharm RIDA®CYCLER and Agilent AriaDx.

### 1.2 Limit of Detection

The SureFast® Fecal Screen 4plex real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

#### **Note: Note:**

Inconsistent mixing ratios\* may cause a loss of sensitivity in the low concentration channel in mixed samples especially with high amplicon concentrations (Cp value < 20).

\* e.g. 99.9 % Enterobacteriaceae and 0.1 % *Escherichia coli*

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Aqua (Art. Nr. F1023) is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at –28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Aqua (Art. Nr. F1023))
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- microcentrifuge with a rotor for the reaction tubes

### 1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

## 1.7 Setup

	<b>Blockcyler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcyler &amp; LightCycler® 480 II</b>
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

## 1.8 Detection channel Set-up

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent AriaDx /Mx</b>	Enterobacteriaceae	FAM	+	
	IAC	HEX	+	
	Enterococci	ROX	+	
	<i>E. coli</i>	Cy5	+	
<b>Bio-Rad CFX96/Dx/Opus</b>	Enterobacteriaceae	FAM	+	
	IAC	VIC/HEX	+	
	Enterococci	ROX	+	
	<i>E. coli</i>	Cy5	+	
<b>R-Biopharm RIDA®CYCLER</b>	Enterobacteriaceae	green	+	
	IAC	yellow	+	
	Enterococci	orange	+	
	<i>E. coli</i>	red	+	
<b>Qiagen Rotor-Gene Q</b>	Enterobacteriaceae	green	+	<b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
	Enterococci	orange	+	
	<i>E. coli</i>	red	+	
<b>Roche LightCycler® 480 II</b>	Enterobacteriaceae	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	533-580	+	
	Enterococci	533-610	+	
	<i>E. coli</i>	618-660	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (Inhibitioncontrol) per reaction.

#### Reactions needed for the qualitative Enterobacteriaceae-, Enterococcen- and *E. coli*- detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

**2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Enterobacteriaceae DNA is detected in the FAM-channel, Enterococci DNA is detected in the ROX-channel and *E. coli* DNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-value  $\leq 2$  compared to the negative control. If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific Enterobacteriaceae, Enterococci or *E. coli* PCR assay.

Result in the respective channel				Interpretation
FAM channel Enterobacteriaceae	ROX channel Enterococci	Cy5 channel <i>E. coli</i>	VIC/HEX channel IAC	
<b>positive</b>	negative	negative	<b>positive/negativ</b>	Enterobacteriaceae-DNA detected
negative	<b>positive</b>	negative	<b>positive/negativ</b>	Enterococci-DNA detected
<b>positive</b>	negative	<b>positive</b>	<b>positive/negativ</b>	<i>E. coli</i> DNA detected
negative	negative	negative	<b>positive</b>	negative no target DNA detected
negative	negative	negative	negative	invalid

### **3 Limitations of the method**

- The presence of PCR inhibitors may cause invalid results.
- A positive test result does not necessary indicate the presence of viable organism. It is indicative for the presence of the target DNA (Enterobacteriaceae DNA).
- Extremely low levels of target below the limit of detection (LoD) may bedetected, but results may not be reproducible.
- 100 % Cross reactivity was observed with DNA extracts from *Aeromonas* and *Aggregatibacter* species .

### **4 Further Information**

#### **4.1 Product Information**

- Detailed information about setup of several real-time PCR devices  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request

#### **4.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).