RIDASCREEN® SET Total

REF R4105

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RIDASCREEN® SET Total

Brief information

RIDASCREEN® SET Total (Art. No. R4105) is an enzyme immunoassay for combined detection of *Staphylococcus aureus* enterotoxins A, B, C, D and E in fluid and solid foods as well as in bacterial cultures.

All reagents required for the enzyme immunoassay are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations. A microtiter plate spectrophotometer is required for determination.

Sample preparation: Extraction according to a simplified method for preparation of samples (only conditionally applicable, see chapter 9. "Sample preparation") by homogenization with buffer and sample preparation

Extraction according to the "Official European Preparation Method"
(dialysis concentration method (see chapter 9.5))

Time requirement: Sample preparation for 10 samplesapprox. 1 h
(simplified preparation)

Sample preparation for 10 samplesapprox. 19 h
(official method, dialysis overnight)

Test implementation (incubation time) 2 h 45 min

Limit of detection

Simplified preparation: Liquid samples 0.25 ng/mL toxin
Solid samples 0.375 ng/g toxin
Supernatants of bacterial cultures 0.25 ng/mL toxin

Dialysis concentration: Liquid samples 0.05 ng/mL toxin
Solid samples 0.05 ng/g toxin

Further information is contained in the validation report.

The specificity of the RIDASCREEN® SET Total test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website <https://food.r-biopharm.com/media/technical-guides/>.

Related products

RIDASCREEN® SET A,B,C,D,E (Art. No. R4101)

1. Intended use

RIDASCREEN® SET Total is a sandwich enzyme immunoassay for the combined detection of Staphylococcus enterotoxins (SET) A, B, C, D and E in fluid and solid foods as well as in bacterial cultures.

2. General information

Besides Salmonella, *Staphylococcus aureus* enterotoxins are the main causative agents of food poisoning. Among the strains of *Staphylococcus aureus*, other Staphylococci species as *S. hyicus* and *S. intermedius* are able to produce one or more heat stable proteins that behave like enterotoxins. Generally, it is assumed that a population of 5×10^5 cells of enterotoxin-producing *Staphylococcus aureus* strains per gram of food is required to lead to an intoxication. However, other studies showed that only 100 - 200 ng of Staphylococcus enterotoxins can lead to symptoms of food poisoning. SET intoxications have been frequently associated with pasta, finished meat products, ham, pies, chicken meat products, fish, fish products, milk, milk products, ice cream, egg products, salads, pastries and cake fillings as well as preparations from these food products. The enterotoxins of the serological groups A, B, C, D and E are very significant.

3. Test principle

RIDASCREEN® SET Total is a reliable test for the detection of the most important toxins (A, B, C, D and E) of *S. aureus*. The surface of the microtiter plate is coated with specific antibodies that selectively bind the enterotoxins contained in the sample. Sample components not bound by the antibodies are then removed in a washing step. Adding marked enterotoxin-specific antibodies as well as enzyme-marked detector molecules forms a sandwich complex (antibody-antigen-antibody-complex). After substrate/chromogen is added, the bound enzyme conjugate converts the chromogen into a blue end product. This color change points to the presence of enterotoxins in the samples. The results can be read photometrically. After addition of the stop solution, which leads to a color change from blue to yellow, the measurement can be made in a microtiter plate spectrophotometer.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including positive and negative control). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Positive control	Red	Ready to use	Red stained	2 mL
Negative control	White	Ready to use	Colorless	2 mL
Wash buffer	Brown	Concentrate	10x	100 mL
Conjugate 1	Red	Ready to use	Green stained	11 mL
Conjugate 2	Black	Ready to use	Blue stained	11 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use	Red stained	13 mL
Stop solution	Yellow	Ready to use		14 mL

5. Reagents required but not provided

5.1 Equipment

- Analytical balance and weighing vessels
- Mixer or equivalent for sample homogenization. For food matrices like whey powder and others that are difficult to homogenize, an ultra turrax is highly recommended to obtain a homogeneous sample.
- 50 mL tubes
- Optional: sterile filter
- 100 µL micropipette

- Incubator 35 - 37 °C (95 - 98.6 °F)
- Multichannel pipette or microwell plate washer
- Microwell plate spectrophotometer (450 / 620 ± 10 nm)

Additional equipment for sample preparation according to the European official screening method (EOSM) Version 5, September 2010:

General remark: it is strictly recommended to use only laboratory vessels (e.g. funnels, beakers, test tubes) of laboratory glass or polypropylene to avoid the absorption of toxins.

- Shaker for beakers (room temperature)
- Refrigerated centrifuge (4 °C), 3130 - 10 000 g, centrifuge tubes
- Dialysis membrane, MWCO: 6 - 8 kD, flat width: 23 ± 2 mm (e.g., Spectra / Por[®]1, ref: 132650, Spectrum)
- Closures, sealing width = 35 mm (e.g. Spectra / Por[®], ref: 132736, Spectrum)
- pH-meter
- 50 mL tubes
- Funnel
- Glass wool
- Laboratory glass trough
- Refrigerator (5 ± 3 °C / 41 ± 5.4 °F) and freezer (≤ -18 °C / ≤ -0.4 °F)
- Vortexer
- Microwell plate spectrophotometer (450 / 620 ± 10 nm)

5.2 Reagents

- PBS buffer, pH 7.4 (0.55 g NaH₂PO₄ x H₂O + 2.85 g Na₂HPO₄ x 2 H₂O + 8.7 g NaCl fill up to 1000 mL with distilled water)
- Distilled water (optional: sterile water)
- n-Heptane (for samples with high fat content)
- Brain heart infusion (BHI) for the pre-enrichment of potentially toxin-forming Staphylococci strains. Ask your local producer / distributor of culture media for the supply with BHI.

Additional reagents for the European official screening method (EOSM) Version 5, September 2010:

- Hydrochloric acid (5N and 1N)
- Sodium hydroxide (5N and 1N)
- Polyethylene glycol 20000 (PEG), quality for synthesis

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to addition in the wells
- Extinction less than or equal to 1.0 ($E_{450/620 \text{ nm}} \leq 1.0$) for the positive control as well as an extinction greater than or equal to 0.2 ($E_{450/620 \text{ nm}} \geq 0.2$) for the negative control

9. Sample preparation

The samples have to be stored at 5 ± 3 °C (41 ± 5.4 °F) until extraction. The samples should be completely thawed before the extraction step.

Important notes: To ensure good precipitation of solids and complete phase separation via centrifugation, increase the centrifugation speed and / or time as needed.

Food proteins may interact with toxins or antibodies of the test kit, in which case they can significantly affect the detection of the toxins. Therefore it is strongly recommended to prepare food samples according to the “*Official European Screening Method of the European Reference Laboratory for Coagulase-positive Staphylococci*” (see chapter 9.5).

The simplified methods for sample preparation described under chapter 9.1 to 9.3 may not be suitable for extracting SET from difficult matrices (e.g. fish, chocolate, pickled vegetables). For these foods, particularly low recovery rates and / or non-specific bindings of matrix proteins to the test antibodies have been observed.

A list of SET recovery rates determined for a wide range of different food matrices is available on request.

9.1 Simplified preparation of milk

- Centrifuge milk samples (10 - 25 mL), especially raw milk samples, in cool condition: 10 min / 3500 g / 10 °C (50 °F)
(if no refrigerated centrifuge is available, pre-cooling of samples is necessary)
- Thoroughly remove the upper cream layer
- Use 100 µL per well in the assay

9.2 Simplified preparation of pasta and rice (cooked), meat, ice cream, processed foods and other foods with a fat content of less than 40 %

- Grind 10 - 25 g of the sample and homogenize with 1.5 mL of PBS buffer (pH 7.4) per g sample (e.g. 10 g sample + 15 mL buffer)
- Shake for 15 min

- Centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- If necessary, remove upper fat layer
- Use 100 µL per well in the assay

9.3 Simplified preparation of foods with a fat content of more than 40 %

- Grind 10 - 25 g of the sample and homogenize with 1.5 mL of PBS buffer (pH 7.4) per g sample (e.g. 10 g sample + 15 ml buffer)
- Shake for 15 min
- Centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- Transfer aqueous phase to another centrifugal vial, add the same volume of n-heptane and mix thoroughly for 5 min
- Centrifuge: 5 min / 3500 g / 10 °C (50 °F)
- Thoroughly remove the upper heptane layer, avoid transferring heptane residues into the wells
- Use 100 µL of the resulting aqueous (lower) phase per well in the assay

9.4 Bacterial cultures

Potentially toxin-forming strains of the species *S. aureus* (as well as *S. hyicus* or *S. intermedius*) have to be pre-cultivated in brain heart infusion (BHI) to ensure optimal formation of enterotoxins.

Important note: Prior to analysis, please make sure that all strains to be tested are present in pure culture.

- Centrifuge supernatants of microbiological fluid cultures for 5 min / at least 3500 g / 10 °C (50 °F)
- Sterile filtration of the supernatant required because unprecipitated or re-suspended cells may disturb the test reaction
- Use 100 µL of the filtrate per well in the assay

Remark: If measured OD values are outside the linear range (about 3.0; please note manufacturer's specifications), the filtrated culture supernatant should be further diluted with PBS buffer and analyzed again.

The cell-free culture supernatant may be stored at -20 °C (-4 °F). Thawed culture supernatants should be analyzed immediately with RIDASCREEN® SET Total and not refrozen.

9.5 Sample preparation according to Standard EN ISO 19020

9.5.1 Preparatory steps before toxin extraction

As staphylococcal enterotoxins can be heterogeneously distributed in a sample, mix the whole sample if possible, or a representative part of it, with a mixer.

Weigh 25 g \pm 0.1 g of the homogenized sample and transfer it into a beaker.

Note 1: For cheese with rind, take about 10 % rind and 90 % cheese.

Note 2: For samples in powder form, reconstitute the sample by weighing 12.5 g of sample and 12.5 g of distilled water or follow the manufacturer's instructions (e.g. milk powder as a 10 % solution).

Note 3: In the case of a suspected staphylococcal food poisoning outbreak, a minimal of 12.5 g of food should be tested.

9.5.2 Extraction of enterotoxins

Add 40 mL of warm distilled or deionized water (38 ± 2 °C / 100.4 ± 7.2 °F) to the test portion and homogenize the mixture by using a turrax or a blender. After use, rinse the mixing system with distilled water and add the rinse water to the sample being analyzed.

Note 1: For a liquid sample, do not add 40 mL of distilled water.

Shake the sample at room temperature for at least 30 min to homogeneously distribute the toxins.

Acidification step: Acidify the mixture with a few drops of hydrochloric acid in order to obtain a **pH between 3.5 and 4.0**.

Note 2: To prevent denaturation of the enterotoxins during acidification, maintain a **pH between 3.5 and 4.0**. A suitable pH-meter must be used for adjustment.

Do not allow pH to fall under 3.0. If the pH drops below 3.0 during acidification, take a fresh 25 g test portion and prepare as described in chapter 9.5.1.

Centrifuge the mixture at least at 3130 x g for 15 min at **4 °C (32.9 °F) or room temperature (20 - 25 °C / 68 - 77 °F)**. Transfer the supernatant into a beaker.

Note 3: It is recommended to rinse vessels and devices with distilled water after each step to recover a maximum of toxin.

Note 4: If the supernatant is not clear enough, centrifuge again as described above.

Note 5: The pH of the supernatant after the first centrifugation has to be < 4.5. If this is not the case, acidify to a **pH between 3.5 and 4.0** and centrifuge again.

Neutralization step: Neutralize the mixture with NaOH solution in order to obtain a **pH between 7.4 and 7.6**. Centrifuge again as described. Remove as much of the neutralized aqueous phase (supernatant) as possible and continue to use.

Note 6: To prevent denaturation of the enterotoxins during neutralization, maintain a **pH between 7.4 and 7.6**. A suitable pH-meter must be used for adjustment.

Do not allow pH to rise above 9.0. If the pH rises above 9.0 during neutralization, take a fresh 25 g test portion and prepare as described in chapter 9.5.1.

9.5.3 Extract concentration via dialysis

For each sample:

- Prepare a 30 % (w/v) PEG solution (30 g PEG / 100 mL distilled water)
- Cut about 50 to 60 cm off a dialysis membrane
- Soak the membrane in distilled water per the manufacturer's instructions (at room temperature for at least 30 minutes)
- Rinse the membrane with distilled water (outside and inside)
- Lock one end of the membrane with a closure, fill it up with the neutralized aqueous phase as prepared in 9.5.2 by using a funnel and a small piece of glass wool to keep out resuspended particles. Lock the other end of the membrane with a second closure.

Note 1: If the sample to analyze contains high amounts of salt or sugar, conduct a dialysis under agitation with 2 L of distilled water two times in one hour.

- Lay the filled dialysis membrane in a laboratory glass trough containing the 30 % (w/v) PEG solution. Allow the extracts to concentrate overnight at 5 ± 3 °C (41 ± 5.4 °F).

Note 2: If the extract is not concentrated enough, increase the dialysis time as needed. It may be necessary to add more PEG powder to the solution.

- Take the dialysis membrane out of the PEG solution and rinse the outside of the membrane with distilled water to remove all traces of PEG.

Remove the concentrated extract using:

- **PBS solution if extract contains milk or milk products.**
- **Deionized water if extract contains no milk or milk products.**

Rinse the inside of the dialysis membrane well to obtain a final concentrated extract mass between 5.0 g to 5.5 g (maximum 5.8 g for sticky extracts).

When removing the concentrate from the dialysis membrane, it is recommended:

- To rub the insides of the dialysis membrane together in order to loosen the toxins from the membrane walls and to recover the maximum amount of enterotoxins.
- To pour small drops of PBS or deionized water on the inside of the membrane multiple times and rinse repeatedly to recover all enterotoxins present.

Carefully transfer the concentrated extract into a glass vial.

Note 3: If the initial weight of the sample to be tested is lower than 25 g (9.5.1, note 3), the weight ratio of the weighed test portion to the concentrated extract should be 5:1.

For outbreaks of food poisoning or for special studies, the test portion mass may differ from the size of 25 g. The ratio of test portion mass to concentrated extract is shown in the following table:

Test portion mass	Concentrated extract
17.5 g - 25.0 g	3.5 g - 5.0 g (ratio 5:1)
12.5 g - 17.5 g	3.5 g (3.9 g max.)

Note 4: If the concentrated extract is analyzed within 48 h, store it at $5 \pm 3 \text{ }^{\circ}\text{C}$ ($41 \pm 5.4 \text{ }^{\circ}\text{F}$), otherwise store it at $\leq -18 \text{ }^{\circ}\text{C}$ ($\leq -0.4 \text{ }^{\circ}\text{F}$). The extract should be completely thawed and homogenized before testing.

10. Test procedure

10.1 Preliminary comments

Bring all reagents as well as prepared samples to room temperature ($20 - 25 \text{ }^{\circ}\text{C}$ / $68 - 77 \text{ }^{\circ}\text{F}$) before use.

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1 + 9) with distilled water (e.g. 10 mL buffer concentrate + 90 mL dist. water). Prior to dilution, dissolve any formed crystals by heating the buffer in a water bath at $37 \text{ }^{\circ}\text{C}$ ($98.6 \text{ }^{\circ}\text{F}$). The diluted and ready-to-use buffer has a shelf life of one week at $2 - 8 \text{ }^{\circ}\text{C}$ ($35 - 46 \text{ }^{\circ}\text{F}$).

Unused reagents should be immediately stored at $2 - 8 \text{ }^{\circ}\text{C}$ ($35 - 46 \text{ }^{\circ}\text{F}$).

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow wells to dry between work steps.

1. Insert a sufficient number of wells into the microwell holder for the controls and all samples. Record control and sample positions.
2. Transfer 100 μL of the controls and samples into separate wells. Use a new pipette tip for each control or sample. Cover the plate and incubate for 60 minutes at $35 - 37 \text{ }^{\circ}\text{C}$ ($95 - 98.6 \text{ }^{\circ}\text{F}$).
3. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Wash each well with 300 μL of wash buffer (10.1). Repeat the washing step four more times.

4. Pipette 100 µL of conjugate 1 into each well. Ensure that the wells are completely empty before proceeding with this step. Cover the plate and incubate for 60 minutes at 35 - 37 °C (95 - 98.6 °F).
5. Dump the liquid out of the wells into a sink. Tap the micro well holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Wash each well with 300 µL of wash buffer (10.1). Repeat the washing step four more times.
6. Pipette 100 µL of conjugate 2 into each well. Ensure that the wells are completely empty before proceeding with this step. Incubate for 30 minutes at 35 - 37 °C (95 - 98.6 °F).
7. Dump the liquid out of the wells into a sink. Tap the micro well holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Wash each well with 300 µL of wash buffer (10.1). Repeat the washing step four more times.
8. Pipette 100 µL of substrate/chromogen into each well. Ensure that the wells are completely empty before proceeding with this step. Incubate for 15 min at 35 - 37 °C (95 - 98.6 °F) in the dark.
9. A color change from red to blue indicates the presence of Staphylococcus enterotoxins in the samples. Color development tends to concentrate along the edge of the wells. Tap the microtiter plate gently to distribute the color evenly before reading.
10. Pipette 100 µL of stop solution into each well. A color change from blue to yellow indicates the presence of Staphylococcus enterotoxins in the samples.
11. Measure the extinction at 450 / 620 nm in a microtiter plate spectrophotometer. The Measurement should be done up to 10 minutes after adding the stop solution.

11. Evaluation

Special software, **RIDASOFT® Win.NET Food&Feed** (Art. No. Z9996FF), is optionally available for evaluating RIDASCREEN® enzyme immunoassays.

11.1 Quality control of the test

- The extinction value for the positive control should be greater than 1.0.
- The extinction value for the negative controls should be less than 0.2.

If these criteria have not been met, the following steps should be checked and corrected before repeating the test:

- Check kit expiry date.
- Ensure sufficient time was allowed for kit components to reach room temperature (20 - 25 °C / 68 - 77 °F).
- A new pipette tip was used for each sample or control to avoid cross-contamination.
- Check sterility of water for preparing wash buffer (use sterile water for preparation).
- Inadequate washing of wells / fewer washing steps than recommended (see chapter 10.2).
- Contaminated pipettes (clean regularly).

11.2 Result interpretation

The cut-off-value for evaluation of results as NEGATIVE or POSITIVE is calculated by adding 0.15 to the OD value of the negative control:

Cut-off-value = extinction value of negative control + 0.15

- A sample is considered POSITIVE if the test is valid (see chapter 11.1) and the extinction of the sample is greater than or equal to the cut-off-value.
- A sample is considered NEGATIVE if the test is valid (see chapter 11.1) and the extinction of the sample is less than the cut-off value-value.

11.3 Identification of enterotoxins

If samples are considered positive with the RIDASCREEN® SET Total, the single toxins can be identified by using RIDASCREEN® SET A,B,C,D,E (Art. No. R4101). Sterile filtration of samples may be necessary in some cases.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2017-08-09	Release version
2020-10-15	10.1 Preliminary comments 11 Evaluation 11.1 Quality control of the test
2022-10-18	8 Indication of instability of reagents 10.2 Test procedure 11.1 Quality control of the test
2026-05-21	Current version 9.5 Sample preparation updated Disclaimer updated

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM)
-  Manufacturer + address

Disclaimer

1. R-Biopharm AG provides a limited statutory warranty (“Gewährleistung”) as mandated and pursuant to German law against defects in title and material. For products with a defined shelf life, the warranty applies for the stated shelf life. For use cycle limited products, it applies until the stated cycle limit is reached, but not beyond the general warranty period. For all other products, the general warranty period is twelve (12) months from the date the risk transfers.
Such warranty is subject to timely notice of the alleged defect, giving sufficient detail to assess the claim. A product’s defectiveness is assessed based on the specifically agreed quality, its suitability for the contractually stated purpose, and whether all agreed-upon accessories, documentation, and instructions were provided.
2. R-Biopharm AG is not liable for issues arising from:
 - a) The failure to read, understand, or follow product instructions and information for use or safety
 - b) The use of unqualified personnel
 - c) The failure to follow industry standards such as Good Laboratory Practices
 - d) The use of unsuitable controls, sample matrices, or procedures
 - e) Other misuse
 - f) unilateral product modifications
 - g) Improper storage by others
 - h) External chemical, electromagnetic, mechanical, or electrolytic influences that exceed the defined parameters.
 - i) Harm caused by external events beyond R-Biopharm's control, including but not limited to burglary, theft, lightning, fire, water, or acts of nature.
3. R-Biopharm AG’s liability is limited to foreseeable, contract-typical losses in case of negligent breaches of obligations that are essential to achieving the agreement’s purpose and that the other party could reasonably rely on (material contract obligations). R Biopharm AG is not liable for losses resulting from non-material contract obligations.
4. ALL OTHER WARRANTIES OR GUARANTEES OF ANY KIND, EXPRESS OR IMPLIED, ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, COURSE OF DEALINGS BETWEEN THE PARTIES OR OTHERWISE. R-Biopharm AG assumes no liability for consequential damages, in particular loss of profit, interruptions, or other indirect damages.
5. The foregoing is not intended to limit liability to the extent such liability cannot be excluded or limited under applicable German law or convention, including cases of intend, willful misconduct, or gross negligence, personal injury or death, express warranties or assumptions of procurement risks, strict liability, or other non-waivable legal mandates.
6. The foregoing liability provisions apply equally to R-Biopharm AG, its directors, officers, employees, and agents.

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are registered trademarks of R-Biopharm AG.