CONGEN

SureFast[®] Speed PREP

Art. No. F1054
100 extractions

User Manual

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1 General Information

1.1 Description

The kit is intended to be used for the isolation of bacteria - DNA from food (enrichments, avulsions and swabs).

To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing.

1.2 Kit components and storage

Kit Code	Reagent / Material	Amount (per box)	Color
L	Lysis Buffer	1 x 20 ml	

The Lysis Buffer should be stored at – 20 bis + 8 °C.

1.3 Additionally required equipment and materials

- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 95°C)
- micro centrifuge (up to 12,000 rpm)
- disposal bags or waste bin
- autoclave
- biological safety cabinet

2 Protocol

2.1 Principle

1. Preparation of the basic material
2. Lysis of the basic material

2.2 Extraction protocol

2.2.1 Preparation of the basic material

Extraction of Enrichment Cultures

If the enrichment culture has been shaken or stirred, let the solids settle out for 5 to 10 min.

Transfer 500 µl of the enrichment or avulsions under sterile conditions into a 2 ml reaction tube (not provided with the kit).

Extraction of Tissues

Transfer 50 mg tissue sample into a 2.0 ml reaction tube (not provided with the kit).

2.2.2 Lysis of the basic material

Add 500 µL of Lysis Buffer (**Code L**) to the 500 µL of the enrichment culture or the 50 mg of the tissue sample. Vortex the tube briefly and put the tube in the heating block.

Incubate for 10 minutes at 95 °C without shaking.

Remove the tube from the heating block and leave the tube at room temperature for 1 minute.

Note: It is necessary to ensure that the sediment is not stirred up and no particles are pipetted in the PCR reaction.

The lysate is ready-to-use for the PCR. If the lysate is not used immediately in the PCR or is intended for a longer storage, transfer 100 µL of the lysate in a new tube (not provided with the kit) and store at -20 °C.

3 Further Information

3.1 Product Information

- Flow chart (Download: www.congen.de/en/company/downloads)
- Validation Report upon request
- Material Safety Data Sheet upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.