
CONGEN

SureFast[®] Foodborne Pathogens 4plex

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 **Content**

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	10
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup	10
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Further Information	14
3.1	Product Information	14
3.2	Technical Support	14

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1 General Information

1.1 Description

The SureFast® Foodborne Pathogens 4plex is a real-time PCR for the direct, qualitative detection and differentiation of specific DNA sequences of *Escherichia coli* virulence factors (*stx1* [subtype a-d] and *stx2* [subtype a-g]), *Listeria monocytogenes* and *Salmonella* spp..

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical verification of instruments was performed on Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Applied Biosystems 7500, Bio-Rad CFX96, R-Biopharm RIDA®CYCLER, Agilent AriaDx and Agilent Mx3005P.

1.2 Limit of Detection

The SureFast® Foodborne Pathogens 4plex real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

Note: Inconsistent mixing ratios may cause a loss of sensitivity in the low concentration channel in mixed samples especially with high amplicon concentrations (Cp value < 20).

1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria (Art. No. F1021) is recommended. To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3). For cultivation of *Escherichia coli* the method of the ISO 13136, for *Listeria monocytogenes* the method of the ISO 11290-1 and for *Salmonella* spp. the method of the ISO 10135 is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Bacteria Art. No. F1021)
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- microcentrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	<i>E. coli stx1/2/2f</i>	FAM	+	
	IAC	HEX	+	
	<i>Listeria monocytogenes</i>	ROX	+	
	<i>Salmonella</i> spp.	Cy5	+	
Agilent AriaDx	<i>E. coli stx1/2/2f</i>	FAM	+	
	IAC	HEX	+	
	<i>Listeria monocytogenes</i>	ROX	+	
	<i>Salmonella</i> spp.	Cy5	+	
Applied Biosystems 7500	<i>E. coli stx1/2/2f</i>	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
	<i>Listeria monocytogenes</i>	ROX	None	
	<i>Salmonella</i> spp.	Cy5	None	
Bio-Rad CFX96	<i>E. coli stx1/2/2f</i>	FAM	+	
	IAC	VIC/HEX	+	
	<i>Listeria monocytogenes</i>	ROX	+	
	<i>Salmonella</i> spp.	Cy5	+	
R-Biopharm RIDA®CYCLER	<i>E. coli stx1/2/2f</i>	green	+	
	IAC	yellow	+	
	<i>Listeria monocytogenes</i>	orange	+	
	<i>Salmonella</i> spp.	red	+	
Roche LightCycler® 480 II	<i>E. coli stx1/2/2f</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	533-580	+	
	<i>Listeria monocytogenes</i>	533-610	+	
	<i>Salmonella</i> spp.	618-660	+	
Roche cobas® z 480 Analyzer	<i>E. coli stx1/2/2f</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	540-580	+	
	<i>Listeria monocytogenes</i>	540-610	+	
	<i>Salmonella</i> spp.	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative *E. coli stx1/2/2f*, *Listeria monocytogenes* and *Salmonella* spp. detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative *E. coli stx1/2/2f*, *Listeria monocytogenes* and *Salmonella* spp. detection in enrichments:

For every enrichment 5 reactions for controls (1x no-template control, 1x extraction control, 1x positive control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Escherichia coli stx1/2/2f DNA is detected in the FAM-channel, *Listeria monocytogenes* DNA is detected in the ROX-channel and *Salmonella* spp. DNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-value ≤ 2 compared to the negative control. If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *E. coli stx1/2/2f*, *Listeria monocytogenes* or *Salmonella* spp. PCR assay.

Result in the respective channel				Interpretation
FAM channel <i>E. coli</i> <i>stx1/2/2f</i>	ROX channel <i>Listeria</i> <i>monocytogenes</i>	Cy5 channel <i>Salmonella</i> spp.	VIC/HEX channel IAC	
positive	negative	negative	positive/negative	<i>E. coli stx1/2/2f</i> DNA detected
negative	positive	negative	positive/negative	<i>Listeria monocytogenes</i> DNA detected
negative	negative	positive	positive/negative	<i>Salmonella</i> spp. DNA detected
negative	negative	negative	positive	negative for <i>E. coli stx1/2/2f</i> , <i>Listeria monocytogenes</i> , <i>Salmonella</i> spp.
negative	negative	negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.