

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

Trade name: **RIDASCREEN® Extractor 2**

Article number: 3055P01

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

Application of the substance / the preparation: In vitro

#### 1.3 Details of the supplier of the safety data sheet

**Manufacturer/Supplier:**

R-Biopharm AG

An der neuen Bergstraße 17

D-64297 Darmstadt

Further information obtainable from: FE e-mail: info@r-biopharm.de

1.4 Emergency telephone number: FE 0049-6151-8102-0

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008



GHS05 corrosion

Eye Dam. 1      H318 Causes serious eye damage.



GHS07

Skin Irrit. 2      H315 Causes skin irritation.

Aquatic Chronic 3      H412 Harmful to aquatic life with long lasting effects.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC



Xi; Irritant

R38-41: Irritating to skin. Risk of serious damage to eyes.

R52/53: Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

#### Information concerning particular hazards for human and environment:

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

#### Classification system:

The classification is according to the latest editions of the EU-lists, and extended by company and literature data.

#### 2.2 Label elements

**Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

**Hazard pictograms** GHS05

**Signal word** Danger

#### Hazard-determining components of labelling:

2-Mercaptoethanol

#### Hazard statements

H315 Causes skin irritation.

H318 Causes serious eye damage.

(Contd. on page 2)

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

**Trade name: RIDASCREEN® Extractor 2**

(Contd. of page 1)

H412 Harmful to aquatic life with long lasting effects.

**Precautionary statements**

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P273 Avoid release to the environment.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER/doctor.

P321 Specific treatment (see on this label).

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

**2.3 Other hazards****Results of PBT and vPvB assessment****PBT:** Not applicable.**vPvB:** Not applicable.

### SECTION 3: Composition/information on ingredients

**3.2 Chemical characterisation: Mixtures****Description:** Mixture: consisting of the following components.**Dangerous components:**

|                                    |  |        |
|------------------------------------|--|--------|
|                                    | 2-Mercaptoethanol  | 2,5-5% |
|                                    | T+ R27<br>T R25<br>C R35<br>N R51/53<br>Acute Tox. 3, H301; Acute Tox. 2, H310<br>Skin Corr. 1A, H314<br>Aquatic Chronic 2, H411 |        |
| CAS: 151-21-3<br>EINECS: 205-788-1 | sodium dodecyl sulphate  | 2,5-5% |
|                                    | Xn R21/22<br>Xi R36/38<br>Acute Tox. 4, H302; Acute Tox. 4, H312; Skin Irrit. 2, H315; Eye Irrit. 2, H319                        |        |

**Additional information:** For the wording of the listed hazard phrases refer to section 16.

### SECTION 4: First aid measures

**4.1 Description of first aid measures****General information:** Immediately remove any clothing soiled by the product.**After inhalation:**

In case of unconsciousness place patient stably in side position for transportation.

**After skin contact:** Immediately wash with water and soap and rinse thoroughly.**After eye contact:**

Rinse opened eye for several minutes under running water. Then consult a doctor.

**After swallowing:** If symptoms persist consult doctor.**4.2 Most important symptoms and effects, both acute and delayed**

No further relevant information available.

**4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

GB

(Contd. on page 3)

# Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

Trade name: RIDASCREEN® Extractor 2

(Contd. of page 2)

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### Suitable extinguishing agents:

CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

### 5.2 Special hazards arising from the substance or mixture

No further relevant information available.

### 5.3 Advice for firefighters

**Protective equipment:** No special measures required.

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Wear protective equipment. Keep unprotected persons away.

### 6.2 Environmental precautions:

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

### 6.3 Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Use neutralising agent.

Dispose contaminated material as waste according to item 13.

### 6.4 Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

Work only in fume cupboard.

**Information about fire - and explosion protection:** No special measures required.

### 7.2 Conditions for safe storage, including any incompatibilities

#### Storage:

**Requirements to be met by storerooms and receptacles:** No special requirements.

**Information about storage in one common storage facility:** Not required.

**Further information about storage conditions:** Keep container tightly sealed.

### 7.3 Specific end use(s)

No further relevant information available.

## SECTION 8: Exposure controls/personal protection

**Additional information about design of technical facilities:** No further data; see item 7.

### 8.1 Control parameters

#### Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

**Additional information:** The lists valid during the making were used as basis.

### 8.2 Exposure controls

#### Personal protective equipment:

#### General protective and hygienic measures:

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

(Contd. on page 4)

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

**Trade name: RIDASCREEN® Extractor 2**

(Contd. of page 3)

Wash hands before breaks and at the end of work.

Avoid contact with the skin.

Avoid contact with the eyes and skin.

**Respiratory protection:** Not required.**Protection of hands:**

Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

**Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

**Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

**For the permanent contact gloves made of the following materials are suitable:**

Nitrile rubber, NBR

Recommended thickness of the material:  $\geq 0.11$  mmValue for the permeation: Level  $\leq 480$ **As protection from splashes gloves made of the following materials are suitable:**

Nitrile rubber, NBR

Recommended thickness of the material:  $\geq 0.11$  mmValue for the permeation: Level  $\leq 480$ **Eye protection:**

Tightly sealed goggles

## SECTION 9: Physical and chemical properties

**9.1 Information on basic physical and chemical properties****General Information****Appearance:**

|                         |                 |
|-------------------------|-----------------|
| <b>Form:</b>            | Fluid           |
| <b>Colour:</b>          | Orange          |
| <b>Odour:</b>           | Characteristic  |
| <b>Odour threshold:</b> | Not determined. |

|                           |   |
|---------------------------|---|
| <b>pH-value at 20 °C:</b> | 6 |
|---------------------------|---|

**Change in condition**

|                                     |               |
|-------------------------------------|---------------|
| <b>Melting point/Melting range:</b> | Undetermined. |
| <b>Boiling point/Boiling range:</b> | 100 °C        |

|                     |                 |
|---------------------|-----------------|
| <b>Flash point:</b> | Not applicable. |
|---------------------|-----------------|

|                                       |                 |
|---------------------------------------|-----------------|
| <b>Flammability (solid, gaseous):</b> | Not applicable. |
|---------------------------------------|-----------------|

(Contd. on page 5)

GB

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

**Trade name: RIDASCREEN® Extractor 2**

(Contd. of page 4)

|   |   |
|---|---|
| <b>Ignition temperature:</b>                    |   |
| <b>Decomposition temperature:</b>               | Not determined.                               |
| <b>Self-igniting:</b>                           | Product is not selfigniting.                  |
| <b>Danger of explosion:</b>                     | Product does not present an explosion hazard. |
| <b>Explosion limits:</b>                        |   |
| <b>Lower:</b>                                   | Not determined.                               |
| <b>Upper:</b>                                   | Not determined.                               |
| <b>Vapour pressure at 20 °C:</b>                | 23 hPa  |
| <b>Density at 20 °C:</b>                        | 1.04397 g/cm <sup>3</sup>                     |
| <b>Relative density</b>                         | Not determined.                               |
| <b>Vapour density</b>                           | Not determined.                               |
| <b>Evaporation rate</b>                         | Not determined.                               |
| <b>Solubility in / Miscibility with water:</b>  | Fully miscible.                               |
| <b>Partition coefficient (n-octanol/water):</b> | Not determined.                               |
| <b>Viscosity:</b>                               |   |
| <b>Dynamic:</b>                                 | Not determined.                               |
| <b>Kinematic:</b>                               | Not determined.                               |
| <b>Solvent content:</b>                         |   |
| <b>Organic solvents:</b>                        | 0.0 %   |
| <b>Water:</b>                                   | 67.8 %  |
| <b>Solids content:</b>                          | 13.1 %  |
| <b>9.2 Other information</b>                    | No further relevant information available.    |

### SECTION 10: Stability and reactivity

**10.1 Reactivity** No further relevant information available.**10.2 Chemical stability****Thermal decomposition / conditions to be avoided:**

No decomposition if used according to specifications.

**10.3 Possibility of hazardous reactions** No dangerous reactions known.**10.4 Conditions to avoid** No further relevant information available.**10.5 Incompatible materials:** No further relevant information available.**10.6 Hazardous decomposition products:** No dangerous decomposition products known.

### SECTION 11: Toxicological information

**11.1 Information on toxicological effects****Acute toxicity** Based on available data, the classification criteria are not met.**LD/LC50 values relevant for classification:****57-13-6 urea**

Oral LD50 8471 mg/kg (rat)

**2-Mercaptoethanol**

Oral LD50 244 mg/kg (rat)

Dermal LD50 150 mg/kg (rabbit)

(Contd. on page 6)

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

Trade name: RIDASCREEN® Extractor 2

(Contd. of page 5)

**151-21-3 sodium dodecyl sulphate**

|      |      |                  |
|------|------|------------------|
| Oral | LD50 | 1288 mg/kg (rat) |
|------|------|------------------|

**Primary irritant effect:****Skin corrosion/irritation**

Causes skin irritation.

**Serious eye damage/irritation**

Causes serious eye damage.

**Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.**CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)****Germ cell mutagenicity** Based on available data, the classification criteria are not met.**Carcinogenicity** Based on available data, the classification criteria are not met.**Reproductive toxicity** Based on available data, the classification criteria are not met.**STOT-single exposure** Based on available data, the classification criteria are not met.**STOT-repeated exposure** Based on available data, the classification criteria are not met.**Aspiration hazard** Based on available data, the classification criteria are not met.

### SECTION 12: Ecological information

**12.1 Toxicity****Aquatic toxicity:** No further relevant information available.**12.2 Persistence and degradability** No further relevant information available.**12.3 Bioaccumulative potential** No further relevant information available.**12.4 Mobility in soil** No further relevant information available.**Ecotoxicological effects:****Remark:** Harmful to fish**Additional ecological information:****General notes:**

Water hazard class 3 (German Regulation) (Self-assessment): extremely hazardous for water

Do not allow product to reach ground water, water course or sewage system, even in small quantities.

Must not reach sewage water or drainage ditch undiluted or unneutralised.

Danger to drinking water if even extremely small quantities leak into the ground.

Harmful to aquatic organisms

**12.5 Results of PBT and vPvB assessment****PBT:** Not applicable.**vPvB:** Not applicable.**12.6 Other adverse effects** No further relevant information available.

### SECTION 13: Disposal considerations

**13.1 Waste treatment methods****Recommendation**

Must not be disposed together with household garbage. Do not allow product to reach sewage system.

Disposal must be made according to official regulations.

**Uncleaned packaging:****Recommendation:** Disposal must be made according to official regulations.**Recommended cleansing agents:** Water, if necessary together with cleansing agents.

GB

(Contd. on page 7)

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

Trade name: RIDASCREEN® Extractor 2

(Contd. of page 6)

### SECTION 14: Transport information

|  |   |
|--|---|
| <b>14.2 UN proper shipping name ADR</b>  | 2810 TOXIC LIQUID, ORGANIC, N.O.S. (THIOGLYCOL, sodium dodecyl sulphate)              |
| <b>14.5 Environmental hazards: Marine pollutant:</b>                           | No  |
| <b>14.6 Special precautions for user</b>                                       | Not applicable.   |
| <b>14.7 Transport in bulk according to Annex II of Marpol and the IBC Code</b> | Not applicable.   |
| <b>UN "Model Regulation":</b>  | UN2810, TOXIC LIQUID, ORGANIC, N.O.S. (THIOGLYCOL, sodium dodecyl sulphate), 6.1, III |

### SECTION 15: Regulatory information

**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

**Directive 2012/18/EU**

**Named dangerous substances - ANNEX I** None of the ingredients is listed.

**National regulations:**

**Waterhazard class:** Water hazard class 3 (Self-assessment): extremely hazardous for water.

**15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

**Relevant phrases**

H301 Toxic if swallowed.  
H302 Harmful if swallowed.  
H310 Fatal in contact with skin.  
H312 Harmful in contact with skin.  
H314 Causes severe skin burns and eye damage.  
H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
H411 Toxic to aquatic life with long lasting effects.

R21/22 Harmful in contact with skin and if swallowed.

R25 Toxic if swallowed.

R27 Very toxic in contact with skin.

R35 Causes severe burns.

R36/38 Irritating to eyes and skin.

R51/53 Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

**Department issuing MSDS:** QS 0049-6151-8102-0

**Contact:** QS 0049-6151-8102-0

**Abbreviations and acronyms:**

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)

ICAO: International Civil Aviation Organisation

(Contd. on page 8)

# Safety data sheet

## according to 1907/2006/EC, Article 31



Printing date 09.12.2015

Version number 10

Revision: 09.12.2015

**Trade name: RIDASCREEN® Extractor 2**

(Contd. of page 7)

ICAO-TI: Technical Instructions by the "International Civil Aviation Organisation" (ICAO)

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 3: Acute toxicity, Hazard Category 3

Acute Tox. 4: Acute toxicity, Hazard Category 4

Acute Tox. 2: Acute toxicity, Hazard Category 2

Skin Corr. 1A: Skin corrosion/irritation, Hazard Category 1A

Skin Irrit. 2: Skin corrosion/irritation, Hazard Category 2

Eye Dam. 1: Serious eye damage/eye irritation, Hazard Category 1

Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2

Aquatic Chronic 2: Hazardous to the aquatic environment - Chronic Hazard, Category 2

Aquatic Chronic 3: Hazardous to the aquatic environment - Chronic Hazard, Category 3

**\* Data compared to the previous version altered.**

GB