


**Reagent Data Sheet**

MV30F02

Version nr: 2

Version date: 07-03-2019

r-biopharm®**Reagent Data Sheet  
Corticosteroid ELISA  
ID: 5081COR**

Component ID	Description	SDS
5000PLSAR	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5081CORCO	Corticosteroid conjugate	X
	Lyophilized	
5081CORS16	Corticosteroid Standard	X
	2 - 1 - 0.5 - 0.25 - 0.125 - 0.0625 - 0 ng/ml	
5081CORAB	Anti-Corticosteroid antibody	X
	Lyophilized	
5000DBMZP	Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet