


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet  
Progesterone ELISA  
ID: 5081PROG**

Component ID	Description	SDS
5000PLSAR	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5081PROGCO	Progesterone conjugate	X
	Lyophilized	
5081PROGS	Progesterone Standard	X
	Lyophilized 5 ng/ml	
5081PROGAB	Progesterone antibody	X
	Lyophilized	
5000DBHZ	Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet