


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet****Penicillin ELISA****ID: 5091PEN**

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5091PENCO	Penicillin conjugate 100x concentrated	X
5091PENS	Penicillin Standard 3x 4 ng/ml lyophilized	X
5091PENAB	Anti-Penicillin antibody 100x concentrated	X
5091CAPDB	Dilution buffer 4x concentrated	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet